



PROCESS FILTRATION PRODUCT GUIDE

mdi
Membrane Technologies

Company Profile



36 Acre Campus



**World Class GMP Compliant
Multilocation Facilities (200,000 sq. ft.)**

Advanced Microdevices (**mdi**) is a leader in innovative membrane technologies. Starting from a single person R&D operation in 1976, **mdi** has developed into a dedicated team of 1000 plus with more than 800,000 products.

The company's core competence is its ability to develop new membrane technologies and innovate existing ones to deliver advantages to the customer for high end purification and separation applications in a wide range of industries such as pharmaceuticals, biopharmaceuticals, biotechnology, food and beverage, hospitals, and immunodiagnostics.

As membranes end up being incorporated into user friendly devices, plastic design and moulding and sealing technologies become an integral part of the chain to deliver value to the customer. Realizing this, **mdi** has grown into a vertically integrated company that helps deliver prototypes rapidly for quicker conversion to products for the market.

mdi products are used for critical applications in pharmaceutical and biopharmaceutical industries, such as sterilization of injectable drugs, sterility testing, sample preparation of drugs that are tested with highly sophisticated instrumentation, and development of new drug entities and formulations. **mdi** also offers world class membranes for making reliable immunoassays for testing of diseases at patient bedside.

mdi products meet or exceed industry standards and many of these are recognized as the best available in the world.

These products are manufactured by highly trained manpower in modern GMP facilities with large ISO class 7 production areas under ISO 9001 certified quality management system and are backed by state of the art QC testing, microbiology, reliability and validation laboratories.

A strong pipeline of new products is constantly being developed in its well equipped R&D labs.

Table of Contents

General	Page No.
➤ MembraneTechnologies	2
➤ Quality Assurance	3
➤ Validation Services	4
➤ Customer Support	5
➤ Filter Selection and Sizing	6
➤ Filter Sizing: Linear Upscaling from R&D to Production Process	7
➤ Easy Connect	8
Filters for Sterilization of Liquids: PES Membrane Filters	9
➤ Filter Selection Chart	10
➤ <i>AseptiCap KL/KS</i> - 25mm, 50mm	11
➤ <i>AseptiCap KL/KS</i> Capsule Filters	13
➤ <i>AseptiCap KS</i> Large Capsule Filters	16
➤ <i>AseptiCap KO</i> Capsule Filters	18
➤ <i>AseptiCap KSO</i> Capsule Filters	19
➤ <i>AseptiCap KSO</i> Large Capsule Filters	21
➤ <i>AseptiSure HS</i> Mini Cartridge Filters	23
➤ <i>AseptiSure HS</i> Cartridge Filters	24
➤ <i>AseptiSure HSR</i> Mini Cartridge Filters	25
➤ <i>AseptiSure HSR</i> Cartridge Filters	26
➤ <i>AseptiSure KS</i> Mini Cartridge Filters	27
➤ <i>AseptiSure KS</i> Cartridge Filters	28
Filters for Sterilization of Liquids: Hydrophilic PVDF Membrane Filters	29
➤ Filter Selection Chart	29
➤ <i>AseptiCap WS</i> Inline Capsule Filters	30
➤ <i>AseptiCap WS</i> Capsule Filters	32
➤ <i>AseptiCap WS</i> Large Capsule Filters	34
➤ <i>AseptiSure WS</i> Mini Cartridge Filters	36
➤ <i>AseptiSure WS</i> Cartridge Filters	37
Filters for Sterilization of Liquids Nylon-66 Membrane Filters	38
➤ Filter Selection Chart	39
➤ <i>AseptiCap NL/NS</i> Capsule Filters	40
➤ <i>AseptiCap NS</i> Large Capsule Filters	42
➤ <i>AseptiSure NS</i> Mini Cartridge Filters	44
➤ <i>AseptiSure NS</i> Cartridge Filters	45
➤ Nylon-66 Membrane Disc Filters	46
Filters for Sterilization of Air/Gas	47
➤ Filter Selection Chart	48
Single Use Gamma Sterilizable PVDF Membrane Filter Devices	
➤ <i>AseptiVent VF-γ</i> 25mm, 37mm, 50mm	49
➤ <i>AseptiVent VF-γ</i> Capsule Filters	51
➤ <i>AseptiVent VF-γ</i> Large Capsule Filters	53
Multiple Use Autoclavable / Steam Sterilizable PTFE Membrane Filter Devices	
➤ <i>AseptiVent TF</i> - 25mm, 37mm, 50mm	55
➤ <i>AseptiVent TF</i> Capsule Filters	57
➤ <i>AseptiVent TF</i> Large Capsule Filters	59
➤ <i>AseptiSure TH</i> Mini Cartridge Filters	61
➤ <i>AseptiSure TH</i> Cartridge Filters	62
➤ <i>AseptiSure TF</i> Mini Cartridge Filters	63
➤ <i>AseptiSure TF</i> Cartridge Filters	64
Filters for Clarification and pre-filtration	65
➤ Filter Selection Chart	66
➤ <i>ClariPro GK</i> Capsule Filters	67
➤ <i>ClariPro GK</i> Large Capsule Filters	68
➤ <i>ClariSure GK</i> Mini Cartridge Filters	69
➤ <i>ClariSure GK</i> Cartridge Filters	70
➤ <i>ClariCap GS</i> Capsule Filters	71
➤ <i>ClariCap GS</i> Large Capsule Filters	72
➤ <i>ClariSure GS</i> Mini Cartridge Filters	73
➤ <i>ClariSure GS</i> Cartridge Filters	74
➤ <i>ClariSure GP</i> Mini Cartridge Filters	75
➤ <i>ClariSure GP</i> Cartridge Filters	76
➤ <i>ClariCap PP</i> Capsule Filters	77
➤ <i>ClariCap PP</i> Large Capsule Filters	78
➤ <i>ClariSure PA</i> Mini Cartridge Filters	79
➤ <i>ClariSure PA</i> Cartridge Filters	80
➤ Microglassfiber Disc Filters	81
Filters for Polishing and Clarification	82
➤ Filter Selection Chart	83
➤ <i>ClariSure PP</i>	84
➤ <i>ClariSure PL</i>	85
➤ <i>ClariSure DP</i>	86
Chemical Compatibility	87
End Connection Availability Chart for Capsule Filters	89
Dimensions: Capsule Filters	90
Adapters and Elastomers Availability Chart for Cartridge Filters	91
Ordering Information	92

mdi offers customized solutions to enhance process efficiency, productivity, product quality, and consistency for wide range applications in biopharmaceuticals, pharmaceuticals, biotechnology, microelectronics, and immunodiagnostics.

Research and Development

A unique multidimensional research and development facility at **mdi**, involving teams specializing in critical areas of membrane technology, biotechnology, electronics, chemistry, and mechanical engineering, continuously strives with an integrated approach to develop advanced, innovative and customized products.

These filters have innovative design inputs to deliver unique performance advantages over competing products in terms of higher retention efficiencies, flow rates, throughputs, and lower filtration losses.

This unparalleled capability to develop custom products and solutions is a continuous source of pride and drives the **mdi** team to push the boundaries of technology and maximize value for the user.



mdi Microfiltration products are well designed with built-in quality assurance.

The careful selection of raw materials, validated production processes and Quality Management System certified by ISO 9001:2015 ensures consistently high quality products. mdi products meet 21 CFR, ASTM, compendia requirements and meet global regulatory expectations.

Facilities

mdi filters are produced and packaged in facilities meeting GMP requirements such as Clean Rooms with Class 10,000 and Class 100 areas for critical processes.

Deep Characterization and Certification

Apart from retention efficiency and other functional parameters such as flowrates, temperature/hydraulic stress etc., mdi filters are deeply characterized for critical areas of concern such as biosafety, bioburden levels, endotoxin levels and extractables.

Quality Control

The filters go through stringent in process and final product testing and quality is ensured by in place QMS.

Traceability

Each sterilizing grade cartridge and capsule filter has unique identification number and is accompanied with individual certificate of quality.



As per regulatory requirements, the pharmaceutical industry has to provide a high level of assurance that the sterile drug product, manufactured through aseptic processing, offers the identity, strength, quality, and purity it purports to have or is represented to possess (Ref. USFDA 21CFR 211.100(a).) Consequently it has become increasingly critical to establish/quantify the impact on the drug due to its interface with various process components under different process conditions.

Sterilizing grade filters are of critical importance in aseptic manufacturing, and it is the drug manufacturer's responsibility to show that the selected filter is able to sterilize the product under the process conditions and it does not affect the purity, quality and strength of the drug product.

Validation Services

mdi ascertain Filter validation services are designed to meet customer specific needs and help achieve regulatory compliance. These include the following:

- **Studies establishing filter integrity test values specific to drug product**
- **Filter fluid interaction studies**
 - Physico-Chemical compatibility studies
 - Extractable/Leachable studies
 - Adsorption studies
- **Microbial retention studies**
- **Throughput Studies**

All of these studies are executed as per pre-approved test methodologies to establish the test conditions and acceptance criteria.

mdi also offers post validation support for regulatory audits.

Validation Guides

mdi filters are validated as per global regulatory requirements. These filters, in compliance with the Regulation Title 21 Code of Federal Regulations (CFR) Part 314.420, have been registered at the U.S. Food and Drug Administration through Drug Master File No. DMF 15554.

Detailed documentation on validation of **mdi** filters for sterilization of fluids (air/gases and liquids) in form of Validation Guides is also available.



mdi technology executives assist in problem solving and process upgradation through experience sharing and developing customized products and systems. Some of these customer support activities are:

» Customized Filtration Solutions

mdi offers customized solutions for complex filtration problems. mdi technology executives will help you in finding solutions to filter difficult to filter fluids & minimizing filtration losses.

» Filtration System Design

Designing an efficient filtration system is an integral part of process optimization for minimizing filtration costs, increasing yields and reducing process time. mdi offers technical support for selection of filter materials by performing throughput studies to optimize filter train and filter sizing.

» Installation and Operational Qualification

mdi provides well documented installation, operational and performance qualification guidelines for all the equipment and systems it offers.

» Regulatory Assistance

mdi provides complete regulatory assistance to it's customers. mdi products and validation services meet global regulatory expectations.

» Technical Seminars

Technical seminars at customer's location are organized to help the users understand the performance characterization of product in use, differentiate between various options available in the market and select the best solutions to suit their requirements. These interactions help create optimized systems and also upgrade current processes in terms of performance and cost effectiveness.



Filter Selection and Sizing

Highly regulated process industries such as pharmaceuticals and biopharmaceuticals work with a very wide variety of process streams/fluids under different process conditions. They continuously face with the challenge of achieving their process objectives efficiently and cost effectively.

Microfiltration is a key process step to achieve critical process objectives that range from sterilization of process fluids to bioburden reduction, polishing and clarification.

These process streams can range from easy to filter water for injection to difficult to filter colloidal solutions, emulsions, liposomal drug delivery systems or large molecule high value therapeutic proteins and vaccine concentrates. Such a wide spectrum of process streams, coupled with different process objectives, is quite a challenge for the process owner.

Selection of the right filters and their sizing to meet various process needs is thus critical to successfully achieve the desired objectives.

For establishing filter type, the following important questions need to be answered with respect to the process:

1. What is the objective of filtration?
 - Sterilization
 - Bioburden reduction
 - Particle removal
 - Clarification
2. What is the fluid to be filtered?
 - Liquid or gases
 - What are the fluid properties such as pH, viscosity, temperature and surface tension
3. What will be the process conditions?
 - Will the filtration system be inline steam sterilized or autoclaved?
 - Will the system be sanitized with chemicals or hot water?
 - What are the sterilization/sanitization conditions?
 - Are the filters going to be used once or multiple times?
 - How many times the system will be sterilized /sanitized?
 - What is the maximum operating temperature?
 - What is the allowable maximum allowable differential pressure?
4. Is it going to be a continuous process or batch filtration?
5. What will be the batch volume for full scale process filtration?
6. What is the maximum allowable filtration time or the minimum desired flow rate?

Once the filter type with respect to MOC, pore size etc. has been established the next step is to establish the filter size.

Filter sizing, although to some extent is dependent on factors such as minimum desired flow rate, fluid viscosity as well as temperature, a critical parameter is the contamination profile of the fluid to be filtered. The nature and quantum of contaminants defines their interaction with the filter, which in turn defines the throughput one can achieve from a given filtration area for the fluid in question. An understanding of this behavior will help define not only the right sized filtration system but also the right combination of pre-filters and final filters to achieve desired/optimum throughput.

mdi offers filter sizing services to product development labs and process owners in full scale manufacturing. These involve small scale throughout studies to establish suitable and cost effective filtration system. Different lab scale filter and pre-filter combinations are used to maximize throughputs. The selected combination is, based on desired batch volumes or throughputs, linearly extrapolated to establish filter size.

For more information please contact our local technology executive or write to us at info@mdimembrane.com

Filter Sizing: Linear Upscaling from R&D to Production Process

Researchers in NDDS and formulation development are concerned about the impact of filter fluid interaction on stability, purity, strength etc. of the drug product. They take a keen interest in filter selection at the formulation development stage itself.

Although preliminary compatibility data support initial filter selection, for stability studies, detailed filter validations are required to provide enough documented evidence to justify specific filter use.

A critical requirement that needs to be addressed at this stage is of scalability from R&D to pilot scale to full scale production processes. Any change in filter MOC for full scale processes will require additional validation.

mdi offers a wide range of filters to provide linear scale up from lab scale to production process. While scaling up the process, the appropriate filter size can be selected by increasing the effective filtration area of filter proportionate to the process fluid volumes.

All materials of construction of core, sleeve, end caps, support layers and housing as well as manufacturing process is identical for all filter devices starting from 5 cm² to 18000 cm² hence process scaling can be facilitated without triggering additional validation studies for given process conditions. **mdi** provides complete documentation for each of the filters thereby reducing the additional validation cost and time.



25 mm, 5 cm²



50 mm, 20 cm²



1", 250 cm²



2", 500 cm²



5", 1000 cm²



8", 2000 cm²



10", 6000 cm²



Easy Connect

Wide Range of End Connections

Pharmaceutical and Biopharmaceutical processes involve transfer of high value fluids through multiple process steps. Making high quality, reliable, flexible and functionally convenient connectivity with filters is of utmost value.

mdi Capsule filters offer a wide range of reliable end connections for functional convenience and customized connectivity.

Validated for Performance

These end connections are manufactured with tight dimension tolerance and are validated for strength and connection integrity under conditions of extreme use.



3/4" Sanitary Flange



1 1/2" Sanitary Flange



1/2" HB



1/2" Single Stepped Hose Barb



1/4" SHB



Quick Connector



3/8" Hose Barb



Female Luer Lock



Male Luer Slip



1/4" MNPT

Customized Connectivity

mdi Capsule filters can also be customized to offer different inlet-outlet combinations to meet the unique connectivity needs in process assemblies. For example, stainless steel components with sanitary flange connections are sometimes required to be connected to single use disposable systems through quick-connectors or hose barb connections.



**1 1/2" Sanitary Flange
to 1/2" Hose Barb**



**1 1/2" Sanitary Flange
to 3/4" Sanitary Flange**



**High Security
1/2" Single Step
Hose Barb Connection**

Filters for Sterilization of Liquids: PES Membrane Filters

mdi produces a wide range of Gamma sterilizable and steam sterilizable PES membrane capsule and cartridge filters to meet filtration requirements of biopharmaceutical and pharmaceutical processing.

These filters meet key process requirements such as high retention efficiency, very high protein recoveries, extremely low extractables, high throughputs, wide chemical compatibility etc.

mdi PES filter devices are available as:

Filter Type	Single Layer	Multiple Layer
Gamma Sterilizable Capsule Filters	<i>AseptiCap KL-γ</i>	<i>AseptiCap KS-γ</i> <i>AseptiCap KSO-γ</i>
Autoclavable Capsule Filters	<i>AseptiCap KL</i>	<i>AseptiCap KS</i>
	<i>AseptiCap KO</i>	<i>AseptiCap KSO</i>
Steam Sterilizable Cartridge Filters	–	<i>AseptiSure KS</i>
High Temperature Resistant Steam Sterilizable Cartridge Filters	–	<i>AseptiSure HS</i> <i>AseptiSure HSR</i>

Applications

Sterile Filtration of:

- Cell culture media
- Cell culture media containing serum
- Media additives
- Final product concentrates
- Buffers
- Adjuvants
- Small Volume Parenterals
- Large Volume Parenterals
- Water for Injection

Quality Assurance

These filter devices are manufactured in Class 10,000 clean rooms under ISO 9001 : 2015 certified quality management systems and are validated to meet compendia and regulatory requirements.

Assurance	
Toxicity	Passes Biological Reactivity Test, In Vivo, as per USP <88> for Class VI plastics
Cytotoxicity	Passes Biological Reactivity Test, In Vitro, USP <87> for Cytotoxicity
Bioburden	Bioburden level is < 1000 cfu/filter device as per ISO 11737-1: 2018
Bacterial Endotoxin	Aqueous extracts exhibit < 0.25 EU/ml as established by Limulus Amebocyte Lysate (LAL) Test as per USP <85>
Non Fiber Releasing	Passes test as per USP and comply with USFDA 21 CFR Part 211.72 for fiber release
Extractables with WFI	Passes test as per USP <661>
Oxidizable Substances	Passes test as per USP <1231>
Particle Shedding	Passes test as per USP <788> for particulate matter in injections
TOC/Conductivity at 25 °C	Meets the WFI requirements of USP <643> for Total Organic Carbon and USP <645> for Water Conductivity after a specified volume of purified water flush
Indirect Food Additive	All Polypropylene components meet the FDA Indirect Food Additive requirements cited in 21 CFR 177.1520
Good Manufacturing Practice	These products are manufactured in a facility which adheres to Good Manufacturing Practices

Filter Selection Chart

Application Area	Key Application Requirements	Gamma Sterilizable Capsule Filters	Steam Sterilizable		
			Capsule Filters	Cartridge Filters	
Biopharmaceuticals					
Media preparation	- Mycoplasma removal (in case of Mammalian Cell Culture)	AseptiCap KS -γ 0.1 μ m PES Membrane Capsule Filter	AseptiCap KS 0.1 μ m PES Membrane Capsule Filter	AseptiSure KS 0.1 μ m PES Membrane Cartridge Filter	AseptiSure HS 0.1 μ m High Temperature PES Membrane Cartridge Filter
	- Microbial retention (in case of Microbial Fermentation)	AseptiCap KS -γ 0.2 μ m PES Membrane Capsule Filter	AseptiCap KS 0.2 μ m PES Membrane Capsule Filter	AseptiSure KS 0.2 μ m PES Membrane Cartridge Filter	AseptiSure HS 0.2 μ m High Temperature PES Membrane Cartridge Filter
Sterile filtration of growth regulators	- Absolute retention - Low protein binding	AseptiCap KS -γ 0.2 μ m PES Membrane Capsule Filter	AseptiCap KS 0.2 μ m PES Membrane Capsule Filter	AseptiSure KS 0.2 μ m PES Membrane Cartridge Filter	AseptiSure HS 0.2 μ m High Temperature PES Membrane Cartridge Filter
Sterile filtration of alkaline/acidic solutions for pH control	- Absolute retention - Compatible with 1-14 pH	AseptiCap KSO-γ 0.2 μ m PES Membrane Capsule Filter	AseptiCap KO/KSO 0.2 μ m PES Membrane Capsule Filter	-	AseptiSure HSR 0.2 μ m High Temperature PES Membrane Cartridge Filter
Cell Harvesting	- Bioburden reduction	AseptiCap KS -γ 0.2 μ m or 0.45 μ m PES Membrane Capsule Filter	AseptiCap KS 0.2 μ m or 0.45 μ m PES Membrane Capsule Filter	AseptiSure KS 0.2 μ m or 0.45 μ m PES Membrane Cartridge Filter	AseptiSure HS 0.2 μ m or 0.45 μ m High Temperature PES Membrane Cartridge Filter
Buffer filtration	- Bioburden reduction - Absolute retention	AseptiCap KS -γ 0.2 μ m or 0.45 μ m PES Membrane Capsule Filter	AseptiCap KS 0.2 μ m or 0.45 μ m PES Membrane Capsule Filter	AseptiSure KS 0.2 μ m or 0.45 μ m PES Membrane Cartridge Filter	AseptiSure HS 0.2 μ m or 0.45 μ m High Temperature PES Membrane Cartridge Filter
Sterile filtration of vaccines and therapeutic proteins	- Absolute retention - Low protein binding - Low holdup volume	AseptiCap KS -γ 0.2 μ m PES Membrane Capsule Filter	AseptiCap KS 0.2 μ m PES Membrane Capsule Filter	-	-
Pharmaceuticals					
Large Volume Parenterals	- Absolute retention - High throughputs	-	-	AseptiSure KS 0.2 μ m PES Membrane Cartridge Filter	AseptiSure HS 0.2 μ m High Temperature PES Membrane Cartridge Filter
Small Volume Parenterals	- Absolute retention - Low Protein Binding - Wide Chemical Compatibility	AseptiCap KS -γ 0.2 μ m PES Membrane Capsule Filter	AseptiCap KS 0.2 μ m PES Membrane Capsule Filter	AseptiSure KS 0.2 μ m PES Membrane Cartridge Filter	AseptiSure HS 0.2 μ m High Temperature PES Membrane Cartridge Filter
WFI	- Absolute retention	AseptiCap KS -γ 0.2 μ m PES Membrane Capsule Filter	AseptiCap KS 0.2 μ m PES Membrane Capsule Filter	AseptiSure KS 0.2 μ m PES Membrane Cartridge Filter	AseptiSure HS 0.2 μ m High Temperature PES Membrane Cartridge Filter

AseptiCap KL/KS Inline Capsule Filters (25mm and 50mm)

Specially designed filters for process development and formulation development labs with identical materials of construction for easy scale up to large capsule and cartridge filters.

50mm is a specially vented device, for use with peristaltic pump, to ensure easy removal of entrapped air in the upstream.

Radiation Sterilizable: AseptiCap KL/KS- γ

Autoclavable: AseptiCap KL/KS




Microbiologically Validated as per ASTM F 838
Complies with USFDA 21 CFR 210.3 (b) (6)
Meets and Exceeds USFDA 21 CFR 177.1520

Specifications


Construction			
Membrane	Hydrophilic PES		
Plastic Parts	Polypropylene		
Final Filter Pore Size	0.1 μm	0.2 μm	0.45 μm
Prefilter Pore Size (In case of AseptiCap KS)	0.2 μm , 0.45 μm	0.8 μm , 0.65 μm , 0.45 μm	0.8 μm , 0.65 μm
Integrity Testing/Retention			
Bubble Point	≥ 26 psi (1.82 Kg/cm ²) with 50% IPA ≥ 65 psi (4.56 Kg/cm ²) with Water	≥ 50 psi (3.52 Kg/cm ²) with Water	≥ 30 psi (2.11 Kg/cm ²) with Water
Microbial Retention	LRV >7 for <i>Acholeplasma laidlawii</i> (ATCC 23206) per cm ²	LRV >7 for <i>Brevundimonas diminuta</i> (ATCC 19146) per cm ²	LRV >7 for <i>Serratia marcescens</i> (ATCC 14756) per cm ²
Size			
Size	25 mm	50 mm	
Effective Filtration Area (Nominal)	5 cm ²	20 cm ²	
Operational			
Max. Operating Temperature	55 °C	60 °C	
Max. Differential Pressure	75 Psi (5 Kg/cm ²) @ 25 °C	42 Psi (3 Kg/cm ²) @ 30 °C	
Hold-up Volume(with air purge)	<50 μL	<200 μL	
Burst Pressure	> 14 Kg/cm ²	> 8 Kg/cm ²	
Sterilization	By Irradiation	AseptiCap KL/KS-γ : Gamma Irradiatable up to 50 kGy. These filters should not be autoclaved or in-line steam sterilized.	
	By Gas	AseptiCap KL/KS : Sterilizable by Ethylene Oxide	
	By Autoclave	AseptiCap KL/KS : Autoclavable at 125°C for 30 minutes, 25 cycles These cannot be In-line steam sterilized	
Shelf Life	2 years after gamma sterilization 3 years after Ethylene Oxide sterilization		
pH Compatibility	Compatible with pH range of 1-10		

Water Flow Rates

AseptiCap KL -γ 25 mm



AseptiCap KL -γ 50 mm



Ordering Information

AseptiCap KL/KS and AseptiCap KL/KS -γ, 25 mm

Type		Size		Pore Size		Inlet/Outlet		Radiation Sterilizable		X	Sterility		Pack Size	
	Code	Dia	Code		Code		Code		Code			Code		Code
AseptiCap KL	IKLX	25 mm	06	0.1 μm	36	Female Luer Lock	M	Yes	R	X	Non Sterile	1	100	04
AseptiCap KS (0.2 μm Upstream)	IKS1			0.2 μm	01	Male Luer Slip	N	No	X		EO Sterile	2		
AseptiCap KS (0.45 μm Upstream)	IKSX			0.45 μm	02	1/8" Hose Barb	H				Gamma Sterile	3		
AseptiCap KS (0.65 μm Upstream)	IKS3					1/4" Hose Barb	B							
AseptiCap KS (0.8 μm Upstream)	IKS5													

Example

IKLX	06	01	MN	R	X	1	04
------	----	----	----	---	---	---	----

AseptiCap KL/KS and AseptiCap KL/KS -γ, 50 mm

Type		Size		Pore Size		Inlet/Outlet		Radiation Sterilizable		X	Sterility		Pack Size	
	Code	Dia	Code		Code		Code		Code			Code		Code
AseptiCap KL	VKLX	50 mm	10	0.1 μm	36	1/4"-3/8" SHB	B	Yes	R	X	Non Sterile	1	10	02
AseptiCap KS (0.2 μm Upstream)	VKS1			0.2 μm	01	3/4" Sanitary Flange	S	No	X		EO Sterile	2		
AseptiCap KS (0.45 μm Upstream)	VKSX			0.45 μm	02						Gamma Sterile	3		
AseptiCap KS (0.65 μm Upstream)	VKS3													
AseptiCap KS (0.8 μm Upstream)	VKS5													

Example

VKSX	10	36	BS	X	X	1	02
------	----	----	----	---	---	---	----

Note: Gamma Sterile filters cannot be Gamma Irradiated again

Example for Non Sterile: VKSX1036BBRX102 **Example for Gamma Sterile:** VKSX1036BBXX302

For End Connection availability and dimensions with different sizes refer Pages 89-90.

AseptiCap KL/KS Small Capsule Filters (1", 2", 5", 8")

Polyethersulfone membrane capsule filters are self contained, ready to use, disposable filtration devices that contain a mini cartridge filter element sealed inside a polypropylene housing. These offer highest packing density of the membrane resulting in a very compact capsule offering long service life.

Radiation Sterilizable: AseptiCap KL/KS-γ

Autoclavable: AseptiCap KL/KS

**Microbiologically Validated
as per ASTM F 838**

**Complies with
USFDA 21 CFR 210.3 (b) (6)**

**Meets and Exceeds
USFDA 21 CFR 177.1520**




Specifications


Construction				
Membrane	Hydrophilic PES			
Support (Drainage) Layers	Polyester			
Plastic Parts	Polypropylene			
Final Filter Pore Size	0.1 µm	0.2 µm	0.45 µm	
Pre-filter Pore Size (In case of AseptiCap KS)	0.2 µm, 0.45 µm	0.8 µm, 0.65 µm, 0.45 µm		0.8 µm, 0.65 µm
Integrity Testing/Retention				
Bubble Point	≥ 26 psi (1.82 Kg/cm ²) with 50% IPA ≥ 65 psi (4.56 Kg/cm ²) with Water	≥ 50 psi (3.52 Kg/cm ²) with Water	≥ 30 psi (2.11 Kg/cm ²) with Water	
Microbial Retention	LRV >7 for <i>Acholeplasma laidlawii</i> (ATCC 23206) per cm ²	LRV >7 for <i>Brevundimonas diminuta</i> (ATCC 19146) per cm ²	LRV >7 for <i>Serratia marcescens</i> (ATCC 14756) per cm ²	
Size				
Size	1"	2"	5"	8"
Effective Filtration Area (Nominal)	250 cm ²	500 cm ²	1000 cm ²	2000 cm ²
Vent and Drain	¼" Hose Barb with Platinum Cured Silicone 'O' ring			
Operational				
Max. Operating Temperature	80 °C @ ≤ 30 psi (2 Kg/cm ²)			
Max. Differential Pressure	60 psi (4 Kg/cm ²) @ 30 °C			
Sterilization	By Irradiation	AseptiCap KL/KS -γ: Gamma Irradiatable up to 50 kGy. These filters should not be autoclaved or in-line steam sterilized.		
	By Gas	AseptiCap KL/KS: Sterilizable by Ethylene Oxide		
	By Autoclave	AseptiCap KL/KS: Autoclavable at 125°C for 30 minutes, 25 cycles These cannot be In-line steam sterilized		
Shelf Life	2 years after gamma sterilization 3 years after Ethylene Oxide sterilization			
pH Compatibility	Compatible with pH range of 1-10			

Water Flow Rates


0.1µm AseptiCap KS, 1" Capsule Filter




0.1µm AseptiCap KS, 2" Capsule Filter




0.1µm AseptiCap KS, 5" Capsule Filter




0.1µm AseptiCap KS, 8" Capsule Filter



0.2µm AseptiCap KS, 1" Capsule Filter



0.2µm AseptiCap KS, 2" Capsule Filter




End Connection Type:


- A:** ¼" Stepped Hose Barb
Q: Single Step ½" Hose Barb
E: 1½" Sanitary Flange
J: Quick Connector
S: ¾" Sanitary Flange
U: Female Luer Lock
D: ½" Hose Barb

Water Flow Rates

0.2µm AseptiCap KS, 5" Capsule Filter



0.2µm AseptiCap KS, 8" Capsule Filter



End Connection Type:

A: ¼" Stepped Hose Barb **Q:** Single Step ½" Hose Barb **E:** 1½" Sanitary Flange **J:** Quick Connector **S:** ¾" Sanitary Flange

Ordering Information

AseptiCap KL/KS and AseptiCap KL/KS -γ

Type		Size		Pore Size		Inlet/Outlet		Radiation Sterilizable		Bell		Sterility		Pack Size	
	Code		Code		Code		Code		Code		Code		Code		Code
AseptiCap KL (Single Layer)	DKLX	1"	51	0.1 µm	36	¼" SHB	A	Yes	R	Yes	B	Non Sterile	1	1	01
AseptiCap KS (0.2 µm Upstream)	DKS1	2"	52	0.2 µm	01	¼" MNPT	B	No	X	No Bell	X	EO Sterile	2		
AseptiCap KS (0.45 µm Upstream)	DKSX	5"	53	0.45 µm	02	½" MNPT	C					Gamma Sterile	3		
AseptiCap KS (0.65 µm Upstream)	DKS3	8"	57			½" Hose Barb	D								
AseptiCap KS (0.8 µm Upstream)	DKS5					1½" Sanitary Flange	E								
						¾" Sanitary Flange	S								
						Quick Connector	J								
						Single Step ½" Hose Barb	Q								
						Female Luer Lock	U								
						Male Luer Slip	W								
						⅜" Hose Barb	N								
						⅝" Hose Barb	I								

Example:

DKSX	57	36	DD	R	X	1	01
------	----	----	----	---	---	---	----

Note: Gamma Sterile filters cannot be Gamma Irradiated again

Example for Non Sterile: DKSX5136DDR01 **Example for Gamma Sterile:** DKSX5136DDX01

For End Connection, bell availability and dimensions with different sizes refer Pages 89-90.

AseptiCap KS Large Capsule Filters (5", 10", 20", 30")

These are large disposable Polyethersulfone membrane capsule filters for high value biopharma manufacturing processes, providing a unique combination of high throughputs and low hold up volumes. These capsule filters offer serial filtration incorporating a large pore size upstream membrane to protect the downstream membrane and do away with the time and expense associated with assembling, cleaning and validating stainless steel housings.



**Microbiologically Validated
as per ASTM F 838**

**Complies with
USFDA 21 CFR 210.3 (b) (6)**

**Meets and Exceeds
USFDA 21 CFR 177.1520**

Radiation Sterilizable: AseptiCap KS- γ



Autoclavable: AseptiCap KS

Specifications



Construction			
Membrane	Hydrophilic PES		
Support (Drainage) Layers	Polyester		
Plastic Parts	Polypropylene		
Final Filter Pore Size	0.1 μm	0.2 μm	0.45 μm
Prefilter Pore Size	0.2 μm , 0.45 μm	0.8 μm , 0.65 μm , 0.45 μm	0.8 μm , 0.65 μm
Integrity Testing/Retention			
Bubble Point	≥ 26 psi (1.82 Kg/cm ²) with 50% IPA ≥ 65 psi (4.56 Kg/cm ²) with Water	≥ 50 psi (3.52 Kg/cm ²) with Water	≥ 30 psi (2.11 Kg/cm ²) with Water
Max. Air Diffusion Flow for 10" Capsule Filters	≤ 29 ml/min @ 50 psi (3.52 Kg/cm ²) with Water	≤ 30 ml/min @ 37 psi (2.6 Kg/cm ²) with Water	≤ 35 ml/min @ 22 psi (1.54 Kg/cm ²) with Water
Microbial Retention	LRV >7 for <i>Acholeplasma laidlawii</i> (ATCC 23206) per cm ²	LRV >7 for <i>Brevundimonas diminuta</i> (ATCC 19146) per cm ²	LRV >7 for <i>Serratia marcescens</i> (ATCC 14756) per cm ²
Size			
Size	5"	10"	20"
Effective Filtration Area (Nominal)	3000 cm ²	6000 cm ²	12000 cm ²
Vent and Drain	$\frac{1}{4}$ " Hose Barb with platinum cured Silicone 'O' ring		
Operational			
Max. Operating Temperature	80 °C @ ≤ 30 psi (2 Kg/cm ²)		
Max. Differential Pressure	60 psi (4 Kg/cm ²) @ 30 °C		
Sterilization	By Irradiation	AseptiCap KL/KS -γ: Gamma Irradiatable up to 50 kGy. These filters should not be autoclaved or in-line steam sterilized.	
	By Gas	AseptiCap KL/KS: Sterilizable by Ethylene Oxide	
	By Autoclave	AseptiCap KL/KS: Autoclavable at 125°C for 30 minutes, 25 cycles These cannot be In-line steam sterilized	
Shelf Life	2 years after gamma sterilization 3 years after Ethylene Oxide sterilization		
pH Compatibility	Compatible with pH range of 1-10		

Water Flow Rates

0.1 µm AseptiCap KS, 10" Large Capsule Filters



0.2 µm AseptiCap KS, 10" Large Capsule Filters



End Connection Type: Q: Single Step ½" Hose Barb E: 1½" Sanitary Flange

Ordering Information

AseptiCap KS and AseptiCap KS-γ

Type		Size		Pore Size		Inlet/Outlet		Radiation Sterilizable		Inline / T-line		Sterility		Pack Size	
	Code		Code		Code		Code		Code		Code		Code		Code
AseptiCap KS (0.2 µm Upstream)	LKS1	5**	53	0.1 µm	36	1½" Sanitary Flange	E	Yes	R	Inline	X	Non Sterile	1	1	01
AseptiCap KS (0.45 µm Upstream)	LKSX	10"	54	0.2 µm	01	Single Step ½" Hose Barb	Q	No	X	T-line**	T	EO Sterile	2		
AseptiCap KS (0.65 µm Upstream)	LKS3	20"	55	0.45 µm	02	¾" Sanitary Flange	S					Gamma Sterile	3		
AseptiCap KS (0.8 µm Upstream)	LKS5	30"	56			3/8" Hose Barb	I								
						1" Hose Barb	Z								

Example:

LKSX	54	02	EE	R	T	1	01
------	----	----	----	---	---	---	----

* Size 5" is available in In-line Capsule Filters Only

** T-line Capsule Filters are available with 1½" Sanitary Flange connections Only

Note: Gamma Sterile filters cannot be Gamma Irradiated again

Example for Non Sterile: LKSX5402EERX101 Example for Gamma Sterile: LKSX5402EEXX301

For End Connection availability and dimensions with different sizes refer Pages 89-90.

AseptiCap KO Small Capsule Filters (1", 2", 5", 8")

AseptiCap KO capsule filters incorporate a low protein binding PES membrane with polypropylene drainage layers to ensure pH compatibility from 1-14 making these ideal for alkaline fluid streams.



Microbiologically Validated as per ASTM F 838
Complies with USFDA 21 CFR 210.3 (b) (6)
Meets and Exceeds USFDA 21 CFR 177.1520

Specifications

Construction				
Membrane	Hydrophilic PES			
Support (Drainage) Layers	Polypropylene			
Plastic Parts	Polypropylene			
Filter Pore Size	0.2 µm	0.45 µm		
Integrity Testing/Retention				
Bubble Point	≥ 50 psi (3.52 Kg/cm ²) with Water	≥ 30 psi (2.11 Kg/cm ²) with Water		
Microbial Retention	LRV >7 for <i>Brevundimonas diminuta</i> (ATCC 19146) per cm ²	LRV >7 for <i>Serratia marcescens</i> (ATCC 14756) per cm ²		
Size				
Size	1"	2"	5"	8"
Effective Filtration Area (Nominal)	250 cm ²	500 cm ²	1000 cm ²	2000 cm ²
Vent and Drain	1/4" Hose Barb with platinum cured Silicone 'O' ring			
Operational				
Max. Operating Temperature	80 °C @ ≤30 psi (2 Kg/cm ²)			
Max. Differential Pressure	60 psi (4 Kg/cm ²) @ 30 °C			
Sterilization	By Gas	Sterilizable by Ethylene Oxide		
	By Autoclave	Autoclavable at 125 °C for 30 minutes, 25 cycles. Cannot be in-line steam sterilized		
Typical Water Flow Rates (0.2 µm, 8")	11 lpm @ 0.70 Kg/cm ² @ 27 °C			
Shelf Life	3 years after Ethylene Oxide sterilization			
pH Compatibility	Compatible with pH range of 1-14			

Ordering Information

Type		Size		Pore Size		Inlet/Outlet		X	X	Sterility		Pack Size	
	Code		Code		Code		Code				Code		Code
AseptiCap KO	DKLO	1"	51	0.2 µm	01	1/4" SHB	A			Non Sterile	1	1	01
		2"	52	0.45 µm	02	1/2" Hose Barb	D			EO Sterile	2		
		5"	53			1 1/2" Sanitary Flange	E						
		8"	57			3/4" Sanitary Flange	S						
						Quick Connector	J						
						Single Step 1/2" Hose Barb	Q						
						Female Luer Lock	U						
						Male Luer Slip	W						
						3/16" Hose Barb	N						
						3/8" Hose Barb	I						

Example:

DKLO	57	01	DD	X	X	1	01
------	----	----	----	---	---	---	----

For End Connection availability and dimensions with different sizes refer Pages 89-90.

AseptiCap KSO Small Capsule Filters (1", 2", 5", 8")

mdi AseptiCap KSO are Polyethersulfone membrane capsule filters offering wide pH (1-14) compatibility. These filters are specially designed for alkaline fluid streams in bio-pharma manufacturing processes, with added advantages of high throughputs and low hold up volumes.

These capsule filters offer serial filtration incorporating a large pore size upstream membrane to protect the downstream membrane for enhanced throughputs.

AseptiCap KSO are validated for use in pharmaceutical and biopharmaceutical applications.

Radiation Sterilizable: AseptiCap KSO- γ

Autoclavable: AseptiCap KSO



**Microbiologically Validated
as per ASTM F 838**

**Complies with
USFDA 21 CFR 210.3 (b) (6)**


**Meets and Exceeds
USFDA 21 CFR 177.1520**

Specifications


Construction				
Membrane	Hydrophilic PES			
Support (Drainage) Layers	Polypropylene			
Plastic Parts	Polypropylene			
Final Filter Pore Size	0.2 μm		0.45 μm	
Pre-filter Pore Size	0.8 μm , 0.45 μm		0.8 μm	
Integrity Testing/Retention				
Bubble Point	≥ 50 psi (3.52 Kg/cm ²) with Water		≥ 30 psi (2.11 Kg/cm ²) with Water	
Microbial Retention	LRV >7 for <i>Brevundimonas diminuta</i> (ATCC 19146) per cm ²		LRV >7 for <i>Serratia marcescens</i> (ATCC 14756) per cm ²	
Size				
Size	1"	2"	5"	8"
Effective Filtration Area (Nominal)	250 cm ²	500 cm ²	1000 cm ²	2000 cm ²
Vent and Drain	1/4" Hose Barb with platinum cured Silicone 'O' ring			
Operational				
Max. Operating Temperature	80 °C @ ≤ 30 psi (2 Kg/cm ²)			
Max. Differential Pressure	60 psi (4 Kg/cm ²) @ 30 °C			
Sterilization	By Irradiation	AseptiCap KSO -γ: Gamma Irradiatable up to 50 kGy. These filters should not be autoclaved or in-line steam sterilized.		
	By Gas	AseptiCap KSO: Sterilizable by Ethylene Oxide		
	By Autoclave	AseptiCap KSO: Autoclavable at 125°C for 30 minutes, 25 cycles These cannot be In-line steam sterilized		
Shelf Life	2 years after Gamma Sterilization 3 years after Ethylene Oxide Sterilization			
pH Compatibility	Compatible with pH range of 1-14			

Water Flow Rates


0.2µm AseptiCap KSO, 1" Capsule Filter




0.2µm AseptiCap KSO, 2" Capsule Filter



0.2µm AseptiCap KSO, 5" Capsule Filter



0.2µm AseptiCap KSO, 8" Capsule Filter



End Connection Type:

A: ¼" Stepped Hose Barb **Q:** Single Step ½" Hose Barb **E:** 1½" Sanitary Flange **J:** Quick Connector **S:** ¾" Sanitary Flange **J:** Quick Connector

Ordering Information

AseptiCap KSO and AseptiCap KSO-γ

Type		Size		Pore Size		Inlet/Outlet		Radiation Sterilizable		X	Sterility		Pack Size		
	Code		Code		Code		Code	Yes	Code	No		Code		Code	
AseptiCap KSO (0.45 µm Upstream)	DKOX	1"	51	0.2 µm	01	¼" SHB	A	Yes	R	No	X	Non Sterile	1	1	01
AseptiCap KSO (0.8 µm Upstream)	DKO5	2"	52	0.45 µm	02	½" Hose Barb	D					EO Sterile	2		
		5"	53			1½" Sanitary Flange	E					Gamma Sterile	3		
		8"	57			¾" Sanitary Flange	S								
						Quick Connector	J								
						Single Step ½" Hose Barb	Q								
						Female Luer Lock	U								
						Male Luer Slip	W								
						¾" Hose Barb	N								
						⅝" Hose Barb	I								

Example:

DKOX	52	01	EE	R	X	1	01
------	----	----	----	---	---	---	----

Note: Gamma Sterile filters cannot be Gamma Irradiated again

Example for Non Sterile: DKOX5202EERX101

Example for Gamma Sterile: DKOX5202EEXX301

For End Connection availability and dimensions with different sizes refer Pages 89-90.

AseptiCap KSO Large Capsule Filters (5", 10", 20", 30")

mdi AseptiCap KSO are Polyethersulfone membrane capsule filters offering wide pH (1-14) compatibility. These filters are specially designed for alkaline fluid streams in biopharma manufacturing processes, with added advantages of high throughputs and low hold up volumes.

These capsule filters offer serial filtration incorporating a large pore size upstream membrane to protect the downstream membrane for enhanced throughputs.

AseptiCap KSO are validated for use in pharmaceutical and biopharmaceutical applications.

Radiation Sterilizable: AseptiCap KSO-γ

Autoclavable: AseptiCap KSO



**Microbiologically Validated
as per ASTM F 838**

**Complies with
USFDA 21 CFR 210.3 (b) (6)**



**Meets and Exceeds
USFDA 21 CFR 177.1520**

Specifications

Construction				
Membrane	Hydrophilic PES			
Support (Drainage) Layers	Polypropylene			
Plastic Parts	Polypropylene			
Final Filter Pore Size	0.2 µm		0.45 µm	
Pre-filter Pore Size	0.8 µm, 0.45µm		0.8 µm	
Integrity Testing/Retention				
Bubble Point	≥ 50 psi (3.52 Kg/cm ²) with Water		≥ 30 psi (2.11 Kg/cm ²) with Water	
Max. Air Diffusion Flows for 10" Capsule Filters	≤ 30 ml/min @ 37 psi (2.6 Kg/cm ²) with Water		≤ 35 ml/min @ 22 psi (1.54 Kg/cm ²) with Water	
Microbial Retention	LRV >7 for <i>Brevundimonas diminuta</i> (ATCC 19146) per cm ²		LRV >7 for <i>Serratia marcescens</i> (ATCC 14756) per cm ²	
Size				
Size	5"	10"	20"	30"
Effective Filtration Area (Nominal)	3000 cm ²	6000 cm ²	12000 cm ²	18000 cm ²
Vent and Drain	¼" Hose Barb with platinum cured Silicone 'O' ring			
Operational				
Max. Operating Temperature	80 °C @ ≤ 30 psi (2 Kg/cm ²)			
Max. Differential Pressure	60 psi (4 Kg/cm ²) @ 30 °C			
Sterilization	By Irradiation	AseptiCap KSO -γ: Gamma Irradiatable up to 50 kGy. These filters should not be autoclaved or in-line steam sterilized.		
	By Gas	AseptiCap KSO: Sterilizable by Ethylene Oxide		
	By Autoclave	AseptiCap KSO: Autoclavable at 125°C for 30 minutes, 25 cycles These cannot be In-line steam sterilized		
Shelf Life	2 years after Gamma Sterilization 3 years after Ethylene Oxide Sterilization			
pH Compatibility	Compatible with pH range of 1-14			

Water Flow Rates

0.2 µm AseptiCap KSO, 10" Large Capsule Filters



End Connection Type: Q: Single Step ½" Hose Barb E: 1½" Sanitary Flange

Ordering Information

AseptiCap KSO and AseptiCap KSO-γ

Type		Size		Pore Size		Inlet/Outlet		Radiation Sterilizable		Inline / T-line		Sterility		Pack Size	
	Code		Code		Code		Code	Yes	Code	Inline	Code		Code		Code
AseptiCap KSO (0.45 µm Upstream)	LKOX	5"	53	0.2 µm	01	1½" Sanitary Flange	E	Yes	R	Inline	X	Non Sterile	1	1	01
AseptiCap KSO (0.8 µm Upstream)	LKO5	10"	54	0.45 µm	02	Single Step ½" Hose Barb	Q	No	X	T-line	T	EO Sterile	2		
		20"	55			¾" Sanitary Flange	S					Gamma Sterile	3		
		30"	56			3/8" Hose Barb	I								
						1" Hose Barb	Z								

Example:

LKOX	54	01	EE	R	T	1	01
------	----	----	----	---	---	---	----

* Size 5" is available in In-line Capsule Filters Only

Note: Gamma Sterile filters cannot be Gamma Irradiated again

Example for Non Sterile: LKOX5402EERX101 **Example for Gamma Sterile:** LKOX5402EEXX301

For End Connection availability and dimensions with different sizes refer Pages 89-90.

AseptiSure HS Mini Cartridge filters

mdi Polyethersulfone (PES) membrane mini cartridge filters type *AseptiSure HS* are high temperature resistant filtration devices. These are designed to withstand high pressure steam sterilization upto 135°C.

AseptiSure HS mini cartridge filters with Polyethersulfone membrane serial layers offer enhanced throughputs, thus ensuring better economics.

These are validated for key performance parameters such as retention efficiency, chemical compatibility, extractables, heat stability and flow rates.



**Microbiologically Validated
as per ASTM F 838**

**Complies with
USFDA 21 CFR 210.3 (b) (6)**

**Meets and Exceeds
USFDA 21 CFR 177.1520**

Specifications

Construction			
Membrane	Hydrophilic PES		
Support (Drainage) Layers	Polyester		
Plastic Parts	Polypropylene		
Final Filter Pore Size	0.1 µm	0.2 µm	0.45 µm
Prefilter Pore Size	0.2 µm, 0.45 µm	0.8 µm, 0.65 µm, 0.45 µm	0.8 µm, 0.65 µm
Integrity Testing/Retention			
Bubble Point	≥ 26 psi (1.82 Kg/cm ²) with 50% IPA ≥ 65 psi (4.56 Kg/cm ²) with Water	> 50psi (3.52Kg/cm ²) with Water	> 30psi (2.11Kg/cm ²) with Water
Microbial Retention	LRV >7 for <i>Acholeplasma laidlawii</i> (ATCC 23206) per cm ²	LRV >7 for <i>Brevundimonas diminuta</i> (ATCC 19146) per cm ²	LRV >7 for <i>Serratia marcescens</i> (ATCC 14756) per cm ²
Size			
Size	2.5"	5"	
Effective Filtration Area (Nominal)	1000 cm ²	2000 cm ²	
Operational			
Water Flow Rate for 0.2µm @ 0.70kg/cm ² @27°C	9 lpm	15 lpm	
Max. Operating Temperature	80 °C @ ≤ 2 Kg/cm ² (30 psi)		
Max. Differential Pressure	50 psi (3.5 Kg/cm ²) @ 25 °C		
Reverse Pressure	≤ 10 psi (0.7 Kg/cm ²) @ 25 °C		
Sterilization	In-line steam sterilizable upto 135 °C for 30 minutes at a maximum differential pressure of 5 psi (0.35 Kg/cm ²), 25 cycles		
pH Compatibility	Compatible with pH range of 1-10		

Ordering Information

Type		Size		Pore Size		Adapter		Elastomer		Sterility		Pack Size	
	Code		Code		Code		Code		Code		Code		Code
<i>AseptiSure HS</i> (0.2 µm Upstream)	CPH1	2.5"	50	0.1 µm	36	4463	E0	Silicone	SS	Non Sterile	1	1	01
<i>AseptiSure HS</i> (0.45 µm Upstream)	CPHX	5"	53	0.2 µm	01	4463B	H0						
<i>AseptiSure HS</i> (0.65 µm Upstream)	CPH3			0.45 µm	02	4440	U0						
<i>AseptiSure HS</i> (0.8 µm Upstream)	CPH5					Seal-K	G0*						
						Seal-O	F0						
						Seal-M	J0						

***G0 adapter is not available with elastomer. Please mention XX in place of elastomer code while ordering**

Example:

CPHX	50	01	E0	SS	1	01
------	----	----	----	----	---	----

For Adapters and Elastomers availability refer Page 91.

AseptiSure HS Standard Cartridge Filters

**Microbiologically Validated
as per ASTM F 838**

**Complies with
USFDA 21 CFR 210.3 (b) (6)**

**Meets and Exceeds
USFDA 21 CFR 177.1520**

AseptiSure HS high temperature resistant, serial layer Polyethersulfone cartridge filters are designed to withstand high pressure differential upto 0.3 Kg/cm² (5 psi) at high steam sterilization temperatures of upto 135 °C.



Specifications

Construction				
Membrane	Hydrophilic PES			
Support (Drainage) Layers	Polyester			
Plastic Parts	Polypropylene			
Final Filter Pore Size	0.1 µm	0.2 µm	0.45 µm	
Prefilter Pore Size	0.2 µm, 0.45 µm	0.8 µm, 0.65 µm, 0.45 µm	0.8 µm, 0.65 µm	
Integrity Testing/Retention				
Bubble Point	≥ 26 psi (1.82 Kg/cm ²) with 50% IPA ≥ 65 psi (4.56 Kg/cm ²) with Water	≥ 50 psi (3.52 Kg/cm ²) with Water	≥ 30 psi (2.11 Kg/cm ²) with Water	
Air Diffusion Flow (10" Cartridge Filter)	≤ 29 ml/min @ 50 psi (3.52 Kg/cm ²) with Water	≤ 30 ml/min @ 37 psi (2.6 Kg/cm ²) with Water	≤ 35 ml/min @ 22 psi (1.54 Kg/cm ²) with Water	
Microbial Retention	LRV >7 for <i>Acholeplasma laidlawii</i> (ATCC 23206) per cm ²	LRV >7 for <i>Brevundimonas diminuta</i> (ATCC 19146) per cm ²	LRV >7 for <i>Serratia marcescens</i> (ATCC 14756) per cm ²	
Size				
Size	5"	10"	20"	30"
Effective Filtration Area (Nominal)	3000 cm ²	6000 cm ²	12000 cm ²	18000 cm ²
Operational				
Typical Water Flow Rates (for 0.2µm @ 0.70 Kg/cm ² @ 27 °C)	25 lpm	45 lpm	81 lpm	113 lpm
Max. Operating Temperature	80 °C @ ≤ 30psi (2 Kg/cm ²)			
Max. Differential Pressure	50 psi (3.5 Kg/cm ²) @ 25 °C			
Reverse Pressure	≤ 10 psi (0.7 Kg/cm ²) @ 25 °C			
Sterilization	In-line steam sterilizable upto 135 °C for 30 minutes at a maximum differential pressure of 5 psi (0.35 Kg/cm ²), 25 cycles			
pH Compatibility	Compatible with pH range of 1-10			

Ordering Information

Type		Size		Pore Size		Adapter		Elastomer		Sterility		Pack Size	
	Code		Code		Code		Code		Code		Code		Code
AseptiSure HS (0.2 µm Upstream)	CPH1	5"	53	0.1 µm	36	7P	A0	Silicone	SS	Non Sterile	1	1	01
AseptiSure HS (0.45 µm Upstream)	CPHX	10"	54	0.2 µm	01	7P without fin	A1	Viton	SV				
AseptiSure HS (0.65 µm Upstream)	CPH3	20"	55	0.45 µm	02	28 with fin	C0	EPDM	SE				
AseptiSure HS (0.8 µm Upstream)	CPH5	30"	56			'O'	D0	FEP Encapsulated Viton	FV				

Example:

CPHX	54	01	A0	SS	1	01
-------------	-----------	-----------	-----------	-----------	----------	-----------

For Adapters and Elastomers availability refer Page 91.

AseptiSure HSR Mini Cartridge Filter

mdi Polyethersulfone (PES) membrane mini cartridge filters type *AseptiSure HSR* are designed to withstand high pressure differential upto 0.3 Kg/cm² (5 psi) at high steam sterilization temperatures of upto 135 °C.

These filters with Polyethersulfone membrane and Polypropylene support layers offer pH compatibility from 1 to 14, and are thus ideal for use with acidic as well as alkaline solutions.

mdi *AseptiSure HSR* Mini Cartridge filters are validated for key performance parameters such as retention efficiency, chemical compatibility, extractables, heat stability and flow rates.



**Microbiologically Validated
as per ASTM F 838**

**Complies with
USFDA 21 CFR 210.3 (b) (6)**

**Meets and Exceeds
USFDA 21 CFR 177.1520**

Specifications

Construction			
Membrane	Hydrophilic PES		
Support (Drainage) Layers	Polypropylene		
Plastic Parts	Polypropylene		
Final Filter Pore Size	0.1 µm	0.2 µm	0.45 µm
Prefilter Pore Size	0.2 µm, 0.45 µm	0.8 µm, 0.65 µm, 0.45 µm	0.8 µm, 0.65 µm
Integrity Testing/Retention			
Bubble Point	≥ 26 psi (1.82 Kg/cm ²) with 50% IPA ≥ 65 psi (4.56 Kg/cm ²) with Water	≥ 50 psi (3.52 Kg/cm ²) with Water	≥ 30 psi (2.11 Kg/cm ²) with Water
Microbial Retention	LRV >7 for <i>Acholeplasma laidlawii</i> (ATCC 23206) per cm ²	LRV >7 for <i>Brevundimonas diminuta</i> (ATCC 19146) per cm ²	LRV >7 for <i>Serratia marcescens</i> (ATCC 14756) per cm ²
Size			
Size	2.5"		5"
Effective Filtration Area (Nominal)	1000 cm ²		2000 cm ²
Operational			
Typical Water Flow Rates (for 0.2µm @ 0.70 Kg/cm ² @ 27 °C)	9 lpm		15 lpm
Max. Operating Temperature	80 °C @ ≤ 30psi (2 Kg/cm ²)		
Max. Differential Pressure	50 psi (3.5 Kg/cm ²) @ 25 °C		
Reverse Pressure	≤ 10 psi (0.7 Kg/cm ²) @ 25 °C		
Sterilization	In-line steam sterilizable upto 135 °C for 30 minutes at a maximum differential pressure of 5 psi (0.35 Kg/cm ²), 25 cycles		
pH Compatibility	Compatible with pH range of 1-14		

Ordering Information

Type		Size		Pore Size		Adapter		Elastomer		Sterility		Pack Size	
	Code		Code		Code		Code		Code		Code		Code
<i>AseptiSure HSR</i> (0.2 µm Upstream)	CHR1	2.5"	50	0.1 µm	36	4463	E0	Silicone	SS	Non Sterile	1	1	01
<i>AseptiSure HSR</i> (0.45 µm Upstream)	CHRX	5"	53	0.2 µm	01	4463B	H0						
<i>AseptiSure HSR</i> (0.65 µm Upstream)	CHR3			0.45 µm	02	4440	U0						
<i>AseptiSure HSR</i> (0.8 µm Upstream)	CHR5					Seal-K	G0*						
						Seal-O	F0						
						Seal-M	J0						

***G0 adapter is not available with elastomer. Please mention XX in place of elastomer code while ordering**

Example:

CHRX	53	01	E0	SS	1	01
-------------	-----------	-----------	-----------	-----------	----------	-----------

For Adapters and Elastomers availability refer Page 91.

AseptiSure HSR Standard Cartridge Filters

mdi Polyethersulfone (PES) membrane cartridge filters type *AseptiSure HSR* are high temperature resistant filtration devices. These are designed to withstand high pressure differential at high steam sterilization temperature upto 135°C. These filters exhibit high mechanical stability, and wide chemical compatibility even with alkaline process fluids.

These filters come with Polyethersulfone membrane serial layers and Polypropylene support layers to offer 1-14 pH compatibility.

mdi *AseptiSure HSR* cartridge filters are validated for key performance parameters such as retention efficiency, chemical compatibility, extractables, heat stability and flow rates.



Microbiologically Validated as per ASTM F 838
Complies with USFDA 21 CFR 210.3 (b) (6)
Meets and Exceeds USFDA 21 CFR 177.1520

Specifications

Construction				
Membrane	Hydrophilic PES			
Support (Drainage) Layers	Polypropylene			
Body and Core	Polypropylene			
Final Filter Pore Size	0.1 µm	0.2 µm	0.45 µm	
Prefilter Pore Size	0.2 µm, 0.45 µm	0.8 µm, 0.65 µm, 0.45 µm		0.8 µm, 0.65 µm
Integrity Testing/Retention				
Bubble Point	≥ 26 psi (1.82 Kg/cm ²) with 50% IPA ≥ 65 psi (4.56 Kg/cm ²) with Water	≥ 50 psi (3.52 Kg/cm ²) with Water	≥ 30 psi (2.11 Kg/cm ²) with Water	
Air Diffusion Flow (10" Cartridge Filter)	≤ 29 ml/min @ 50 psi (3.52 Kg/cm ²) with Water	≤ 30 ml/min @ 37 psi (2.6 Kg/cm ²) with Water	≤ 35 ml/min @ 22 psi (1.54 Kg/cm ²) with Water	
Microbial Retention	LRV >7 for <i>Acholeplasma laidlawii</i> (ATCC 23206) per cm ²	LRV >7 for <i>Brevundimonas diminuta</i> (ATCC 19146) per cm ²	LRV >7 for <i>Serratia marcescens</i> (ATCC 14756) per cm ²	
Size				
Size	5"	10"	20"	30"
Effective Filtration Area (Nominal)	3000 cm ²	6000 cm ²	12000 cm ²	18000 cm ²
Operational				
Typical Water Flow Rates (for 0.2µm @ 0.70 Kg/cm ² @ 27 °C)	25 lpm	45 lpm	81 lpm	113 lpm
Max. Operating Temperature	80 °C @ ≤ 30 psi (2 Kg/cm ²)			
Max. Differential Pressure	50 psi (3.5 Kg/cm ²) @ 25 °C			
Reverse Pressure	≤ 10 psi (0.7 Kg/cm ²) @ 25 °C			
Sterilization	In-line steam sterilizable upto 135 °C for 30 minutes at a maximum differential pressure of 5 psi (0.35 Kg/cm ²), 25 cycles			
pH Compatibility	Compatible with pH range of 1-14			

Ordering Information

Type		Size		Pore Size		Adapter		Elastomer		Sterility		Pack Size	
	Code		Code		Code		Code		Code		Code		Code
<i>AseptiSure HSR</i> (0.2 µm Upstream)	CHR1	5"	53	0.1 µm	36	7P	A0	Silicone	SS	Non Sterile	1	1	01
<i>AseptiSure HSR</i> (0.45 µm Upstream)	CHRX	10"	54	0.2 µm	01	7P without fin	A1	Viton	SV				
<i>AseptiSure HSR</i> (0.65 µm Upstream)	CHR3	20"	55	0.45 µm	02	28 with fin	C0	EPDM	SE				
<i>AseptiSure HSR</i> (0.8 µm Upstream)	CHR5	30"	56			'O'	D0	FEP Encapsulated Viton	FV				

Example:

CHRX	54	01	A0	SS	1	01
-------------	-----------	-----------	-----------	-----------	----------	-----------

For Adapters and Elastomers availability refer Page 91.

AseptiSure KS Mini Cartridge filters

mdi Polyethersulfone (PES) membrane mini cartridge filters type *AseptiSure KS* are serial filtration devices with a larger pore size upstream PES membrane layer to protect the downstream final PES membrane layer from premature clogging and to give enhanced throughputs, thus resulting in better economics.

mdi *AseptiSure KS* filters are validated for key performance parameters such as retention efficiency, chemical compatibility, extractables, heat stability and flow rates. These are available in a variety of pore sizes to suit specific microfiltration needs in critical and specialized process applications.

**Microbiologically Validated
as per ASTM F 838**
**Complies with
USFDA 21 CFR 210.3 (b) (6)**
**Meets and Exceeds
USFDA 21 CFR 177.1520**



Specifications

Construction			
Membrane	Hydrophilic PES		
Support (Drainage) Layers	Polyester		
Plastic Parts	Polypropylene		
Final Filter Pore Size	0.1 µm	0.2 µm	0.45 µm
Pre-filter Pore Size	0.2 µm, 0.45 µm	0.8 µm, 0.65 µm, 0.45 µm	0.8 µm, 0.65 µm
Integrity Testing/Retention			
Bubble Point	≥ 26 psi (1.82 Kg/cm ²) with 50% IPA ≥ 65 psi (4.56 Kg/cm ²) with Water	> 50 psi (3.52Kg/cm ²) with Water	> 30 psi (2.11Kg/cm ²) with Water
Microbial Retention	LRV >7 for <i>Acholeplasma laidlawii</i> (ATCC 23206) per cm ²	LRV >7 for <i>Brevundimonas diminuta</i> (ATCC 19146) per cm ²	LRV >7 for <i>Serratia marcescens</i> (ATCC 14756) per cm ²
Size			
Size	2.5"		5"
Effective Filtration Area (Nominal)	1000 cm ²		2000 cm ²
Operational			
Water Flow Rate for 0.2µm @ 0.70kg/cm ² @27°C	9 lpm		15 lpm
Max. Operating Temperature	80 °C @ ≤ 2 Kg/cm ² (30 psi)		
Max. Differential Pressure	50 psi (3.5 Kg/cm ²) @ 25 °C		
Reverse Pressure	≤ 10 psi (0.7 Kg/cm ²) @ 25 °C		
Sterilization	In-line steam sterilizable at 121°C for 30 minutes at a maximum differential pressure of 3 psi (0.21 kg/cm ²) , 25 cycles		
pH Compatibility	Compatible with pH range of 1-10		

Ordering Information

Type		Size		Pore Size		Adapter		Elastomer		Sterility		Pack Size	
	Code		Code		Code		Code		Code		Code		Code
<i>AseptiSure KS</i> (0.2 µm Upstream)	CPK1	2.5"	50	0.1 µm	36	4463	E0	Silicone	SS	Non Sterile	1	1	01
<i>AseptiSure KS</i> (0.45 µm Upstream)	CPKX	5"	53	0.2 µm	01	4463B	H0						
<i>AseptiSure KS</i> (0.65 µm Upstream)	CPK3			0.45 µm	02	4440	U0						
<i>AseptiSure KS</i> (0.8 µm Upstream)	CPK5					Seal-K	G0*						
						Seal-O	F0						
						Seal-M	J0	*G0 adapter is not available with elastomer. Please mention XX in place of elastomer code while ordering					

Example:

CPKX	50	01	E0	SS	1	01
------	----	----	----	----	---	----

For Adapters and Elastomers availability refer Page 91.

AseptiSure KS Standard Cartridge Filters

AseptiSure KS serial filtration Polyethersulfone cartridge filters incorporate a large pore size upstream membrane layer to protect the downstream terminal filtration membrane layer.



**Microbiologically Validated
as per ASTM F 838**

**Complies with
USFDA 21 CFR 210.3 (b) (6)**

**Meets and Exceeds
USFDA 21 CFR 177.1520**

Specifications

Construction				
Membrane	Hydrophilic PES			
Support (Drainage) Layers	Polyester			
Plastic Parts	Polypropylene			
Final Filter Pore Size	0.1 µm	0.2 µm	0.45 µm	
Prefilter Pore Size	0.2 µm, 0.45 µm	0.8 µm, 0.65 µm, 0.45 µm		0.8 µm, 0.65 µm
Integrity Testing/Retention				
Bubble Point	≥ 26 psi (1.82 Kg/cm ²) with 50% IPA ≥ 65 psi (4.56 Kg/cm ²) with Water	≥ 50 psi (3.52 Kg/cm ²) with Water	≥ 30 psi (2.11 Kg/cm ²) with Water	
Air Diffusion Flow (10" Cartridge Filter)	≤ 29 ml/min @ 50 psi (3.52 Kg/cm ²) with Water	≤ 30 ml/min @ 37 psi (2.6 Kg/cm ²) with Water	≤ 35 ml/min @ 22 psi (1.54 Kg/cm ²) with Water	
Microbial Retention	LRV >7 for <i>Acholeplasma laidlawii</i> (ATCC 23206) per cm ²	LRV >7 for <i>Brevundimonas diminuta</i> (ATCC 19146) per cm ²	LRV >7 for <i>Serratia marcescens</i> (ATCC 14756) per cm ²	
Size				
Size	5"	10"	20"	30"
Effective Filtration Area (Nominal)	3000 cm ²	6000 cm ²	12000 cm ²	18000 cm ²
Operational				
Typical Water Flow Rates (for 0.2µm @ 0.70 Kg/cm ² @ 27 °C)	25 lpm	45 lpm	81 lpm	113 lpm
Max. Operating Temperature	80 °C @ ≤ 30 psi (2 Kg/cm ²)			
Max. Differential Pressure	50 psi (3.5 Kg/cm ²) @ 25 °C			
Reverse Pressure	≤ 10 psi (0.7 Kg/cm ²) @ 25 °C			
Sterilization	In-line steam sterilizable at 121 °C for 30 minutes at a maximum differential pressure of 3 psi (0.21 kg/cm ²) , 25 cycles			
pH Compatibility	Compatible with pH range of 1-10			

Ordering Information

Type		Size		Pore Size		Adapter		Elastomer		Sterility		Pack Size	
	Code		Code		Code		Code		Code		Code		Code
AseptiSure KS (0.2 µm Upstream)	CPK1	5"	53	0.1 µm	36	7P	A0	Silicone	SS	Non Sterile	1	1	01
AseptiSure KS (0.45 µm Upstream)	CPKX	10"	54	0.2 µm	01	7P without fin	A1	Viton	SV				
AseptiSure KS (0.65 µm Upstream)	CPK3	20"	55	0.45 µm	02	28 with fin	C0	EPDM	SE				
AseptiSure KS (0.8 µm Upstream)	CPK5	30"	56			'O'	D0	FEP Encapsulated Viton	FV				

Example:

CPKX	53	01	E0	SS	1	01
------	----	----	----	----	---	----

For Adapters and Elastomers availability refer Page 91.

Filters for Sterilization of Liquids:

PVDF Membrane Filters

mdi produces a wide range of Gamma sterilizable and steam sterilizable Hydrophilic PVDF membrane capsule and cartridge filters to meet filtration requirements of biopharmaceutical and pharmaceutical processing.

These filters meet key process requirements such as high retention efficiency, very high protein recoveries, extremely low extractables, high throughputs, wide chemical compatibility etc.

Applications

Sterile Filtration of:

- Antibodies
- Protein Solutions
- Buffers
- Cell culture media
- Small volume parenterals
- Vaccine concentrates

Types Available

Gamma Sterilizable Capsule Filters	<i>AseptiCap WS -γ</i>
Autoclavable Capsule Filters	<i>AseptiCap WS</i>
High Temperature Resistant Steam Sterilizable Cartridge Filters	<i>AseptiSure WS</i>

Quality Assurance

These filter devices are manufactured in Class 10,000 clean rooms under ISO 9001 : 2015 certified quality management systems and are validated to meet compendia and regulatory requirements.

Assurance	
Toxicity	Passes Biological Reactivity Test, In Vivo, as per USP <88> for Class VI plastics
Cytotoxicity	Passes Biological Reactivity Test, In Vitro, as per USP <87> for Cytotoxicity
Bioburden	Bioburden level is < 1000 cfu/filter device as per ISO 11737-1: 2018
Bacterial Endotoxin	Aqueous extracts exhibit < 0.25 EU/ml as established by Limulus Amebocyte Lysate (LAL) Test as per USP <85>
Non Fiber Releasing	Passes test as per USP and comply with USFDA 21 CFR Part 210.3 (b)(6) for fiber release
Extractables with WFI	Passes test as per USP <661>
Oxidizable Substances	Passes test as per USP <1231>
Particle Shedding	Passes test as per USP <788> for particulate matter in injections
TOC/Conductivity at 25 °C	Meets the WFI requirements of USP <643> for Total Organic Carbon and USP <645> for Water Conductivity after a specified volume of purified water flush
Indirect Food Additive	All Polypropylene components meet the FDA Indirect Food Additive requirements cited in 21 CFR 177.1520
Good Manufacturing Practice	These products are manufactured in a facility which adheres to Good Manufacturing Practices

Filter Selection Chart

Application Area	Key Application Requirements	Gamma Sterilizable	Steam Sterilizable	
		Capsule Filters	Capsule Filters	Cartridge Filters
Biopharmaceuticals				
Media preparation	- Mycoplasma removal (in case of mammalian Cell Culture)	<i>AseptiCap WS -γ</i> 0.1 μm PVDF Membrane Capsule Filter	<i>AseptiCap WS</i> 0.1 μm PVDF Membrane Capsule Filter	<i>AseptiSure WS</i> 0.1 μm High Temperature PVDF Membrane Cartridge Filter
	- Microbial retention (in case of microbial Fermentation)	<i>AseptiCap WS -γ</i> 0.2 μm PVDF Membrane Capsule Filter	<i>AseptiCap WS</i> 0.2 μm PVDF Membrane Capsule Filter	<i>AseptiSure WS</i> 0.2 μm High Temperature PVDF Membrane Cartridge Filter
Sterile filtration of vaccines and therapeutic proteins	- Absolute retention - Low protein binding - Low holdup volume	<i>AseptiCap WS -γ</i> 0.2 μm PVDF Membrane Capsule Filter	<i>AseptiCap WS</i> 0.2 μm PVDF Membrane Capsule Filter	-

AseptiCap WS Inline Capsule Filters (25 mm and 50 mm)

These are low protein binding hydrophilic PVDF gamma sterilizable membrane inline capsule filters, designed for sterile filtration of very small fluid volumes in formulation and process development labs.

Radiation Sterilizable: AseptiCap WS- γ

Autoclavable: AseptiCap WS



**Microbiologically Validated
as per ASTM F 838**

**Complies with
USFDA 21 CFR 210.3 (b) (6)**

**Meets and Exceeds
USFDA 21 CFR 177.1520**

Specifications

Construction			
Membrane	Hydrophilic PVDF		
Plastic Parts	Polypropylene		
Final Filter Pore Size	0.1 μm	0.2 μm	0.45 μm
Prefilter Pore Size	0.2 μm and 0.45 μm	0.8 μm and 0.45 μm	0.8 μm
Integrity Testing/Retention			
Bubble Point	≥ 28 psi (1.96 Kg/cm ²) with 50% IPA ≥ 70 psi (4.92 Kg/cm ²) with Water	≥ 50 psi (3.52 Kg/cm ²) with Water	≥ 30 psi (2.11 Kg/cm ²) with Water
Microbial Retention	LRV >7 for <i>Acholeplasma laidlawii</i> (ATCC 23206) per cm ²	LRV >7 for <i>Brevundimonas diminuta</i> (ATCC 19146) per cm ²	LRV >7 for <i>Serratia marcescens</i> (ATCC 14756) per cm ²
Size			
Size	25 mm	50 mm	
Effective Filtration Area (Nominal)	5 cm ²	20 cm ²	
Operational			
Max. Operating Temperature	55 °C	60 °C	
Max. Differential Pressure	75 Psi (5 Kg/cm ²) @ 25 °C	42 Psi (3 Kg/cm ²) @ 30 °C	
Hold-up Volume(with air purge)	<50 μL	<200 μL	
Burst Pressure	> 14 Kg/cm ²	> 8 Kg/cm ²	
Sterilization	By Irradiation	AseptiCap WS-γ : Gamma Irradiatable up to 50 kGy. These filters should not be autoclaved or in-line steam sterilized.	
	By Gas	AseptiCap WS : Sterilizable by Ethylene Oxide	
	By Autoclave	AseptiCap WS : Autoclavable at 125°C for 30 minutes, 2 cycles These cannot be In-line steam sterilized	
Shelf Life	2 years after gamma sterilization 3 years after Ethylene Oxide sterilization		

Ordering Information

AseptiCap WS and AseptiCap WS-γ, 25 mm

Type		Size		Pore Size		Inlet/Outlet		Radiation Sterilizable		X	Sterility		Pack Size	
	Code	Dia	Code		Code		Code		Code			Code		Code
AseptiCap WS (0.2 μm Upstream)	IWS1	25 mm	06	0.1 μm	36	Female Luer Lock	M	Yes	R		Non Sterile	1	100	04
AseptiCap WS (0.45 μm Upstream)	IWSX			0.2 μm	01	Male Luer Slip	N	No	X		EO Sterile	2		
AseptiCap WS (0.8 μm Upstream)	IWS5			0.45 μm	02	1/8" Hose Barb	H				Gamma Sterile	3		
						1/4" Hose Barb	B							

Example

IWSX	06	01	MN	R	X	1	04
------	----	----	----	---	---	---	----

AseptiCap WS and AseptiCap WS-γ, 50 mm

Type		Size		Pore Size		Inlet/Outlet		Radiation Sterilizable		X	Sterility		Pack Size	
	Code	Dia	Code		Code		Code		Code			Code		Code
AseptiCap WS (0.2 μm Upstream)	VWS1	50 mm	10	0.1 μm	36	1/4" SHB	B	Yes	R		Non Sterile	1	10	02
AseptiCap WS (0.45 μm Upstream)	VWSX			0.2 μm	01	3/4" Sanitary Flange	S	No	X		EO Sterile	2		
AseptiCap WS (0.8 μm Upstream)	VWS5			0.45 μm	02						Gamma Sterile	3		

Example

VWSX	10	36	BS	X	X	1	02
------	----	----	----	---	---	---	----

Note: Gamma Sterile filters cannot be Gamma Irradiated again

Example for Non Sterile: VWSX1036BBRX102 **Example for Gamma Sterile:** VWSX1036BBXX302

For End Connection availability and dimensions with different sizes refer Pages 89-90.

AseptiCap WS Small Capsule Filters (1", 2", 5", 8")

Microbiologically Validated as per ASTM F 838
Complies with USFDA 21 CFR 210.3 (b) (6)
Meets and Exceeds USFDA 21 CFR 177.1520

mdi AseptiCap WS are low protein binding hydrophilic PVDF membrane capsule filters offering serial filtration incorporating a larger pore size upstream membrane to protect the downstream membrane for enhanced throughput.

These capsule filters are validated to meet compendia and regulatory requirements and are well characterized. They meet key process requirements such as absolute retention efficiency, extremely low extractables, high throughputs, wide chemical compatibility and other important characteristics.

Radiation Sterilizable : AseptiCap WS- γ

Autoclavable : AseptiCap WS




Specifications


Construction				
Membrane	Hydrophilic PVDF			
Support (Drainage) Layers	Polyester			
Plastic Parts	Polypropylene			
Final Filter Pore Size	0.1 μm	0.2 μm	0.45 μm	
Prefilter Pore Size	0.2 μm , 0.45 μm	0.8 μm and 0.45 μm	0.8 μm	
Integrity Testing/Retention				
Bubble Point	≥ 28 psi (1.96 Kg/cm ²) with 50% IPA ≥ 70 psi (4.92 Kg/cm ²) with Water	≥ 50 psi (3.52 Kg/cm ²) with Water	≥ 30 psi (2.11 Kg/cm ²) with Water	
Microbial Retention	LRV >7 for <i>Acholeplasma laidlawii</i> (ATCC 23206) per cm ²	LRV >7 for <i>Brevundimonas diminuta</i> (ATCC 19146) per cm ²	LRV >7 for <i>Serratia marcescens</i> (ATCC 14756) per cm ²	
Size				
Size	1"	2"	5"	8"
Effective Filtration Area (Nominal)	250 cm ²	500 cm ²	1000 cm ²	2000 cm ²
Vent and Drain	$\frac{1}{4}$ " Hose Barb with Platinum Cured Silicone 'O' ring			
Operational				
Max. Operating Temperature	80 °C @ ≤ 30 psi (2 Kg/cm ²)			
Max. Differential Pressure	60 psi (4 Kg/cm ²) @ 30 °C			
Sterilization	By Irradiation	AseptiCap WS - γ : Gamma Irradiatable up to 50 kGy. These filters should not be autoclaved or in-line steam sterilized.		
	By Gas	AseptiCap WS : Sterilizable by Ethylene Oxide		
	By Autoclave	AseptiCap WS : Autoclavable at 125°C for 30 minutes, 2 cycles These cannot be In-line steam sterilized		
Shelf Life	2 years after gamma sterilization 3 years after Ethylene Oxide sterilization			

Water Flow Rates


0.2 µm AseptiCap WS, 1" Capsule Filters




0.2 µm AseptiCap WS, 2" Capsule Filters



0.2 µm AseptiCap WS, 5" Capsule Filters



0.2 µm AseptiCap WS, 8" Capsule Filters



End Connection Type:

A: ¼" Stepped Hose Barb **D:** ½" Hose Barb **E:** 1½" Sanitary Flange **J:** Quick Connector **U:** Female Luer Lock **Q:** Single Step ½" Hose Barb
S: ¾" Sanitary Flange

Ordering Information

AseptiCap WS and AseptiCap WS-γ

Type		Size		Pore Size		Inlet/Outlet		Radiation Sterilizable		Bell		Sterility		Pack Size	
	Code		Code		Code		Code		Code		Code		Code		Code
AseptiCap WS (0.2 µm Upstream)	DWS1	1"	51	0.1 µm	36	¼" SHB	A	Yes	R	Yes	B	Non Sterile	1	1	01
AseptiCap WS (0.45 µm Upstream)	DWSX	2"	52	0.2 µm	01	¼" MNPT	B	No*	X	No Bell	X	EO Sterile	2		
AseptiCap WS (0.8 µm Upstream)	DWS5	5"	53	0.45 µm	02	½" MNPT	C					Gamma Sterile	3		
		8"	57			½" Hose Barb	D								
						1½" Sanitary Flange	E								
						¾" Sanitary Flange	S								
						Quick Connector	J								
						Single Step ½" Hose Barb	Q								
						Female Luer Lock	U								
						Male Luer Slip	W								
						⅜" Hose Barb	N								
						⅝" Hose Barb	I								

Example:

DWSX	57	36	DD	R	X	1	01
------	----	----	----	---	---	---	----

*Gamma Sterile filters cannot be Gamma Irradiated again

Example for Non Sterile: DWSX5136DDRX101 Example for Gamma Sterile: DWSX5136DDXX301

For End Connection, bell availability and dimensions with different sizes refer Pages 89-90.

AseptiCap WS Large Capsule Filters (5", 10", 20", 30")

mdi AseptiCap WS are low protein binding hydrophilic PVDF membrane capsule filters offering serial filtration incorporating a larger pore size upstream membrane to protect the downstream membrane for enhanced throughput.

These capsule filters are validated to meet compendia and regulatory requirements and are well characterized. They meet key process requirements such as absolute retention efficiency, extremely low extractables, high throughputs, wide chemical compatibility and other important characteristics.

Radiation Sterilizable : AseptiCap WS- γ

Autoclavable : AseptiCap WS




Microbiologically Validated as per ASTM F 838
Complies with USFDA 21 CFR 210.3 (b) (6)
Meets and Exceeds USFDA 21 CFR 177.1520

Specifications


Construction				
Membrane	Hydrophilic PVDF			
Support (Drainage) Layers	Polyester			
Plastic Parts	Polypropylene			
Final Filter Pore Size	0.1 μm	0.2 μm	0.45 μm	
Prefilter Pore Size	0.2 μm , 0.45 μm	0.8 μm and 0.45 μm	0.8 μm	
Integrity Testing/Retention				
Bubble Point	≥ 28 psi (1.96 Kg/cm ²) with 50% IPA ≥ 70 psi (4.92 Kg/cm ²) with Water	≥ 50 psi (3.52 Kg/cm ²) with Water	≥ 30 psi (2.11 Kg/cm ²) with Water	
Max. Air Diffusion Flow for 10" Capsule Filters	≤ 30 ml/min @ 50 psi (3.52 Kg/cm ²) with Water	≤ 30 ml/min @ 37 psi (2.6 Kg/cm ²) with Water	≤ 30 ml/min @ 22 psi (1.54 Kg/cm ²) with Water	
Microbial Retention	LRV >7 for <i>Acholeplasma laidlawii</i> (ATCC 23206) per cm ²	LRV >7 for <i>Brevundimonas diminuta</i> (ATCC 19146) per cm ²	LRV >7 for <i>Serratia marcescens</i> (ATCC 14756) per cm ²	
Size				
Size	5"	10"	20"	30"
Effective Filtration Area (Nominal)	3000 cm ²	6000 cm ²	12000 cm ²	18000 cm ²
Vent and Drain	$\frac{1}{4}$ " Hose Barb with platinum cured Silicone 'O' ring			
Operational				
Max. Operating Temperature	80 °C @ ≤ 30 psi (2 Kg/cm ²)			
Max. Differential Pressure	60 psi (4 Kg/cm ²) @ 30 °C			
Sterilization	By Irradiation	AseptiCap WS-γ : Gamma Irradiatable up to 50 kGy. These filters should not be autoclaved or in-line steam sterilized.		
	By Gas	AseptiCap WS : Sterilizable by Ethylene Oxide		
	By Autoclave	AseptiCap WS : Autoclavable at 125°C for 30 minutes, 2 cycles These cannot be In-line steam sterilized		
Shelf Life	2 years after gamma sterilization 3 years after Ethylene Oxide sterilization			

Water Flow Rates


0.2 µm AseptiCap WS, 5" Large Capsule Filters




0.2 µm AseptiCap WS, 10" Large Capsule Filters



0.2 µm AseptiCap WS, 20" Large Capsule Filters



0.2 µm AseptiCap WS, 30" Large Capsule Filters



End Connection Type: E: 1½" Sanitary Flange

Ordering Information

AseptiCap WS and AseptiCap WS-γ

Type		Size		Pore Size		Inlet/Outlet		Radiation Sterilizable		Inline / T-line		Sterility		Pack Size	
	Code		Code		Code		Code		Code		Code		Code		Code
AseptiCap WS (0.2 µm Upstream)	LWS1	5"	53	0.1 µm	36	1½" Sanitary Flange	E	Yes	R	Inline	X	Non Sterile	1	1	01
AseptiCap WS (0.45 µm Upstream)	LWSX	10"	54	0.2 µm	01	Single Step ½" Hose Barb	Q	No	X	T-line**	T	EO Sterile	2		
AseptiCap WS (0.8 µm Upstream)	LWS5	20"	55	0.45 µm	02	¾" Sanitary Flange	S					Gamma Sterile	3		
		30"	56			3/8" Hose Barb	I								
						1" Hose Barb	Z								

Example:

LWSX	54	01	EE	R	T	1	01
------	----	----	----	---	---	---	----

* Size 5" is available in In-line Capsule Filters Only

** T-line Capsule Filters are available with 1½" Sanitary Flange connections Only

Note: Gamma Sterile filters cannot be Gamma Irradiated again

Example for Non Sterile: LWSX5402EERX101 Example for Gamma Sterile: LWSX5402EEXX301

For End Connection availability and dimensions with different sizes refer Pages 89-90.

AseptiSure WS Mini Cartridge filters

**Microbiologically Validated
as per ASTM F 838**
**Complies with
USFDA 21 CFR 210.3 (b) (6)**
**Meets and Exceeds
USFDA 21 CFR 177.1520**

AseptiSure WS are low protein binding hydrophilic PVDF membrane mini cartridge filters offering serial filtration incorporating a large pore size upstream membrane to protect the downstream membrane for enhanced throughput.

These cartridge filters are validated to meet compendia and regulatory requirements and are well characterized. They meet key process requirements such as absolute retention efficiency, extremely low extractables, high throughputs, wide chemical compatibility and other important characteristics.



Specifications

Construction			
Membrane	Hydrophilic PVDF		
Support (Drainage) Layers	Polyester		
Plastic Parts	Polypropylene		
Final Filter Pore Size	0.1 µm	0.2 µm	0.45 µm
Prefilter Pore Size	0.2 µm, 0.45 µm	0.8 µm and 0.45 µm	0.8 µm
Integrity Testing/Retention			
Bubble Point	≥ 28 psi (1.96 Kg/cm ²) with 50% IPA ≥ 70 psi (4.92 Kg/cm ²) with Water	> 50psi (3.52Kg/cm ²) with Water	> 30psi (2.11Kg/cm ²) with Water
Microbial Retention	LRV >7 for <i>Acholeplasma laidlawii</i> (ATCC 23206) per cm ²	LRV >7 for <i>Brevundimonas diminuta</i> (ATCC 19146) per cm ²	LRV >7 for <i>Serratia marcescens</i> (ATCC 14756) per cm ²
Size			
Size	2.5"		5"
Effective Filtration Area (Nominal)	1000 cm ²		2000 cm ²
Operational			
Typical Water Flow Rate for 0.2µm @ 0.70kg/cm ² @27°C	7.0 lpm		15.7 lpm
Max. Operating Temperature	80 °C @ ≤ 2 Kg/cm ² (30 psi)		
Max. Differential Pressure	50 psi (3.5 Kg/cm ²) @ 25 °C		
Reverse Pressure	≤ 10 psi (0.7 Kg/cm ²) @ 25 °C		
Sterilization	In-line steam sterilizable upto 135 °C for 30 minutes at a maximum differential pressure of 5 psi (0.35 Kg/cm ²), 3 cycles		

Ordering Information

Type		Size		Pore Size		Adapter		Elastomer		Sterility		Pack Size	
	Code		Code		Code		Code		Code		Code		Code
AseptiSure WS (0.2 µm Upstream)	CWH1	2.5"	50	0.1 µm	36	4463	E0	Silicone	SS	Non Sterile	1	1	01
AseptiSure WS (0.45 µm Upstream)	CWHX	5"	53	0.2 µm	01	4463B	H0						
AseptiSure WS (0.8 µm Upstream)	CWH5			0.45 µm	02	4440	U0						
						Seal-K	G0*			*G0 adapter is not available with elastomer. Please mention XX in place of elastomer code while ordering			
						Seal-O	F0						
						Seal-M	J0						

Example:

CWHX	50	01	E0	SS	1	01
------	----	----	----	----	---	----

For Adapters and Elastomers availability refer Page 91.

AseptiSure WS Standard Cartridge Filters

**Microbially Validated
as per ASTM F 838**

**Complies with
USFDA 21 CFR 210.3 (b) (6)**

**Meets and Exceeds
USFDA 21 CFR 177.1520**

AseptiSure WS are low protein binding hydrophilic PVDF membrane cartridge filters offering serial filtration incorporating a large pore size upstream membrane to protect the downstream membrane for enhanced throughput.



Specifications

Construction				
Membrane	Hydrophilic PVDF			
Support (Drainage) Layers	Polyester			
Plastic Parts	Polypropylene			
Final Filter Pore Size	0.1 µm	0.2 µm	0.45 µm	
Prefilter Pore Size	0.2 µm, 0.45 µm	0.8 µm and 0.45 µm	0.8 µm	
Integrity Testing/Retention				
Bubble Point	≥ 28 psi (1.96 Kg/cm ²) with 50% IPA ≥ 70 psi (4.92 Kg/cm ²) with Water	≥ 50 psi (3.52 Kg/cm ²) with Water	≥ 30 psi (2.11 Kg/cm ²) with Water	
Air Diffusion Flow (10" Cartridge Filter)	≤ 30 ml/min @ 50 psi (3.52 Kg/cm ²) with Water	≤ 30 ml/min @ 37 psi (2.6 Kg/cm ²) with Water	≤ 30 ml/min @ 22 psi (1.54 Kg/cm ²) with Water	
Microbial Retention	LRV >7 for <i>Acholeplasma laidlawii</i> (ATCC 23206) per cm ²	LRV >7 for <i>Brevundimonas diminuta</i> (ATCC 19146) per cm ²	LRV >7 for <i>Serratia marcescens</i> (ATCC 14756) per cm ²	
Size				
Size	5"	10"	20"	30"
Effective Filtration Area (Nominal)	3000 cm ²	6000 cm ²	12000 cm ²	18000 cm ²
Operational				
Typical Water Flow Rates (for 0.2µm @ 0.70 Kg/cm ² @ 27 °C)	18 lpm	33 lpm	60 lpm	84 lpm
Max. Operating Temperature	80 °C @ ≤ 30psi (2 Kg/cm ²)			
Max. Differential Pressure	50 psi (3.5 Kg/cm ²) @ 25 °C			
Reverse Pressure	≤ 10 psi (0.7 Kg/cm ²) @ 25 °C			
Sterilization	In-line steam sterilizable upto 135 °C for 30 minutes at a maximum differential pressure of 5 psi (0.35 Kg/cm ²), 3 cycles			

Ordering Information

Type		Size		Pore Size		Adapter		Elastomer		Sterility		Pack Size	
	Code		Code		Code		Code		Code		Code		Code
AseptiSure WS (0.2 µm Upstream)	CWH1	5"	53	0.1 µm	36	7P	A0	Silicone	SS	Non Sterile	1	1	01
AseptiSure WS (0.45 µm Upstream)	CWHX	10"	54	0.2 µm	01	7P without fin	A1	Viton	SV				
AseptiSure WS (0.8 µm Upstream)	CWH5	20"	55	0.45 µm	02	28 with fin	C0	EPDM	SE				
		30"	56			'O'	D0	FEP Encapsulated Viton	FV				

Example:

CWHX	54	01	A0	SS	1	01
------	----	----	----	----	---	----

For Adapters and Elastomers availability refer Page 91.

Filters for Sterilization of Liquids: Nylon 66 Membrane Filters

mdi Nylon 66 membrane filters are sterilizing grade filters offering absolute retention and wide chemical compatibility.

mdi Nylon 66 filter devices are available as:

Filter Type	Single Layer	Multiple Layer
Capsule Filters	<i>AseptiCap NL</i>	<i>AseptiCap NS</i>
Cartridge Filters	–	<i>AseptiSure NS</i>
Membrane Disc Filters	–	NN

Applications

- Sterilization of compatible solvents and chemicals
- Sterilization of disinfectants in pharmaceutical process and lab areas
- Filtration of hospital disinfectants
- Filtration of rinse water for endoscopes and other hospital equipment and surfaces
- Sterilizing filtration in pharmaceutical for aqueous and non aqueous solutions

Quality Assurance

These filter devices are manufactured in Class 10,000 clean rooms under ISO 9001 : 2015 certified quality management systems and are validated to meet compendia and regulatory requirements.

Assurance	
Toxicity	Passes Bioreactivity test, In Vivo, as per USP <88> for Class VI plastics
Bioburden	Bioburden level is < 1000 cfu/filter device as per ISO 11737-1: 2018
Bacterial Endotoxin	Aqueous extracts exhibit < 0.25 EU/ml as established by Limulus Amebocyte Lysate (LAL) Test as per USP <85>
Non Fiber Releasing	Passes test as per USP and comply with USFDA 21 CFR Part 210.3 (b)(6) for fiber release
Extractables with WFI	Passes test as per USP <661>
Oxidizable Substances	Within limits as specified in USP <1231>
Particle Shedding	Passes test as per USP <788> for particulate matter in injections
TOC/Conductivity at 25 °C	Meets the WFI requirements of USP <643> for Total Organic Carbon and USP <645> for Water Conductivity after a specified volume of purified water flush
Indirect Food Additive	All Polypropylene components meet the FDA Indirect Food Additive requirements cited in 21 CFR 177.1520
Good Manufacturing Practice	These products are manufactured in a facility which adheres to Good Manufacturing Practices

Filter Selection Chart

Application Area	Key Application Requirements	Capsule Filters	Cartridge Filters	Disc Filters
Sterilization of compatible solvents and chemicals	<ul style="list-style-type: none"> - Absolute retention - Wide chemical compatibility 	AseptiCap NL/NS 0.2 µm Nylon 66 Membrane Capsule Filter	AseptiSure NS 0.2 µm Nylon 66 Membrane Cartridge Filter	0.2 µm Nylon 66 Membrane Disc Filters
Sterilization of disinfectants in pharmaceutical labs and process areas	<ul style="list-style-type: none"> - Absolute retention - Wide chemical compatibility 	AseptiCap NL/NS 0.2 µm Nylon 66 Membrane Capsule Filter	AseptiSure NS 0.2 µm Nylon 66 Membrane Cartridge Filter	0.2 µm Nylon 66 Membrane Disc Filters
Filtration of hospital disinfectants	<ul style="list-style-type: none"> - Absolute retention - Wide chemical compatibility 	AseptiCap NL/NS 0.2 µm Nylon 66 Membrane Capsule Filter	AseptiSure NS 0.2 µm Nylon 66 Membrane Cartridge Filter	0.2 µm Nylon 66 Membrane Disc Filters
Filtration of rinse water for endoscopes and other hospital equipments	<ul style="list-style-type: none"> - Absolute retention 	AseptiCap NL/NS 0.2 µm Nylon 66 Membrane Capsule Filter	AseptiSure NS 0.2 µm Nylon 66 Membrane Cartridge Filter	-
Sterilizing filtration of aqueous and non aqueous solutions	<ul style="list-style-type: none"> - Absolute retention - Wide chemical compatibility 	AseptiCap NL/NS 0.2 µm Nylon 66 Membrane Capsule Filter	AseptiSure NS 0.2 µm Nylon 66 Membrane Cartridge Filter	0.2 µm Nylon 66 Membrane Disc Filters

AseptiCap NL/NS Small Capsule Filters (1", 2", 5", 8")

**Microbiologically Validated
as per ASTM F 838**

**Complies with
USFDA 21 CFR 210 (b) (6)**

**Meets and Exceeds
USFDA 21 CFR 177.1520**

AseptiCap NL/NS Nylon 66 membranes capsule filters are sterilizing grade filters offering absolute retention, wide chemical compatibility, and very low hold up volumes.

Single Layered: AseptiCap NL
Multiple Layered: AseptiCap NS




Specifications


Construction					
Membrane		Nylon 66			
Support (Drainage) Layers		Polyester			
Plastic Parts		Polypropylene			
Final Filter Pore Size		0.2 µm		0.45 µm	
Prefilter Pore Size (In case of AseptiCap NS)		0.8 µm, 0.45 µm		0.8 µm	
Integrity Testing/Retention					
Bubble Point (with water)		≥ 50 psi (3.52 Kg/cm ²)		≥ 30 psi (2.1 Kg/cm ²)	
Microbial Retention (LRV >7 for)		<i>Brevundimonas diminuta</i> (ATCC 19146) per cm ²		<i>Serratia marcescens</i> (ATCC 14756) per cm ²	
Size					
Size		1"	2"	5"	8"
Effective Filtration Area (Nominal)	AseptiCap NL	250 cm ²	900 cm ²	1800 cm ²	2700 cm ²
	AseptiCap NS	200 cm ²	700 cm ²	1400 cm ²	2100 cm ²
Vent and Drain		¼" Hose Barb with Silicone 'O' ring			
Operational					
Max. Operating Temperature		80 °C @ ≤ 2 Kg/cm ² (30 psi)			
Max. Differential Pressure		4 Kg/cm ² (60 psi) @ 30 °C			
Sterilization	By Gas	Sterilizable by Ethylene Oxide			
	By Autoclave	Autoclavable at 125 °C for 30 minutes. Cannot be In-line steam sterilized			

Water Flow Rates


0.2µm AseptiCap NS, 1" Capsule Filter




0.2µm AseptiCap NS, 2" Capsule Filter



0.2µm AseptiCap NS, 5" Capsule Filter



0.2µm AseptiCap NS, 8" Capsule Filter



End Connection Type:

A: ¼" Stepped Hose Barb **Q:** Single Step ½" Hose Barb **E:** 1½" Sanitary Flange **S:** ¾" Sanitary Flange **D:** ½" Hose Barb

Ordering Information

Type		Size		Pore Size		Inlet/Outlet		X	Bell		Sterility		Pack Size	
	Code		Code		Code		Code			Code		Code		Code
AseptiCap NL	DNLX	1"	51	0.2 µm	01	¼" SHB	A		Yes	B	Non Sterile	1	1	01
AseptiCap NS (0.45 µm Upstream)	DNSX	2"	52	0.45 µm	02	¼" MNPT	B		No Bell	X	EO Sterile	2		
AseptiCap NS (0.8 µm Upstream)	DNS5	5"	53			½" MNPT	C							
		8"	57			½" Hose Barb	D							
						1½" Sanitary Flange	E							
						¾" Sanitary Flange	S							
						Quick Connector	J							
						Single Step ½" Hose Barb	Q							
						Female Luer Lock	U							
						Male Luer Slip	W							
						⅜" Hose Barb	N							
						⅝" Hose Barb	I							

Example:

DNSX	52	01	DD	X	X	1	01
------	----	----	----	---	---	---	----

For End Connection, bell availability and dimensions with different sizes refer Pages 89-90.

AseptiCap NS Large Capsule Filters (5", 10", 20", 30")

AseptiCap NS Nylon 66 large membranes capsule filters are double layered sterilizing grade filters offering absolute retention, wide chemical compatibility, and serial filtration for enhanced throughputs.

**Microbially Validated
as per ASTM F 838**

**Complies with
USFDA 21 CFR 210.3 (b) (6)**

**Meets and Exceeds
USFDA 21 CFR 177.1520**





Specifications



Construction				
Membrane	Nylon 66			
Support (Drainage) Layers	Polyester			
Plastic Parts	Polypropylene			
Final Filter Pore Size	0.2 µm		0.45 µm	
Prefilter Pore Size	0.8 µm, 0.45 µm		0.8 µm	
Integrity Testing/Retention				
Bubble Point (with water)	≥ 50 psi (3.52 Kg/cm ²)		≥ 30 psi (2.1 Kg/cm ²)	
Air Diffusion Flow for 10" Capsule Filters (with water)	≤ 30 ml/min @ 37 psi (2.60 kg/cm ²)		≤ 30 ml/min @ 22 psi (1.54 kg/cm ²)	
Microbial Retention (LRV >7 for)	<i>Brevundimonas diminuta</i> (ATCC 19146) per cm ²		<i>Serratia marcescens</i> (ATCC 14756) per cm ²	
Size				
Size	5"	10"	20"	30"
Effective Filtration Area (Nominal)	3000 cm ²	6000 cm ²	12000 cm ²	18000 cm ²
Vent and Drain	¼" Hose Barb with Silicone 'O' ring			
Operational				
Max. Operating Temperature	80 °C @ ≤ 30 psi (2 Kg/cm ²)			
Max. Differential Pressure	60 psi (4 Kg/cm ²) @ 30 °C			
Sterilization	By Gas	Sterilizable by Ethylene Oxide		
	By Autoclave	Autoclavable at 125 °C for 30 minutes. Cannot be In-line steam sterilized.		

Water Flow Rates

0.2 µm AseptiCap NS, 5" Large Capsule Filters



0.2 µm AseptiCap NS, 10" Large Capsule Filters



End Connection Type: Q: Single Step 1/2" Hose Barb E: 1 1/2" Sanitary Flange

Ordering Information

Type		Size		Pore Size		Inlet/Outlet		X	Inline / T-line		Sterility		Pack Size	
	Code		Code		Code		Code			Code		Code		Code
AseptiCap NS (0.45 µm Upstream)	LNSX	5"	53	0.2 µm	01	1 1/2" Sanitary Flange	E		Inline	X	Non Sterile	1	1	01
		10"	54	0.45 µm	02	Single Step 1/2" Hose Barb	Q		T-line	T	EO Sterile	2		
AseptiCap NS (0.8 µm Upstream)	LNSS	20"	55			3/4" Sanitary Flange	S							
		30"	56			3/8" Hose Barb	I							
						1" Hose Barb	Z							

Example:

LNSX	54	01	QQ	X	T	1	01
------	----	----	----	---	---	---	----

* Size 5" is available in In-line Capsule Filters Only

For End Connection availability and dimensions with different sizes refer Pages 89-90.

AseptiSure NS Mini Cartridge filters

**Microbiologically Validated
as per ASTM F 838**

**Complies with
USFDA 21 CFR 210.3 (b) (6)**

**Meets and Exceeds
USFDA 21 CFR 177.1520**

mdi AseptiSure NS Nylon 66 membrane mini cartridge filters are sterilizing grade filters offering absolute retention and wide chemical compatibility. These filters offer serial filtration for enhanced throughput. The upstream layer is of larger pore size to protect the downstream final filtration layer.


Specifications

Construction		
Membrane	Nylon 66	
Support (Drainage) Layers	Polyester	
Plastic Parts	Polypropylene	
Final Filter Pore Size	0.2 µm	0.45 µm
Prefilter Pore Size	0.8 µm, 0.45 µm	0.8 µm
Integrity Testing/Retention		
Bubble Point (with water)	≥ 50 psi (3.52 Kg/cm ²)	≥ 30 psi (2.1 Kg/cm ²)
Microbial Retention (LRV >7 for)	<i>Brevundimonas diminuta</i> (ATCC 19146) per cm ²	<i>Serratia marcescens</i> (ATCC 14756) per cm ²
Size		
Size	2.5"	5"
Effective Filtration Area (Nominal)	1000 cm ²	2000 cm ²
Operational		
Max. Operating Temperature	80 °C @ ≤ 2 Kg/cm ² (30 psi)	
Max. Differential Pressure	< 3.5Kg/cm ² (50 psi) @ 25°C	
Reverse Pressure	< 0.7 Kg/cm ² (10 psi) @ 25 °C	
Sterilization	Autoclavable/In-line Steam Sterilizable at 121 °C for 30 minutes @ a maximum differential pressure of 3psi (0.21 kg/cm ²)	



Water Flow Rates

0.2µm AseptiSure NS Mini Cartridge Filters



Ordering Information

Type		Size		Pore Size		Adapter		Elastomer		Sterility		Pack Size	
	Code		Code		Code		Code		Code		Code		Code
AseptiSure NS (0.45 µm Upstream)	CPNX	2.5"	50	0.2 µm	01	4463	E0	Silicone	SS	Non Sterile	1	1	01
AseptiSure NS (0.8 µm Upstream)	CPN5	5"	53	0.45 µm	02	4463B	H0						
						4440	U0						
						Seal-K	G0*						
						Seal-O	F0						
						Seal-M	J0						

*G0 adapter is not available with elastomer. Please mention XX in place of elastomer code while ordering

Example:

CPNX	50	01	E0	SS	1	01
------	----	----	----	----	---	----

For Adapters and Elastomers availability refer Page 91.

AseptiSure NS Standard Cartridge Filters

AseptiSure NS Nylon 66 membrane cartridge filters are sterilizing grade filters offering absolute retention and wide chemical compatibility. These filters offer serial filtration for enhanced throughput. The upstream layer is of larger pore size to protect the downstream final filtration layer.

**Microbiologically Validated
as per ASTM F 838**

**Complies with
USFDA 21 CFR 210.3 (b) (6)**

**Meets and Exceeds
USFDA 21 CFR 177.1520**


Specifications

Construction				
Membrane	Nylon 66			
Support (Drainage) Layers	Polyester			
Plastic Parts	Polypropylene			
Final Filter Pore Size	0.2 µm	0.45 µm		
Prefilter Pore Size	0.8 µm, 0.45 µm	0.8 µm		
Integrity Testing/Retention				
Air Diffusion Flow (with water wetted) (10" Cartridge Filter)	≤ 30 ml/min @ 37 psi (2.60 kg/cm ²)	≤ 30 ml/min @ 22 psi (1.54 kg/cm ²)		
Microbial Retention (LRV >7 for)	<i>Brevundimonas diminuta</i> (ATCC 19146) per cm ²	<i>Serratia marcescens</i> (ATCC 14756) per cm ²		
Size				
Size	5"	10"	20"	30"
Effective Filtration Area (Nominal)	3000 cm ²	6000 cm ²	12000 cm ²	18000 cm ²
Operational				
Max. Operating Temperature	80 °C @ ≤ 30 psi (2 Kg/cm ²)			
Max. Differential Pressure	50 psi (3.5Kg/cm ²) @ 25 °C			
Reverse Pressure	< 10 psi (0.7 Kg/cm ²) @ 25 °C			
Sterilization	Autoclavable/In-line Steam Sterilizable at 121 °C for 30 minutes @ a maximum differential pressure of 3psi (0.21 kg/cm ²)			



Water Flow Rates

AseptiSure NS, 10" Cartridge Filters



Ordering Information

Type		Size		Pore Size		Adapter		Elastomer		Sterility		Pack Size	
	Code		Code		Code		Code		Code		Code		Code
AseptiSure NS (0.45 µm Upstream)	CPNX	5"	53	0.2 µm	01	7P	A0	Silicone	SS	Non Sterile	1	1	01
AseptiSure NS (0.8 µm Upstream)	CPN5	10"	54	0.45 µm	02	7P without fin	A1	Viton	SV				
		20"	55			28 with Fin	C0	EPDM	SE				
		30"	56			'O'	D0	FEP Encapsulated Viton	FV				

Example:

CPNX	54	01	A0	SS	1	01
------	----	----	----	----	---	----

For Adapters and Elastomers availability refer Page 91.

Nylon 66 Membrane Disc Filters - Type NN

**Microbiologically Validated
as per ASTM F 838**

**Complies with
USFDA 21 CFR 210.3 (b) (6)**

NN membrane disc filters are double layered, hydrophilic, non-media migrating, biologically inert, plain white absolute membrane filters offering wide chemical compatibility.



Specifications

Construction				
Membrane	Nylon 66			
Pore Size	0.2 µm, 0.45 µm, 0.8 µm, 1.2 µm			
Size	90 mm, 142 mm, 293 mm			
Integrity Testing/Retention				
Bubble Point (with Water)	0.2 µm	≥ 50 psi (3.51 Kg/cm ²) with water		
	0.45 µm	≥ 32 psi (2.25 Kg/cm ²) with water		
Microbial Retention	0.2 µm	LRV > 7 for <i>Brevundimonas diminuta</i>		
	0.45 µm	LRV > 7 for <i>Serratia marcescens</i>		
Operational				
Max. Operating Temperature	80 °C continuous			
Max. Differential Pressure	5 Kg/cm ²			
Water Flow Rates (27 °C @ 0.70 Kg/cm ²)	0.2 µm	0.45 µm	0.8 µm	1.2 µm
	14 ml/min/cm ²	37 ml/min/cm ²	120 ml/min/cm ²	180 ml/min/cm ²
Sterilization	Autoclavable at 121 °C for 30 minutes			

Ordering Information

Type		Size		Pore Size		XX	XX	Sterility		Pack Size	
Code		Code		Code				Code		Code	
NN	NNXX	90 mm	14	0.2 µm	01			Non Sterile	1	50	03
		142 mm	16	0.45 µm	02						
		293 mm	19	0.8 µm	03						
				1.2 µm	10						
↓ ↓ ↓ ↓ ↓ ↓ ↓ ↓ ↓ ↓ ↓ ↓											
Example:	NNXX	14	01	XX	XX	1	03				

Filters for Sterilization of Air/Gases

PVDF Membrane Filter Devices

mdi Gamma Sterilizable *AseptiVent VF-γ* are hydrophobic PVDF membrane single use capsule filters with a wide range of end connections and different sizes for linear scalability to use with disposable single use assemblies for biopharmaceutical processes.

These filters are validated for microbial retention with liquid bacterial challenge test as per ASTM F838 to ensure reliable performance under worst case conditions.

PTFE Membrane Filter Devices

mdi produces a wide range of PTFE membrane capsule and cartridge filters to meet filtration requirements of biopharmaceutical and pharmaceutical processing.

These filters are validated for microbial retention with liquid bacterial challenge test as per ASTM F838 to provide a high degree of sterility assurance for critical applications such as bioreactor/fermentor venting etc. As they offer wide chemical compatibility with organic solvents, these are ideal for manufacture of sterile API.

Quality Assurance

These filter devices are manufactured in Class 10,000 clean rooms under ISO 9001: 2015 certified quality management systems and are validated to meet compendia and regulatory requirements.

Filter Type	Single Use	Multiple Use
Gamma Sterilizable Capsule Filters	<i>AseptiVent VF-γ</i>	–
Autoclavable Capsule Filters	–	<i>AseptiVent TF</i>
Steam Sterilizable Cartridge Filters	–	<i>AseptiSure TF</i>
High Temperature Resistant Steam Sterilizable Cartridge Filters	–	<i>AseptiSure TH</i>

Applications

- Sterile air sparging in fermentors and bioreactors
- Sterile venting of cell factories, bioreactors and fermentors
- Fermentor exhaust
- Sterilization of environmental air in isolators
- Venting of sterile collection vessels
- Cleaning sterile surfaces
- WFI tank venting
- Nitrogen blanketing
- Sterile filtration of API and solvents
- Dry powder injectable filling
- Sterile air for dryers and micronizers

Assurance	
Toxicity	Passes Bioreactivity test, In Vivo, as per USP <88> for Class VI plastics
Bioburden	Bioburden level is < 1000 cfu/filter device as per ISO 11737-1: 2018
Bacterial Endotoxin	Aqueous extracts exhibit < 0.5 EU/ml as established by Limulus Amebocyte Lysate (LAL) Test as per USP <85>
Non Fiber Releasing	Passes test as per USP and comply with USFDA 21 CFR Part 210.3 (b)(6) for fiber release
Oxidizable Substances	Within limits as specified in USP <1231>
Particle Shedding	Passes test as per USP <788> for particulate matter in injections
Indirect Food Additive	All Polypropylene components meet the FDA Indirect Food Additive requirements cited in 21 CFR 177.1520
Good Manufacturing Practice	These products are manufactured in a facility which adheres to Good Manufacturing Practices

Filter Selection Chart

Application Area	Key Application Requirements	Gamma Sterilizable Capsule Filters	Steam Sterilizable	
			Capsule Filters	Cartridge Filters
Sterile venting for fermentors and bioreactors	<ul style="list-style-type: none"> - Absolute retention - High hydrophobicity - High flow rates 	AseptiVent VF-γ 0.2 μm PVDF Membrane Capsule Filters	AseptiVent TF 0.2 μm PTFE Membrane Capsule Filter	AseptiSure TF/TH 0.2 μm PTFE Membrane Cartridge Filter
Sterile air sparging in fermentors and bioreactors	<ul style="list-style-type: none"> - Absolute retention - High hydrophobicity - High flow rates 	AseptiVent VF-γ 0.2 μm PVDF Membrane Capsule Filters	AseptiVent TF 0.2 μm PTFE Membrane Capsule Filter	AseptiSure TF/TH 0.2 μm PTFE Membrane Cartridge Filter
Sterile air for cell factories	<ul style="list-style-type: none"> - Absolute retention - High hydrophobicity 	AseptiVent VF-γ 0.2 μm PVDF Membrane Capsule Filters	AseptiVent TF 0.2 μm PTFE Membrane Capsule Filter	–
Venting of small bioreactors	<ul style="list-style-type: none"> - Absolute retention - High hydrophobicity 	AseptiVent VF-γ 0.2 μm PVDF Membrane Capsule Filters	AseptiVent TF 0.2 μm PTFE Membrane Capsule Filter	–
Fermentor exhaust	<ul style="list-style-type: none"> - Absolute retention - High hydrophobicity - High flow rates 	–	AseptiVent TF 0.2 μm PTFE Membrane Capsule Filter	AseptiSure TF/TH 0.2 μm PTFE Membrane Cartridge Filter
Venting of sterile collection vessels	<ul style="list-style-type: none"> - Absolute retention - High hydrophobicity - High flow rates 	AseptiVent VF-γ 0.2 μm PVDF Membrane Capsule Filters	AseptiVent TF 0.2 μm PTFE Membrane Capsule Filter	AseptiSure TF/TH 0.2 μm PTFE Membrane Cartridge Filter
Nitrogen blanketing in sterile API	<ul style="list-style-type: none"> - Absolute retention - High flow rates 	–	AseptiVent TF 0.2 μm PTFE Membrane Capsule Filter	AseptiSure TF/TH 0.2 μm PTFE Membrane Cartridge Filter
Cleaning sterile surfaces	<ul style="list-style-type: none"> - Absolute retention - High flow rates 	–	AseptiVent TF 0.2 μm PTFE Membrane Capsule Filter	AseptiSure TF/TH 0.2 μm PTFE Membrane Cartridge Filter
Dry powder injectable filling	<ul style="list-style-type: none"> - Absolute retention - High flow rates 	–	AseptiVent TF 0.2 μm PTFE Membrane Capsule Filter	AseptiSure TF/TH 0.2 μm PTFE Membrane Cartridge Filter
WFI tank venting	<ul style="list-style-type: none"> - Absolute retention - High hydrophobicity - High flow rates 	–	–	AseptiSure TF/TH 0.2 μm PTFE Membrane Cartridge Filter
Sterile filtration of API/Solvents	<ul style="list-style-type: none"> - Absolute retention - High flow rates 	–	–	AseptiSure TF/TH 0.2 μm PTFE Membrane Cartridge Filter
Sterile air for dryers and micronizers	<ul style="list-style-type: none"> - Absolute retention - High flow rates 	–	AseptiVent TF 0.2 μm PTFE Membrane Capsule Filter	AseptiSure TF/TH 0.2 μm PTFE Membrane Cartridge Filter

AseptiVent VF-γ Inline Capsule Filters (25mm, 37mm, 50mm)

0.2 μm **AseptiVent VF-γ**, Gamma Sterilizable PVDF membrane vent filters are validated for microbial retention with liquid bacterial challenge test as per ASTM F838 to provide a high degree of sterility assurance for critical applications such as small bioreactors, sterile tank venting, bottle venting, barrier filter for vacuum pump etc.

**Microbially Validated
as per ASTM F 838**

**Complies with
USFDA 21 CFR 210.3 (b) (6)**

**Meets and Exceeds
USFDA 21 CFR 177.1520**




Specifications


Construction			
Final Filter Pore Size	0.2 μm		
Membrane	Hydrophobic PVDF		
Plastic Parts	Polypropylene		
Integrity Testing/Retention			
Bubble Point	≥ 18 psi (1.27 Kg/cm ²) with 50% IPA/ Water Solution		
Microbial Retention	LRV >7 for <i>Brevundimonas diminuta</i> (ATCC 19146) per cm ²		
Size			
Size	25 mm	37 mm	50 mm
Effective Filtration Area (Nominal)	5 cm ²	10 cm ²	20 cm ²
Operational			
Max. Operating Temperature	60 °C		
Max. Differential Pressure	1.5 Kg/cm ² (21 psi) @ 30° C		
Burst Pressure	> 14 Kg/cm ²	> 8 Kg/cm ²	> 8 Kg/cm ²
Sterilization By Irradiation	Gamma Irradiatable up to 50 kGy. These filters should not be autoclaved or in-line steam sterilized.		
Shelf Life	2 years after Gamma Sterilization		

Air Flow Rates

0.2 µm AseptiVent VF-γ, 37 mm Capsule Filters



0.2 µm AseptiVent VF-γ, 50 mm Capsule Filters



End Connection Type: B: 1/4" Stepped Hose Barb S: 3/4" Sanitary Flange

Ordering Information

AseptiVent VF -γ, 25 mm

Type		Size		Pore Size		Inlet/Outlet		Radiation Sterilizable		X	Sterility		Pack Size	
	Code		Code		Code		Code		Code			Code		Code
AseptiVent VF	IVFX	25 mm	06	0.2 µm	01	Female Luer Lock	M	Yes	R		Non Sterile	1	100	04
						Male Luer Slip	N	No*	X		Gamma Sterile	3		
						1/8" Hose Barb	H							
						1/4" Hose Barb	B							
Example:	IVFX	06	01	MN				R		X	1		04	

*Gamma Sterile filters cannot be Gamma Irradiated again

AseptiVent VF -γ, 37 mm and 50 mm

Type		Size		Pore Size		Inlet/Outlet		Radiation Sterilizable		X	Sterility		Pack Size	
	Code		Code		Code		Code		Code			Code		Code
AseptiVent VF	IVFX	37 mm	08	0.2 µm	01	1/4" SHB	B	Yes	R		Non Sterile	1	10	02
		50 mm	10			3/4" Sanitary Flange*	S	No**	X		Gamma Sterile	3		
Example:	IVFX	08	01	BB				R		X	1		02	

*3/4" Sanitary Flange connection is available only in 50 mm filter

**Gamma Sterile filters cannot be Gamma Irradiated again

For End Connection availability and dimensions with different sizes refer Pages 89-90.

AseptiVent VF-γ Small Capsule Filters (1", 2", 5", 8")

AseptiVent VF-γ PVDF membrane vent capsule filters are validated for Microbial Retention with liquid bacterial challenge test as per ASTM F838 to provide a high degree of sterility assurance for critical applications.

Microbially Validated as per ASTM F 838
Complies with USFDA 21 CFR 210.3 (b) (6)
Meets and Exceeds USFDA 21 CFR 177.1520




Specifications


Construction				
Pore Size	0.2 μm			
Membrane	Hydrophobic PVDF			
Support (Drainage) Layers	Polyester			
Plastic Parts	Polypropylene			
Integrity Testing/Retention				
Bubble Point	≥ 18 psi (1.26 Kg/cm ²) with 50% IPA/Water Solution			
Microbial Retention	LRV >7 for <i>Brevundimonas diminuta</i> (ATCC 19146) per cm ²			
Size				
Size	1"	2"	5"	8"
Effective Filtration Area (Nominal)	250 cm ²	500 cm ²	1000 cm ²	2000 cm ²
Vent and Drain	¼" Hose Barb with Silicone 'O' ring			
Operational				
Max. Operating Temperature	80 °C @ ≤ 30 psi (2 Kg/cm ²)			
Max. Differential Pressure	60 psi (4 Kg/cm ²) @ 30 °C			
Sterilization By Irradiation	Gamma Irradiatable up to 50 kGy. These filters should not be autoclaved or in-line steam sterilized.			
Shelf Life	2 years after Gamma Sterilization			

Air Flow Rates


0.2 µm AseptiVent VF-γ, 1" Capsule Filters, EE Connection




0.2 µm AseptiVent VF-γ, 2" Capsule Filters, EE Connection



0.2 µm AseptiVent VF-γ, 5" Capsule Filters, EE Connection



0.2 µm AseptiVent VF-γ, 8" Capsule Filters, EE Connection



End Connection Type E: 1½" Sanitary Flange

Ordering Information

Type		Size		Pore Size		Inlet/Outlet		Radiation Sterilizable		X	Sterility		Pack Size	
	Code		Code		Code		Code		Code				Code	
AseptiVent VF	DVLX	1"	51	0.2 µm	01	¼" SHB	A	Yes	R		Non Sterile	1	1	01
		2"	52			½" Hose Barb	D	No*	X		Gamma Sterile	3		
		5"	53			1½" Sanitary Flange	E							
		8"	57			¾" Sanitary Flange	S							
						Quick Connector	J							
						Single Step ½" Hose Barb	Q							
						Female Luer Lock	U							
						Male Luer Slip	W							
						⅜" Hose Barb	N							
						⅝" Hose Barb	I							

Example:

DVLX	57	01	EE	R	X	1	01
------	----	----	----	---	---	---	----

* Gamma Sterile filters cannot be Gamma Irradiated again

For End Connection availability and dimensions with different sizes refer Pages 89-90.

AseptiVent VF-γ Large Capsule Filters (5", 10", 20", 30")

AseptiVent VF-γ PVDF membrane vent filters are validated for microbial retention with liquid bacterial challenge test as per ASTM F838 to provide a high degree of sterility assurance for critical applications.



**Microbially Validated
as per ASTM F 838**

**Complies with
USFDA 21 CFR 210.3 (b) (6)**


**Meets and Exceeds
USFDA 21 CFR 177.1520**

Specifications


Construction				
Pore Size	0.2 µm			
Membrane	Hydrophobic PVDF			
Support (Drainage) Layers	Polyester			
Integrity Testing/Retention				
Bubble Point	≥ 18 psi (1.26 Kg/cm ²) with 50% IPA/Water Solution			
Microbial Retention	LRV >7 for <i>Brevundimonas diminuta</i> (ATCC 19146) per cm ²			
Size				
Size	5"	10"	20"	30"
Effective Filtration Area (Nominal)	3000 cm ²	6000 cm ²	12000 cm ²	18000 cm ²
Vent and Drain	¼" Hose Barb with Silicone 'O' ring			
Operational				
Typical Air Flow Rate	67 Nm ³ /h @ ΔP = 2 psi (15 psi inlet)			
Max. Operating Temperature	80° C @ 2Kg/cm ² (30psi)			
Max. Differential Pressure	4Kg/cm ² (60psi) @ 30° C			
Minimum Acceptable Bubble Point (with 50% IPA)	≥ 1.26 Kg/cm ² (18 psi)			
Sterilization By Irradiation	Gamma Irradiatable up to 50 kGy. These filters should not be autoclaved or in-line steam sterilized.			

Air Flow Rates

0.2 µm AseptiVent VF-γ, 5" Capsule Filters, EE Connection



0.2 µm AseptiVent VF-γ, 10" Capsule Filters, EE Connection



End Connection Type: E: 1½" Sanitary Flange

Ordering Information

Type		Size		Pore Size		Inlet/Outlet		Radiation Sterilizable		Inline / T-line		Sterility		Pack Size	
	Code		Code		Code				Code		Code		Code		Code
AseptiVent VF	LVLX	5"	53	0.2 µm	01	1½" Sanitary Flange	E	Yes	R	Inline	X	Non Sterile	1	1	01
		10"	54			Single Step ½" Hose Barb	Q	No*	X	T-line	T	Gamma Sterile	3		
		20"	55			¾" Sanitary Flange	S								
		30"	56			¾" Hose Barb	I								
						1" Hose Barb	Z								

Example:

LVLX	54	01	EE	R	X	1	01
------	----	----	----	---	---	---	----

* Gamma Sterile filters cannot be Gamma Irradiated again

**Size 5" is available in In-line Capsule Filters Only

For End Connection availability and dimensions with different sizes refer Pages 89-90.

AseptiVent TF Inline Capsule Filters (25 mm, 37 mm, 50 mm)

AseptiVent TF Disposable Inline PTFE gas filters are convenient pre-fabricated devices used for sterilization of gases and as a bacterial air vent in various pharmaceutical and biopharmaceutical processes.

Microbially Validated as per ASTM F 838
Complies with USFDA 21 CFR 210.3 (b) (6)
Meets and Exceeds USFDA 21 CFR 177.1520




Specifications


Construction			
Membrane	Hydrophobic PTFE		
Plastic Parts	Polypropylene		
Final Filter Pore Size	0.2 µm	0.45 µm	
Integrity Testing/Retention			
Bubble Point	≥ 22 psi (1.54 Kg/cm ²) with 70% IPA/Water Solution	≥ 10 psi (0.7 Kg/cm ²) with 70% IPA/Water Solution	
Microbial Bacterial Retention	LRV >7 for <i>Brevundimonas diminuta</i> (ATCC 19146) per cm ²	LRV >7 for <i>Serratia marcescens</i> ATCC 14756) per cm ²	
Size			
Size	25 mm	37 mm	50 mm
Effective Filtration Area (Nominal)	5 cm ²	10 cm ²	20 cm ²
Operational			
Max. Operating Temperature	60 °C		
Max. Differential Pressure	42 psi (3 Kg/cm ²) @ 30 °C		
Burst Pressure	> 14 Kg/cm ²	> 8 Kg/cm ²	> 8 Kg/cm ²
Sterilization	By Gas	Sterilizable by Ethylene Oxide	
	By Autoclave	Autoclavable at 125 °C for 30 minutes, 30 cycles. Cannot be In-line steam sterilized	
Shelf Life	3 years after Ethylene Oxide Sterilization		

Air Flow Rates


0.2 µm AseptiVent TF, 25 mm Capsule Filters



0.2 µm AseptiVent TF, 37 mm Capsule Filters



0.2 µm AseptiVent TF, 50 mm Capsule Filters



End Connection Type:

B: ¼" Stepped Hose Barb C: 1/8" MNPT S: ¾" Sanitary Flange

Ordering Information

AseptiVent TF- 25 mm

Type		Size		Pore Size		Inlet/Outlet		XX	Sterility		Pack Size	
	Code		Code		Code		Code			Code		Code
AseptiVent TF	ITFX	25 mm	06	0.2 µm	01	Female Luer Lock	M		Non Sterile	1	100	04
				0.45 µm	02	Male Luer Slip	N		EO Sterile	2		
						½" Hose Barb	H					
						¼" Hose Barb	B					
Example:												
	ITFX		06		01		MN			1		04

AseptiVent TF- 37 mm, 50 mm

Type		Size		Pore Size		Inlet/Outlet		XX	Sterility		Pack Size	
	Code		Code		Code		Code			Code		Code
AseptiVent TF	ITFX	37 mm*	08	0.2 µm	01	¼" SHB	B		Non Sterile	1	10	02
		50 mm	10	0.45 µm	02	½" MNPT	C		EO Sterile	2		
						¾" Sanitary Flange	S					
Example:												
	ITFX		08		01		BB			1		02

* Note: AseptiVent TF- 37 mm is available with BB connection only

For End Connection availability and dimensions with different sizes refer Pages 89-90.

AseptiVent TF Small Capsule Filters (1", 2", 5", 8")

AseptiVent TF capsule filters employ hydrophobic PTFE membrane offering absolute retention and very wide chemical compatibility making these useful for sterile filtration of air/gases as well as aggressive solvents.

Microbially Validated as per ASTM F 838
Complies with USFDA 21 CFR 210.3 (b) (6)
Meets and Exceeds USFDA 21 CFR 177.1520




Specifications


Construction				
Membrane	Hydrophobic PTFE			
Support (Drainage) Layers	Polypropylene			
Plastic Parts	Polypropylene			
Final Filter Pore Size	0.2 µm		0.45 µm	
Integrity Testing/Retention				
Bubble Point (with 70% IPA Wetted)	≥ 22 psi (1.55 Kg/cm ²)		≥ 10 psi (0.7 Kg/cm ²)	
Microbial Retention (LRV >7 for)	<i>Brevundimonas diminuta</i> (ATCC 19146) per cm ²		<i>Serratia marcescens</i> (ATCC 14756) per cm ²	
Size				
Size	1"	2"	5"	8"
Effective Filtration Area (Nominal)	250 cm ²	500 cm ²	1000 cm ²	2000 cm ²
Vent and Drain	¼" Hose Barb with Silicone 'O' ring			
Operational				
Max. Operating Temperature	80 °C @ ≤30 psi (2 Kg/cm ²)			
Max. Differential Pressure	60 psi (4 Kg/cm ²) @ 30 °C			
Sterilization	By Gas	Sterilizable by Ethylene Oxide		
	By Autoclave	Autoclavable at 125 °C for 30 minutes, 50 cycles. Cannot be In-line steam sterilized		
Shelf Life	3 years after Ethylene Oxide sterilization			

Air Flow Rates


0.2µm AseptiVent TF, 1" Capsule Filter




0.2µm AseptiVent TF, 2" Capsule Filter



0.2µm AseptiVent TF, 5" Capsule Filter



0.2µm AseptiVent TF, 8" Capsule Filter



End Connection Type **A:** ¼" Stepped Hose Barb **E:** 1½" Sanitary Flange **Q:** Single Step ½" Hose Barb

Ordering Information

Type		Size		Pore Size		Inlet/Outlet		X	X	Sterility		Pack Size	
	Code		Code		Code		Code				Code		Code
AseptiVent TF	DTLX	1"	51	0.2 µm	01	¼" SHB	A			Non Sterile	1	1	01
		2"	52	0.45 µm	02	¼" MNPT	B			EO Sterile	2		
		5"	53			½" MNPT	C						
		8"	57			½" Hose Barb	D						
						1½" Sanitary Flange	E						
						¾" Sanitary Flange	S						
						Quick Connector	J						
						Single Step ½" Hose Barb	Q						
						Female Luer Lock	U						
						Male Luer Slip	W						
						¾" Hose Barb	N						
						⅜" Hose Barb	I						

Example:

DTLX	53	01	DD	X	X	1	01
------	----	----	----	---	---	---	----

For End Connection availability and dimensions with different sizes refer Pages 89-90.

AseptiVent TF Large Capsule Filters (5", 10", 20", 30")

AseptiVent TF PTFE large capsule filters offer absolute retention and wide chemical compatibility for sterile filtration of air/gases as well as aggressive solvents in large volume.



**Microbially Validated
as per ASTM F 838**

**Complies with
USFDA 21 CFR 210.3 (b) (6)**


**Meets and Exceeds
USFDA 21 CFR 177.1520**

Specifications

Construction				
Membrane	Hydrophobic PTFE			
Support (Drainage) Layers	Polypropylene			
Plastic Parts	Polypropylene			
Pore Size	0.2 µm		0.45 µm	
Integrity Testing/Retention				
Bubble Point (with 70% IPA Wetted)	≥22psi(1.55Kg/cm ²)		≥ 10 psi (0.7 Kg/cm ²)	
Microbial Retention (LRV >7 for)	<i>Brevundimonas diminuta</i> (ATCC 19146) per cm ²		<i>Serratia marcescens</i> (ATCC 14756) per cm ²	
Size				
Size	5"	10"	20"	30"
Effective Filtration Area (Nominal)	3000 cm ²	6000 cm ²	12000 cm ²	18000 cm ²
Vent and Drain	¼" Hose Barb with Silicone 'O' ring			
Operational				
Max. Operating Temperature	80 °C @ ≤ 30 psi (2 Kg/cm ²)			
Max. Differential Pressure	60 psi (4 Kg/cm ²) @ 30 °C			
Sterilization	By Gas	Sterilizable by Ethylene Oxide		
	By Autoclave	Autoclavable at 125 °C for 30 minutes, 30 cycles. Cannot be In-line steam sterilized		
Shelf Life	3 years after Ethylene Oxide sterilization			

Air Flow Rates

0.2 µm AseptiVent TF, 10" Large Capsule Filters



End Connection Type E: 1½" Sanitary Flange

Ordering Information

Type		Size		Pore Size		Inlet/Outlet		X	Inline / T-line		Sterility		Pack Size	
	Code		Code		Code		Code			Code		Code		Code
AseptiVent TF	LTLX	5"*	53	0.2 µm	01	1½" Sanitary Flange	E		Inline	X	Non Sterile	1	1	01
		10"	54	0.45 µm	02	Single Step ½" Hose Barb	Q		T-line	T	EO Sterile	2		
		20"	55			¾" Sanitary Flange	S							
		30"	56			⅜" Hose Barb	I							
						1" Hose Barb	Z							

Example:

LTLX	54	01	EE	X	X	1	01
------	----	----	----	---	---	---	----

* Size 5" is available in In-line Capsule Filters Only

For End Connection availability and dimensions with different sizes refer Pages 89-90.

AseptiSure TH Mini Cartridge filters

mdi AseptiSure TH PTFE membrane mini cartridge filters are specially designed high temperature resistant PTFE filters which are steam sterilizable at upto 135°C. These filters are validated with liquid microbial challenge test as per ASTM F 838 to offer absolute retention even under high moisture conditions.

These are also validated for other key performance parameters such as chemical compatibility, extractables, heat stability, flow rates, blow through and ability to withstand accidental reverse pressure. These are available in a variety of pore sizes to suit specific microfiltration needs in critical and specialized process applications for air as well as liquid.

**Microbiologically Validated
as per ASTM F 838**

**Complies with
USFDA 21 CFR 210.3 (b) (6)**

**Meets and Exceeds
USFDA 21 CFR 177.1520**


Specifications

Construction			
Pore Size		0.2 µm	0.45 µm
Membrane	Hydrophobic PTFE		
Support (Drainage) Layers	Polypropylene		
Plastic Parts	Polypropylene		
Integrity Testing/Retention			
Bubble Point (70% IPA/Water)		22psi (1.54 kg/cm ²)	10psi (0.7 kg/cm ²)
Water Intrusion Rate	2.5"	≤ 0.3 ml/min @ 2.0kg/cm ²	–
	5"	≤ 0.6 ml/min @ 2.0kg/cm ²	–
Microbial Bacterial Retention (LRV >7 for)		<i>Brevundimonas diminuta</i> (ATCC 19146) per cm ²	<i>Serratia marcescens</i> (ATCC 14756) per cm ²
Size			
Size		2.5"	5"
Effective Filtration Area (Nominal)		1000 cm ²	2000 cm ²
Operational			
Max. Operating Temperature	95 °C @ ≤ 2 Kg/cm ² (30 psi)		
Max. Differential Pressure	3.5 Kg/cm ² (50 psi) @ 25 °C		
Reverse Pressure	≤ 0.7 Kg/cm ² (10 psi) @ 25 °C		
Sterilization	Autoclavable/In-line steam sterilizable at 135 °C @ maximum differential pressure of 5 psi (0.35 kg/cm ²) for 30 minutes, 80 cycles		



Air Flow Rates

0.2 µm AseptiSure TH Mini Cartridge Filters



Ordering Information

Type		Size		Pore Size		Adapter		Elastomer		Sterility		Pack Size	
	Code		Code		Code		Code		Code		Code		Code
AseptiSure TH	CPH	2.5"	50	0.2 µm	01	4463	E0	Silicone	SS	Non Sterile	1	1	01
		5"	53	0.45 µm	02	4463B	H0						
						4440	U0						
						Seal-K	G0*						
						Seal-O	F0						
						Seal-M	J0						

*G0 adapter is not available with elastomer. Please mention XX in place of elastomer code while ordering

Example:

CPH	50	01	E0	SS	1	01
-----	----	----	----	----	---	----

For Adapters and Elastomers availability refer Page 91.

AseptiSure TH Standard Cartridge Filters

**Microbiologically Validated
as per ASTM F 838**

**Complies with
USFDA 21 CFR 210.3 (b) (6)**

**Meets and Exceeds
USFDA 21 CFR 177.1520**

AseptiSure TH are high temperature resistant PTFE membrane cartridge filters which are steam sterilizable at 135°C. These filters are validated with liquid microbial challenge test as per ASTM F 838 to offer absolute retention even under high moisture conditions.


Specifications

Construction				
Final Filter Pore Size	0.2 µm	0.45 µm		
Membrane	Hydrophobic PTFE			
Support (Drainage) Layers	Polypropylene			
Plastic Parts	Polypropylene			
Integrity Testing/Retention				
Air Diffusion Flow (with 70% IPA Wetted) (10" Cartridge Filter)	≤ 45 ml/min @ 16 psi (1.12 Kg/cm ²)	≤ 45 ml/min @ 8 psi (0.56 Kg/cm ²)		
Microbial Retention (LRV >7 for)	<i>Brevundimonas diminuta</i> (ATCC 19146) per cm ²	<i>Serratia marcescens</i> (ATCC 14756) per cm ²		
Size				
Size	5"	10"	20"	30"
Effective Filtration Area (Nominal)	3000 cm ²	6000 cm ²	12000 cm ²	18000 cm ²
Operational				
Max. Operating Temperature	95 °C @ ≤ 30 psi (2 Kg/cm ²)			
Max. Differential Pressure	50 psi (3.5 Kg/cm ²) @ 25 °C			
Reverse Pressure	≤ 10 psi (0.7 Kg/cm ²) @ 25 °C			
Sterilization	Autoclavable/In-line steam sterilizable at 135 °C @ maximum differential pressure of 5 psi (0.35 Kg/cm ²) for 30minutes, 80 cycles .			



Air Flow Rates

0.2 µm AseptiSure TH Cartridge Filters



Ordering Information

Type		Size		Pore Size		Adapter		Elastomer		Sterility		Pack Size	
	Code		Code		Code		Code		Code		Code		Code
AseptiSure TH	CPTH	5"	53	0.2 µm	01	7P	A0	Silicone	SS	Non Sterile	1	1	01
		10"	54	0.45 µm	02	7P without fin	A1	Viton	SV				
		20"	55			28 with fin	C0	EPDM	SE				
		30"	56			'O'	D0	FEP Encapsulated Viton	FV				

Example:

CPTH	56	01	A0	SS	1	01
------	----	----	----	----	---	----

For Adapters and Elastomers availability refer Page 91.

AseptiSure TF Mini Cartridge filters

**Microbiologically Validated
as per ASTM F 838**
**Complies with
USFDA 21 CFR 210.3 (b) (6)**
**Meets and Exceeds
USFDA 21 CFR 177.1520**

mdi AseptiSure TF PTFE membrane mini cartridge filters are hydrophobic filters offering absolute retention. These filters are designed for sterilizing filtration of gases. The high quality of membrane and design of cartridge assures long life and ability to withstand adverse process conditions experienced during use.

mdi AseptiSure TF filters are validated for key performance parameters such as retention efficiency, chemical compatibility, extractables, heat stability, flow rates and blow through.


Specifications

Construction		
Final Filter Pore Size	0.2 µm	0.45 µm
Membrane	Hydrophobic PTFE	
Support (Drainage) Layers	Polypropylene	
Plastic Parts	Polypropylene	
Integrity Testing/Retention		
Bubble Point	22psi (1.52 Bar) with 70% IPA/Water Solution	10psi (0.69 Bar) with 70% IPA/Water Solution
Water Intrusion Rate	2.5"	≤ 0.3 ml/min @ 2.0kg/cm ²
	5"	≤ 0.6 ml/min @ 2.0kg/cm ²
Microbial Bacterial Retention (LRV >7 for)	<i>Brevundimonas diminuta</i> (ATCC 19146) per cm ²	<i>Serratia marcescens</i> (ATCC 14756) per cm ²
Size		
Size	2.5"	5"
Effective Filtration Area (Nominal)	1000 cm ²	2000 cm ²
Operational		
Max. Operating Temperature	80 °C @ ≤ 2 Kg/cm ² (30 psi)	
Max. Differential Pressure	3.5 Kg/cm ² (50 psi) @ 25 °C	
Reverse Pressure	≤ 0.7 Kg/cm ² (10 psi) @ 25 °C	
Sterilization	Autoclavable/In-line steam sterilizable at 121°C @ maximum differential pressure of 3 psi (0.21 kg/cm ²) for 30 minutes, 100 cycles	



Air Flow Rates

0.2 µm AseptiSure TF Mini Cartridge Filters



Ordering Information

Type		Size		Pore Size		Adapter		Elastomer		Sterility		Pack Size	
	Code		Code		Code		Code		Code		Code		Code
AseptiSure TF	CPTF	2.5"	50	0.2 µm	01	4463	E0	Silicone	SS	Non Sterile	1	1	01
		5"	53	0.45 µm	02	4463B	H0						
						4440	U0						
						Seal-K	G0*						
						Seal-O	F0						
						Seal-M	J0						

*G0 adapter is not available with elastomer. Please mention XX in place of elastomer code while ordering

Example:

CPTF	50	01	E0	SS	1	01
------	----	----	----	----	---	----

For Adapters and Elastomers availability refer Page 91.

AseptiSure TF Standard Cartridge Filters

**Microbially Validated
as per ASTM F 838-05**

**Complies with
USFDA 21 CFR 210.3 (b) (6)**

**Meets and Exceeds
USFDA 21 CFR 177.1520**

AseptiSure TF cartridge filters employ hydrophobic PTFE membrane offering absolute retention, wide chemical compatibility, and are validated with liquid bacterial challenge test.


Specifications

Construction				
Final Filter Pore Size	0.2 µm	0.45 µm		
Membrane	Hydrophobic PTFE			
Support (Drainage) Layers	Polypropylene			
Body and Core	Polypropylene			
Integrity Testing/Retention				
Air Diffusion Flow (with 70% IPA Wetted) (10" Cartridge Filter)	≤ 45 ml/min @ 16 psi (1.12 Kg/cm ²)	≤ 45 ml/min @ 8 psi (0.56 Kg/cm ²)		
Microbial Retention (LRV >7 for)	<i>Brevundimonas diminuta</i> (ATCC 19146) per cm ²	<i>Serratia marcescens</i> (ATCC 14756) per cm ²		
Size				
Size	5"	10"	20"	30"
Effective Filtration Area (Nominal)	3000 cm ²	6000 cm ²	12000 cm ²	18000 cm ²
Operational				
Max. Operating Temperature	80 °C @ ≤ 30 psi (2 Kg/cm ²)			
Max. Differential Pressure	50 psi (3.5 Kg/cm ²) @ 25 °C			
Reverse Pressure	≤ 10 psi (0.7 Kg/cm ²) @ 25 °C			
Sterilization	Autoclavable/In-line steam sterilizable at 121°C @ maximum differential pressure of 3 psi (0.21 kg/cm ²) for 30 minutes, 100 cycles			



Air Flow Rates

0.2 µm AseptiSure TF Cartridge Filters



Ordering Information

Type		Size		Pore Size		Adapter		Elastomer		Sterility		Pack Size	
	Code		Code		Code		Code		Code		Code		Code
AseptiSure TF	CPTF	5"	53	0.2 µm	01	7P	A0	Silicone	SS	Non Sterile	1	1	01
		10"	54	0.45 µm	02	7P without fin	A1	Viton	SV				
		20"	55			28 with fin	C0	EPDM	SE				
		30"	56			'O'	D0	FEP Encapsulated Viton	FV				

Example:

CPTF	54	01	A0	FV	1	01
------	----	----	----	----	---	----

For Adapters and Elastomers availability refer Page 91.

Filters for Clarification and Pre-filtration

mdi offers a range of pre-filters designed to protect the terminal sterilizing grade membrane filters and maximize throughputs.

These are biologically inert filters with wide chemical compatibility, offering very high retention efficiency and are available in cartridge filter and capsule filter formats, in different sizes, pore sizes and end connections to suit different needs.

These filter devices are available as:

Filter Type	Product Name
Polyethersulfone Membrane Capsule Filters with Microglassfiber Upstream	<i>ClariPro GK</i>
Polyethersulfone Membrane Cartridge Filters with Microglassfiber Upstream	<i>ClariSure GK</i>
Microglassfiber Capsule Filter	<i>ClariCap GS</i>
Microglassfiber Cartridge Filters	<i>ClariSure GS</i> <i>ClariSure GP</i>
Polypropylene Capsule Filters	<i>ClariCap PP</i>
Polypropylene Cartridge Filters	<i>ClariSure PA</i>

Applications

- Precipitate removal post viral inactivation
- Pre-filtration of cell culture media
- Pre-filtration of serum and other viscous biologicals
- Pre-filtration of serum solutions
- Clarification of cell harvest supernatant
- Pre-filtration of protein solutions
- Pre-filtration of high value difficult to filter drug solutions
- Pre-filtration of large volume parenterals
- Pre-filtration of difficult to filter SVP
- Polishing of turbid solutions
- Pre-filtration of fermentor air

Quality Assurance

These filter devices are manufactured in Class 10,000 clean rooms under ISO 9001:2015 certified quality management systems and are deeply validated to meet compendia and regulatory requirements.

Assurance	
Toxicity	Passes Bioreactivity test, In Vivo, as per USP <88> for Class VI plastics
Non Fiber Releasing	Passes test as per USP and comply with USFDA 21 CFR Part 210.3 (b)(6) for fiber release
Extractables with WFI	Passes test as per USP <661>
Oxidizable Substances	Within limits as specified in USP <1231>
Particle Shedding	Passes test as per USP <788> for particulate matter in injection
Indirect Food Additive	All Polypropylene components meet the FDA Indirect Food Additive requirements cited in 21 CFR 177.1520
Good Manufacturing Practice	These products are manufactured in a facility which adheres to Good Manufacturing Practices

Filter Selection Chart

Application Area	Key Application Requirements	Capsule Filters		Cartridge Filters	
Biopharmaceuticals					
Precipitate removal post viral inactivation	<ul style="list-style-type: none"> - High retention efficiency - High throughput 	ClariPro GK PES Membrane Capsule Filters with Microglassfiber upstream	ClariCap GS Microglassfiber Capsule Filters	ClariSure GK PES Membrane Cartridge Filters with Microglassfiber upstream	ClariSure GS Microglassfiber Cartridge Filters
Pre-filtration of cell culture media	<ul style="list-style-type: none"> - High retention efficiency - High throughput 	ClariPro GK PES Membrane Capsule Filters with Microglassfiber upstream	ClariCap GS Microglassfiber Capsule Filters	ClariSure GK PES Membrane Cartridge Filters with Microglassfiber upstream	ClariSure GS Microglassfiber Cartridge Filters
Prefiltration of serum and other viscous biologicals	<ul style="list-style-type: none"> - High retention efficiency - High throughput 	ClariPro GK PES Membrane Capsule Filters with Microglassfiber upstream	ClariCap GS Microglassfiber Capsule Filters	ClariSure GK PES Membrane Cartridge Filters with Microglassfiber upstream	ClariSure GS Microglassfiber Cartridge Filters
Pre-filtration of serum solutions	<ul style="list-style-type: none"> - High retention efficiency - High throughput 	ClariPro GK PES Membrane Capsule Filters with Microglassfiber upstream	ClariCap GS Microglassfiber Capsule Filters	ClariSure GK PES Membrane Cartridge Filters with Microglassfiber upstream	ClariSure GS Microglassfiber Cartridge Filters
Clarification of cell harvest supernatants	<ul style="list-style-type: none"> - High retention efficiency - High throughput 	ClariPro GK PES Membrane Capsule Filters with Microglassfiber upstream	ClariCap GS Microglassfiber Capsule Filters	ClariSure GK PES Membrane Cartridge Filters with Microglassfiber upstream	ClariSure GS Microglassfiber Cartridge Filters
Pre-filtration of fermentor air	<ul style="list-style-type: none"> - High retention efficiency 	-	ClariCap PP Polypropylene Capaule Filters	ClariSure PA Pleated Polypropylene Cartridge Filters	-
Pre-filtration of proteinaceous liquids	<ul style="list-style-type: none"> - Low hold up volume - High throughput 	ClariPro GK PES Membrane Capsule Filters with Microglassfiber upstream	ClariCap PP Polypropylene Capaule Filters	-	-
Pharmaceuticals					
Pre-filtration of high value difficult to filter drug solutions	<ul style="list-style-type: none"> - High retention efficiency - High throughput 	ClariPro GK PES Membrane Capsule Filters with Microglassfiber upstream	ClariCap GS Microglassfiber Capsule Filters	ClariSure GK PES Membrane Cartridge Filters with Microglassfiber upstream	ClariSure GS Microglassfiber Cartridge Filters
Pre-filtration of large volume parenterals	<ul style="list-style-type: none"> - High retention efficiency - High throughput 	-	-	ClariSure PA Pleated Polypropylene Cartridge Filters	ClariSure GP Microglassfiber Cartridge Filters
Pre-filtration of difficult to filter SVP like Oxytetracycline	<ul style="list-style-type: none"> - High retention efficiency - High throughput 	-	-	ClariSure GK PES Membrane Cartridge Filters with Microglassfiber upstream	ClariSure GP Microglassfiber Cartridge Filters
Polishing of turbid solutions	<ul style="list-style-type: none"> - Very high retention efficiency for colloidal particles - High throughput 	ClariPro GK PES Membrane Capsule Filters with Microglassfiber upstream	ClariCap GS Microglassfiber Capsule Filters	ClariSure GK PES Membrane Cartridge Filters with Microglassfiber upstream	ClariSure GS Microglassfiber Cartridge Filters

ClariPro GK Small Capsule Filters (1", 2", 5", 8")

Complies with
USFDA 21 CFR 210.3 (b)(6)
Meets and Exceeds
USFDA 21 CFR 177.1520

ClariPro GK hydrophilic PES membrane capsule filters are ready to use, disposable filtration devices. These filters are specially designed filters incorporating a microglassfiber upstream layer and a downstream PES membrane layer and are used as pre-filters in biopharmaceuticals process development as well as manufacturing processes for difficult to filter solutions.

Radiation Sterilizable: ClariPro GK- γ
Autoclavable: ClariPro GK

Specifications

Construction	
Pore Size	0.1 μ m, 0.2 μ m, 0.5 μ m
Membrane	Hydrophilic PES
Pre-filter	Microglassfiber
Support (Drainage) Layers	Polyester
Plastic Parts	Polypropylene


Size				
Size	1"	2"	5"	8"
Effective Filtration Area (Nominal)	150 cm ²	500 cm ²	1000 cm ²	1500 cm ²
Vent and Drain	1/4" Hose Barb with Silicone 'O' ring			

Operational	
Max. Operating Temperature	80 °C @ \leq 30 psi (2 Kg/cm ²)
Max. Differential Pressure	60 psi (4 Kg/cm ²) @ 30 °C
Sterilization	By Irradiation ClariPro GK- γ : Gamma Irradiatable up to 50 kGy. These filters should not be autoclaved or in-line steam sterilized.
	By Gas ClariPro GK: Sterilizable by Ethylene Oxide
	ClariPro GK: Autoclavable at 125 °C for 30 minutes, 25 cycle.
These cannot be In-line steam sterilized	



Water Flow Rates

0.5 μ m ClariPro GK Capsule Filters



E: 1 1/2" Sanitary Flange Connections

Ordering Information

Type		Size		Pore Size		Inlet/Outlet		Radiation Sterilizable		X	Sterility		Pack Size	
	Code		Code		Code		Code		Code				Code	
ClariPro GK	DGKX	1"	51	0.1 μ m	36	1/4" SHB	A	Yes	R		Non Sterile	1	1	01
		2"	52	0.2 μ m	01	1/2" Hose Barb	D	No*	X		EO Sterile	2		
		5"	53	0.5 μ m	04	1 1/2" Sanitary Flange	E				Gamma Sterile	3		
		8"	57			3/4" Sanitary Flange	S							
						Quick Connector	J							
						Single Step 1/2" Hose Barb	Q							
						Female Luer Lock	U							
						Male Luer Slip	W							
						3/16" Hose Barb	N							
						3/8" Hose Barb	I							

Example:

DGKX	57	04	DD	R	X	1	01
------	----	----	----	---	---	---	----

*Gamma Sterile capsule filters cannot be Gamma Irradiated again

For End Connection availability and dimensions with different sizes refer Pages 89-90.

ClariPro GK Large Capsule Filters (5", 10", 20", 30")

Complies with
USFDA 21 CFR 210.3 (b)(6)

Meets and Exceeds
USFDA 21 CFR 177.1520

ClariPro GK hydrophilic PES membrane large capsule filters are ready to use, disposable filtration devices with a microglassfiber upstream layer and a downstream PES membrane final layer.

The microglassfiber layer offers high dirt holding capacity and efficient retention of colloidal particles to give clear downstream in case of difficult to filter turbid solutions. The downstream PES membrane offers absolute retention and low protein binding along with high flow rates.

Radiation Sterilizable: ClariPro GK-γ

Autoclavable: ClariPro GK

Specifications

Construction	
Pore Size	0.1 μm, 0.2 μm, 0.5 μm
Membrane	Hydrophilic PES
Pre-filter	Microglassfiber
Support (Drainage) Layers	Polyester
Plastic Parts	Polypropylene

Size				
Size	5"	10"	20"	30"
Effective Filtration Area (Nominal)	2500 cm ²	5000 cm ²	10000 cm ²	15000 cm ²

Vent and Drain ¼" Hose Barb with Silicone 'O' ring

Operational	
Max. Operating Temperature	80 °C @ ≤ 30 psi (2 Kg/cm ²)
Max. Differential Pressure	60 psi (4 Kg/cm ²) @ 30 °C

ClariPro GK-γ: Gamma Irradiatable up to 50 kGy. These filters should not be autoclaved or in-line steam sterilized.

Sterilization By Gas ClariPro GK: Sterilizable by Ethylene Oxide


By Autoclave ClariPro GK: Autoclavable at 125 °C for 30 minutes, 25 cycles

These cannot be In-line steam sterilized



Water Flow Rates

0.5 μm ClariPro GK Large Capsule Filters



E: 1½" Sanitary Flange Connections

Ordering Information

Type		Size		Pore Size		Inlet/Outlet		Radiation Sterilizable		Inline / T-line		Sterility		Pack Size	
	Code		Code		Code		Code		Code		Code		Code		Code
ClariProGK	LGKX	5"	53	0.1 μm	36	1½" Sanitary Flange	E	Yes	R	Inline	X	Non Sterile	1	1	01
		10"	54	0.2 μm	01	Single Step ½" Hose Barb	Q	No*	X	T-line	T	EO Sterile	2		
		20"	55	0.5 μm	04	¾" Sanitary Flange	S					Gamma Sterile	3		
		30"	56			⅝" Hose Barb	I								
						1" Hose Barb	Z								

Example:

LGKX	54	01	EE	R	T	1	01
------	----	----	----	---	---	---	----

* Gamma Sterile capsule filters cannot be Gamma Irradiated again

** Size 5" is available in In-line Capsule Filters Only

For End Connection availability and dimensions with different sizes refer Pages 89-90.

ClariSure GK Mini Cartridge Filters

Complies with
USFDA 21 CFR 210.3 (b)(6)

Meets and Exceeds
USFDA 21 CFR 177.1520

ClariSure GK hydrophilic PES membrane mini cartridge filters offer a microglassfiber upstream layer with a PES membrane final filter to combine high dirt holding capacities with efficient retention of colloidal particles to give clear downstream in case of difficult to filter turbid solutions.


Specifications

Construction	
Pore Size	0.1 µm, 0.2 µm, 0.5 µm
Membrane	Hydrophilic PES
Pre-filter	Microglassfiber
Support (Drainage) Layers	Polyester
Plastic Parts	Polypropylene
Size	
Size	2.5" 5"
Effective Filtration Area (Nominal)	1000 cm ² 1500 cm ²
Operational	
Max. Operating Temperature	80 °C @ ≤ 30 psi (2 Kg/cm ²)
Max. Differential Pressure	50 psi (3.5 Kg/cm ²) @ 25 °C
Sterilization	Autoclavable/In-line steam sterilizable at 121°C for 30 minutes, 25 cycles



Water Flow Rates

0.5 µm ClariSure GK Mini Cartridge Filters



Ordering Information

Type		Size		Pore Size		Adapter		Elastomer		Sterility		Pack Size	
	Code		Code		Code		Code		Code		Code		Code
ClariSure GK	CGKX	2.5"	50	0.1 µm	36	4463	E0	Silicone	SS	Non Sterile	1	1	01
		5"	53	0.2 µm	01	4463B	H0						
				0.5 µm	04	4440	U0						
						Seal-K	G0*						
						Seal-O	F0						
						Seal-M	J0						

*G0 adapter is not available with elastomer. Please mention **XX** in place of elastomer code while ordering

Example:

CGKX	50	04	E0	SS	1	01
------	----	----	----	----	---	----

For Adapters and Elastomers availability refer Page 91.

ClariSure GK Standard Cartridge Filters

Complies with
USFDA 21 CFR 210.3 (b)(6)

Meets and Exceeds
USFDA 21 CFR 177.1520

ClariSure GK hydrophilic PES membrane cartridge filters offer a microglassfiber upstream layer with a PES membrane final filter to combine high dirt holding capacities with efficient retention of colloidal particles to give clear downstream in case of difficult to filter turbid solutions.




Specifications

Construction	
Pore Size	0.1 µm, 0.2 µm, 0.5 µm
Membrane	Hydrophilic PES
Pre-filter	Microglassfiber
Support (Drainage) Layers	Polyester
Plastic Parts	Polypropylene
Size	
Size	5" 10" 20" 30"
Effective Filtration Area (Nominal)	2500 cm ² 5000 cm ² 10000 cm ² 15000 cm ²
Operational	
Max. Operating Temperature	80 °C @ ≤ 30 psi (2 Kg/cm ²)
Max. Differential Pressure	50 psi (3.5 Kg/cm ²) @ 25 °C
Sterilization	Autoclavable/In-line steam sterilizable at 121°C for 30 minutes, 25 cycles

Water Flow Rates

0.5 µm ClariSure GK Cartridge Filters



Ordering Information

Type		Size		Pore Size		Adapter		Elastomer		Sterility		Pack Size	
	Code		Code		Code		Code		Code		Code		Code
ClariSure GK	CGKX	5"	53	0.1 µm	36	7P	A0	Silicone	SS	Non Sterile	1	1	01
		10"	54	0.2 µm	01	7P without fin	A1	Viton	SV				
		20"	55	0.5 µm	04	28 with fin	C0	EPDM	SE				
		30"	56			'O'	D0	FEP Encapsulated Viton	FV				

Example:

CGKX	56	36	A0	SS	1	01
------	----	----	----	----	---	----

For Adapters and Elastomers availability refer Page 91.

ClariCap GS Small Capsule Filters (1", 2", 5", 8")

Complies with
USFDA 21 CFR 210.3 (b)(6)

Meets and Exceeds
USFDA 21 CFR 177.1520

ClariCap GS capsule filters employ microglassfiber filter media for efficient retention of colloidal particles to give clear downstream in case of difficult to filter turbid solutions. These specially designed filtration devices are non media migrating with a heat calendered polypropylene layer in the downstream.


Specifications

Construction				
Pore Size	0.7 µm, 1 µm, 1.5 µm, 2 µm			
Filter Media	Microglassfiber			
Plastic Parts	Polypropylene			
Size				
Size	1"	2"	5"	8"
Effective Filtration Area (Nominal)	150 cm ²	400 cm ²	800 cm ²	1200 cm ²
Vent and Drain	¼" Hose Barb with Silicone 'O' ring			
Operational				
Max. Operating Temperature	80 °C @ ≤ 30 psi (2 Kg/cm ²)			
Max. Differential Pressure	60 psi (4 Kg/cm ²) @ 30 °C			
Sterilization	Autoclavable at 125 °C for 30 minutes, 30 cycles. Cannot be In-line steam sterilized			



Water Flow Rates

0.7 µm ClariCapGS Capsule Filters



End Connection : D: ½" Hose Barb

Ordering Information

Type		Size		Pore Size		Inlet/Outlet		X	Bell		Sterility		Pack Size	
	Code		Code		Code		Code			Code		Code		Code
ClariCap GS	DGSX	1"	51	0.7 µm	41	¼" SHB	A		Yes	B	Non Sterile	1	1	01
		2"	52	1 µm	05	½" Hose Barb	D		No Bell	X				
		5"	53	1.5 µm	14	1½" Sanitary Flange	E							
		8"	57	2 µm	15	¾" Sanitary Flange	S							
						Quick Connector	J							
						Single Step ½" Hose Barb	Q							
						Female Luer Lock	U							
						Male Luer Slip	W							
						⅜" Hose Barb	N							
						⅝" Hose Barb	I							

Example:

DGSX	57	05	DD	X	X	1	01
------	----	----	----	---	---	---	----

For End Connection availability, Bell and dimensions with different sizes refer Pages 89-90.

ClariCap GS Large Capsule Filters (5", 10", 20", 30")

**Complies with
USFDA 21 CFR 210.3 (b)(6)**

**Meets and Exceeds
USFDA 21 CFR 177.1520**

ClariCap GS large capsule filters are multilayered, high throughput filters, specially designed for difficult to filter solutions.

These are high efficiency pre-filters combining the unique abilities of microglassfiber filter media to retain colloidal particles and heat calendered polypropylene filter media to ensure non media migration.


Specifications

Construction				
Pore Size	0.7 µm, 1 µm, 1.5 µm, 2 µm			
Filter Media	Microglassfiber			
Plastic Parts	Polypropylene			
Size				
Size	5"	10"	20"	30"
Effective Filtration Area (Nominal)	1700 cm ²	3400 cm ²	6800 cm ²	10200 cm ²
Vent and Drain	¼" Hose Barb with Silicon 'O' ring			
Operational				
Max. Operating Temperature	80 °C @ ≤ 30 psi (2 Kg/cm ²)			
Max. Differential Pressure	60 psi (4 Kg/cm ²) @ 30 °C			
Sterilization	Autoclavable at 125 °C for 30 minutes, 30 cycles. Cannot be In-line steam sterilized			



Water Flow Rates

0.7µm ClariCap GS Large Capsule Filters



E: 1½" Sanitary Flange Connections

Ordering Information

Type		Size		Pore Size		Inlet/Outlet		X	Inline / T-line		Sterility		Pack Size	
	Code		Code		Code		Code			Code		Code		Code
ClariCap GS	LGSX	5"*	53	0.7 µm	41	1½" Sanitary Flange	E			Code	Non Sterile	1	1	01
		10"	54	1 µm	05	Single Step ½" Hose Barb	Q		Inline	X				
		20"	55	1.5 µm	14	¾" Sanitary Flange	S		T-line	T				
		30"	56	2 µm	15	3/8" Hose Barb	I							
						1" Hose Barb	Z							

Example:

LGSX	54	41	EE	X	T	1	01
------	----	----	----	---	---	---	----

*Size 5" is available in Inline capsule filters only

For End Connection availability and dimensions with different sizes refer Page 89-90.

ClariSure GS Mini Cartridge Filters

Complies with
USFDA 21 CFR 210.3 (b)(6)

Meets and Exceeds
USFDA 21 CFR 177.1520

ClariSure GS mini cartridge filters are multilayered, high throughput filters, specially designed for difficult to filter solutions.

It is a very high efficiency pre-filter combining the unique abilities of microglassfiber filter media to retain colloidal particles and heat calendered polypropylene filter media to ensure non media migration.




Specifications

Construction	
Pore Size	0.7 µm, 1 µm
Filter Media	Microglassfiber
Support (Drainage) Layers	Polyester
Plastic Parts	Polypropylene
Size	
Size	2.5" 5"
Effective Filtration Area (Nominal)	800 cm ² 1200 cm ²
Operational	
Max. Operating Temperature	80 °C @ ≤ 30 psi (2 Kg/cm ²)
Max. Differential Pressure	50 psi (3.5 Kg/cm ²) @ 25 °C
Sterilization	Autoclavable/In-line steam sterilizable at 121°C for 30 minutes, 30 cycles

Water Flow Rates

0.7 µm ClariSure GS Mini Cartridge Filters



Ordering Information

Type		Size		Pore Size		Adapter		Elastomer		Sterility		Pack Size	
	Code		Code		Code		Code		Code		Code		Code
ClariSure GS	CPGS	2.5"	50	0.7 µm	41	4463	E0	Silicone	SS	Non Sterile	1	1	01
		5"	53	1 µm	05	4463B	H0						
						4440	U0						
						Seal-K	G0*						
						Seal-O	F0						
						Seal-M	J0						

*G0 adapter is not available with elastomer. Please mention XX in place of elastomer code while ordering

Example:

CPGS	50	41	E0	SS	1	01
------	----	----	----	----	---	----

For Adapters and Elastomers availability refer Page 91.

ClariSure GS Standard Cartridge Filters

Complies with
USFDA 21 CFR 210.3 (b)(6)

Meets and Exceeds
USFDA 21 CFR 177.1520

ClariSure GS cartridge filters are multilayered, high throughput filters, specially designed for difficult to filter solutions.

It is a very high efficiency pre-filter combining the unique abilities of microglassfiber filter media to retain colloidal particles and heat calendered polypropylene filter media to ensure non media migration.


Specifications

Construction				
Pore Size	0.7 µm, 1 µm			
Filter Media	Microglassfiber			
Support (Drainage) Layers	Polyester			
Plastic Parts	Polypropylene			
Size				
Size	5"	10"	20"	30"
Effective Filtration Area (Nominal)	1700 cm ²	3400 cm ²	6800 cm ²	10200 cm ²
Operational				
Max. Operating Temperature	80 °C @ ≤ 30 psi (2 Kg/cm ²)			
Max. Differential Pressure	50 psi (3.5 Kg/cm ²) @ 30 °C			
Sterilization	Autoclavable/In-line steam sterilizable at 121 °C for 30 minutes, 30 cycles			



Water Flow Rates

0.7 µm ClariSure GS Cartridge Filters



Ordering Information

Type		Size		Pore Size		Adapter		Elastomer		Sterility		Pack Size	
	Code		Code		Code		Code		Code		Code		Code
ClariSure GS	CPGS	5"	53	0.7 µm	41	7P	A0	Silicone	SS	Non Sterile	1	1	01
		10"	54	1 µm	05	7P without fin	A1	Viton	SV				
		20"	55			28 with fin	C0	EPDM	SE				
		30"	56			'O'	D0	FEP Encapsulated Viton	FV				

Example:

CPGS	55	41	A0	SV	1	01
------	----	----	----	----	---	----

For Adapters and Elastomers availability refer Page 91.

ClariSure GP Mini Cartridge Filters

Complies with
USFDA 21 CFR 210.3 (b)(6)

Meets and Exceeds
USFDA 21 CFR 177.1520

ClariSure GP mini cartridge filters are specially designed multilayered precision filtration devices. These cartridge filters act as throughput enhancers, specially in case of difficult to filter solutions.

A microglassfiber upstream layer retains very fine colloidal particles and a downstream polypropylene layer checks any kind of media migration.

Specifications

Construction		
Pore Size	0.5 µm, 1 µm, 2 µm	
Filter Media	Microglassfiber	
Support (Drainage) Layers	Polypropylene	
Plastic Parts	Polypropylene	
Size		
Size	2.5"	5"
Effective Filtration Area (Nominal)	800 cm ²	1600 cm ²
Operational		
Max. Operating Temperature	80 °C @ ≤ 30 psi (2 Kg/cm ²)	
Max. Differential Pressure	50 psi (3.5 Kg/cm ²) @ 25 °C	
Sterilization	Autoclavable/In-line steam sterilizable at 121°C for 30 minutes, 30 cycles	



Ordering Information

Type		Size		Pore Size		Adapter		Elastomer		Sterility		Pack Size	
	Code		Code		Code		Code		Code		Code		Code
ClariSure GP	CPGP	2.5"	50	0.5 µm	04	4463	E0	Silicone	SS	Non Sterile	1	1	01
		5"	53	1 µm	05	4463B	H0						
				2 µm	15	4440	U0						
						Seal-K	G0*	*G0 adapter is not available with elastomer. Please mention XX in place of elastomer code while ordering					
						Seal-O	F0						
						Seal-M	J0						

Example:

CPGP	50	05	E0	SS	1	01
------	----	----	----	----	---	----

For Adapters and Elastomers availability refer Page 91.

ClariSure GP Standard Cartridge Filters

**Complies with
USFDA 21 CFR 210.3 (b)(6)**

**Meets and Exceeds
USFDA 21 CFR 177.1520**

ClariSure GP cartridge filters are specially designed multilayered precision filtration devices. These cartridge filters act as throughput enhancers, specially in case of difficult to filter solutions.

A microglassfiber upstream layer retains very fine colloidal particles and a downstream polypropylene layer checks any kind of media migration.




Specifications

Construction	
Pore Size	0.5 μm, 1 μm, 2 μm
Filter Media	Microglassfiber
Support (Drainage) Layers	Polypropylene
Plastic Parts	Polypropylene
Size	
Size	5" 10" 20" 30"
Effective Filtration Area (Nominal)	2500 cm ² 5000 cm ² 10000 cm ² 15000 cm ²
Operational	
Max. Operating Temperature	80 °C @ ≤ 30 psi (2 Kg/cm ²)
Max. Differential Pressure	50 psi (3.5 Kg/cm ²) @ 25 °C
Sterilization	Autoclavable/In-line steam sterilizable at 121°C for 30 minutes, 30 cycles

Water Flow Rates

ClariSure GP, 10" Cartridge Filters



Ordering Information

Type		Size		Pore Size		Adapter		Elastomer		Sterility		Pack Size	
	Code		Code		Code		Code		Code		Code		Code
ClariSure GP	CPGP	5"	53	0.5 μm	04	7P	A0	Silicone	SS	Non Sterile	1	1	01
		10"	54	1 μm	05	7P without fin	A1	Viton	SV				
		20"	55	2 μm	15	28 with fin	C0	EPDM	SE				
		30"	56			'O'	D0	FEP Encapsulated Viton	FV				

Example:

CPGP	55	04	A0	SS	1	01
------	----	----	----	----	---	----

For Adapters and Elastomers availability refer Page 91.

ClariCap PP Small Capsule Filters (1", 2", 5", 8")

ClariCap PP capsule filters employ high retention efficiency polypropylene filter media with wide chemical compatibility for pre-filtration and polishing applications.

Complies with
USFDA 21 CFR 210.3 (b)(6)

Meets and Exceeds
USFDA 21 CFR 177.1520


Specifications

Construction				
Pore Size	0.5 µm, 1 µm, 2.5 µm, 5 µm, 10 µm, 20 µm			
Filter Media	Polypropylene			
Support (Drainage) Layers	Polypropylene			
Plastic Parts	Polypropylene			
Size				
Size	1"	2"	5"	8"
Effective Filtration Area (Nominal)	250 cm ²	500 cm ²	800 cm ²	1600 cm ²
Vent and Drain	¼" Hose Barb with Silicone 'O' ring			
Operational				
Max. Operating Temperature	80 °C @ ≤ 30 psi (2 Kg/cm ²)			
Max. Differential Pressure	60 psi (4 Kg/cm ²) @ 30 °C			
Sterilization	Autoclavable at 125 °C for 30 minutes, 30 cycles. Cannot be In-line steam sterilized			



Water Flow Rates

ClariCap PP, 5" Capsule Filters



Ordering Information

Type		Size		Pore Size		Inlet/Outlet		X	X	Sterility		Pack Size	
	Code		Code		Code		Code				Code		Code
ClariCap PP	DOLX	1"	51	0.5 µm	04	¼" SHB	A			Non Sterile	1	1	01
		2"	52	1 µm	05	¼" MNPT	B						
		5"	53	2.5 µm	06	½" MNPT	C						
		8"	57	5 µm	07	½" Hose Barb	D						
				10 µm	08	1½" Sanitary Flange	E						
				20 µm	11	¾" Sanitary Flange	S						
						Quick Connector	J						
						Single Step ½" Hose Barb	Q						
						Female Luer Lock	U						
						Male Luer Slip	W						
						¾" Hose Barb	N						
						¾" Hose Barb	I						

Example:

DOLX	53	05	CC	X	X	1	01
------	----	----	----	---	---	---	----

For End Connection availability, Bell and dimensions with different sizes refer Pages 89-90.

ClariCap PP Large Capsule Filters (5", 10", 20", 30")

**Complies with
USFDA 21 CFR 210.3 (b)(6)**

**Meets and Exceeds
USFDA 21 CFR 177.1520**

ClariCap PP large capsule filters employ high retention efficiency Polypropylene filter media for wide chemical compatibility, efficient pre-filtration for clarification and polishing applications.




Specifications

Construction	
Pore Size	0.5 µm, 1 µm, 2.5 µm, 5 µm, 10 µm, 20 µm
Filter Media	Polypropylene
Support (Drainage) Layers	Polypropylene
Plastic Parts	Polypropylene
Size	
Size	5" 10" 20" 30"
Effective Filtration Area (Nominal)	2500 cm ² 5000 cm ² 10000 cm ² 15000 cm ²
Vent and Drain	¼" Hose Barb with Silicone 'O' ring
Operational	
Max. Operating Temperature	80 °C @ ≤ 30 psi (2 Kg/cm ²)
Max. Differential Pressure	60 psi (4 Kg/cm ²) @ 30 °C
Sterilization	Autoclavable at 125 °C for 30 minutes, 30 cycles. Cannot be In-line steam sterilized

Water Flow Rates

ClariCap PP 10", Large Capsule Filters



End Connection Type:

E: 1½" Sanitary Flange

Ordering Information

Type		Size		Pore Size		Inlet/Outlet		X	Inline / T-line		Sterility		Pack Size	
	Code		Code		Code		Code			Code		Code		Code
ClariCap PP	LOLX	5"	53	0.5 µm	04	1½" Sanitary Flange	E				Non Sterile	1	1	01
		10"	54	1 µm	05	Single Step ½" Hose Barb	Q		Inline	X				
		20"	55	2.5 µm	06	¾" Sanitary Flange	S		T-line	T				
		30"	56	5 µm	07	⅝" Hose Barb	I							
				10 µm	08	1" Hose Barb	Z							
				20 µm	11									

Example:

LOLX	54	06	QQ	X	T	1	01
------	----	----	----	---	---	---	----

* Size 5" is available in In-line Capsule Filters Only

For End Connection availability and dimensions with different sizes refer Pages 89-90.

ClariSure PA Mini Cartridge Filters

Complies with
USFDA 21 CFR 210.3 (b)(6)

Meets and Exceeds
USFDA 21 CFR 177.1520

mdi ClariSure PA mini cartridge filters are 100% polypropylene, very high retention efficiency pleated cartridge filters offering large filtration area.

These filters have heat stable construction and are used as pre-filters to sterilizing membrane cartridge filters.




Specifications


Construction		
Pore Size	0.5 μm, 1 μm, 1.2 μm, 2.5 μm, 5 μm, 8 μm, 10 μm	
Filter Media	Polypropylene	
Support (Drainage) Layers	Polypropylene	
Plastic Parts	Polypropylene	
Size		
Size	2.5"	5"
Effective Filtration Area (Nominal)	800 cm ²	1600 cm ²
Operational		
Max. Operating Temperature	80 °C @ ≤ 30 psi (2 Kg/cm ²)	
Max. Differential Pressure	50 psi (3.5 Kg/cm ²) @ 25 °C	
Sterilization	Autoclavable/In-line steam sterilizable at 121°C for 30 minutes, 100 cycles	

Water Flow Rates

ClariSure PA, 2.5" Mini Cartridge Filters



ClariSure PA, 5" Mini Cartridge Filters



Ordering Information

Type		Size		Pore Size		Adapter		Elastomer		Sterility		Pack Size	
	Code		Code		Code		Code		Code		Code		Code
ClariSure PA	CPPA	2.5"	50	0.5 μm	04	4463	E0	Silicone	SS	Non Sterile	1	1	01
		5"	53	1 μm	05	4463B	H0						
				1.2 μm	10	4440	U0						
				2.5 μm	06	Seal-K	G0*						
				5 μm	07	Seal-O	F0						
				8 μm	17	Seal-M	J0						
				10 μm	08								

*G0 adapter does not come with elastomer. Please mention XX in place of elastomer while ordering

Example:

CPPA	50	05	E0	SS	1	01
------	----	----	----	----	---	----

For Adapters and Elastomers availability refer Page 91.

ClariSure PA Standard Cartridge Filters

Complies with
USFDA 21 CFR 210.3 (b)(6)

Meets and Exceeds
USFDA 21 CFR 177.1520

ClariSure PA cartridge filters are 100% polypropylene, very high retention efficiency pleated cartridge filters offering large filtration area.

These filters have heat stable construction and are used as pre-filters to sterilizing membrane cartridge filters.


Specifications

Construction				
Pore Size	0.5 µm, 1 µm, 1.2 µm, 2.5 µm, 5 µm, 8 µm, 10 µm			
Filter Media	Polypropylene			
Support (Drainage) Layers	Polypropylene			
Plastic Parts	Polypropylene			
Size				
Size	5"	10"	20"	30"
Effective Filtration Area (Nominal)	2500 cm ²	5000 cm ²	10000 cm ²	15000 cm ²
Operational				
Max. Operating Temperature	80 °C @ ≤ 30 psi (2 Kg/cm ²)			
Max. Differential Pressure	50 psi (3.5 Kg/cm ²) @ 25 °C			
Sterilization	Autoclavable/In-line steam sterilizable at 121 °C for 30 minutes, 100 cycles			



Water Flow Rates

ClariSure PA, 10" Cartridge Filters



Ordering Information

Type		Size		Pore Size		Adapter		Elastomer		Sterility		Pack Size	
	Code		Code		Code		Code		Code		Code		Code
ClariSure PA	CPPA	5"	53	0.5 µm	04	7P	A0	Silicone	SS	Non Sterile	1	1	01
		10"	54	1 µm	05	7P without fin	A1	Viton	SV				
		20"	55	1.2 µm	10	28 with fin	C0	EPDM	SE				
		30"	56	2.5 µm	06	'O'	D0	FEP Encapsulated Viton	FV				
				5 µm	07								
				8 µm	17								
				10 µm	08								

Example:

CPPA	56	05	A0	SS	1	01
------	----	----	----	----	---	----

For Adapters and Elastomers availability refer Page 91.

Microglassfiber Disc Filters

Microglassfiber Disc Filters Type - GF2

GF2 filters are high dirt holding microglassfiber disc filters specially designed for pre-filtration of solutions with high dirt load.

Pore Size: 1.5 µm

Special Features

- High flow rates
- High dirt holding capacity

Fine Microglassfiber Disc Filters - Type GFS

GFS filters are high retention efficiency fine microglassfiber disc filters.

Pore Size: 2 µm

Special Features

- High retention efficiency

Positively Charged Microglassfiber Disc Filters- Type GFSZ


GFSZ filters are positively charged high retention efficiency microglassfiber disc filters.

Pore Size: 2 µm

Special Features

- Very high retention efficiency for negatively charged particles

Ordering Information



Filters for Polishing and Clarification

mdi offers a range of cartridge filters for polishing and clarification applications requiring absolute removal of particulate contaminants.

These are non media migrating, biologically and chemically inert, large area cartridge filters, offering very high (99.999%) to moderate high (99.9%) retention efficiencies to suit different applications.

Types Available

>> **ClariSure PP**- Pleated Polypropylene cartridge filters

>> **ClariSure PL**- Pleated Profile cartridge filters

>> **ClariSure DP**- Pleated Depth cartridge filters

Applications

- Filtration of organic solvents in non sterile API
- Polishing filtration of Non Sterile API
- Filtration of precipitating agents
- Filtration of wash solvents for final non sterile API
- Filtration of air to dryers and micronizers
- Final wash water for Ampoule/Vial washing and Bung washing
- Bottle washing in oral formulations
- Filtration of feed water for RO plants

Quality Assurance

These filter devices are manufactured in Class 10,000 clean rooms under ISO 9001:2015 certified quality management systems and are validated to meet compendia and regulatory requirements.

Assurance	
Toxicity	Passes Bioreactivity test, In Vivo, as per USP <88> for Class VI plastics
Non Fiber Releasing	Passes test as per USP and comply with USFDA 21 CFR Part 210.3 (b)(6) for fiber release
Extractables with WFI	Passes test as per USP <661>
Oxidizable Substances	Within limits as specified in USP <1231>
Particle Shedding	Passes test as per USP <788> for particulate matter in injections
Indirect Food Additive	All Polypropylene components meet the FDA Indirect Food Additive requirements cited in 21 CFR 177.1520
Good Manufacturing Practice	These products are manufactured in a facility which adheres to Good Manufacturing Practices

Filter Selection Chart

Application Area	Key Application Requirements	Cartridge Filters		
Filtration of organic solvents in non sterile API	<ul style="list-style-type: none"> - High retention efficiency - Wide chemical compatibility 	ClariSure PP Pleated Polypropylene Cartridge Filters	ClariSure PL Pleated Profile Cartridge Filters	ClariSure DP Pleated Depth Cartridge Filters
Polishing filtration of non sterile API	<ul style="list-style-type: none"> - Very high retention efficiency - Wide chemical compatibility 	ClariSure PP Pleated Polypropylene Cartridge Filters	ClariSure PL Pleated Profile Cartridge Filters	-
Filtration of precipitating agents	<ul style="list-style-type: none"> - Very high retention efficiency - Wide chemical compatibility 	ClariSure PP Pleated Polypropylene Cartridge Filters	ClariSure PL Pleated Profile Cartridge Filters	-
Filtration of wash solvents for final non sterile API	<ul style="list-style-type: none"> - Very high retention efficiency - Wide chemical compatibility 	ClariSure PP Pleated Polypropylene Cartridge Filters	ClariSure PL Pleated Profile Cartridge Filters	-
Filtration of air to dryers and micronizers	<ul style="list-style-type: none"> - Very high retention efficiency 	ClariSure PP Pleated Polypropylene Cartridge Filters	ClariSure PL Pleated Profile Cartridge Filters	-
Final wash water for ampoule/vial washing and bung washing	<ul style="list-style-type: none"> - Very high retention efficiency - Wide chemical compatibility 	ClariSure PP Pleated Polypropylene Cartridge Filters	ClariSure PL Pleated Profile Cartridge Filters	-
Bottle washing in oral formulations	<ul style="list-style-type: none"> - Very high retention efficiency - Wide chemical compatibility 	ClariSure PP Pleated Polypropylene Cartridge Filters	ClariSure PL Pleated Profile Cartridge Filters	-
Filtration of feed water for RO plants	<ul style="list-style-type: none"> - High dirt holding capacity 	-	-	ClariSure DP Pleated Depth Cartridge Filters

ClariSure PP Standard Cartridge Filters

**Complies with
USFDA 21 CFR 210.3 (b)(6)**

**Meets and Exceeds
USFDA 21 CFR 177.1520**

ClariSure PP cartridge filters are 100% polypropylene construction precision filtration devices. These are very high retention efficiency (99.999%) pleated cartridge filters offering large filtration area.

These filters are used as terminal filters in applications requiring highly efficient particulate removal.




Specifications

Construction	
Pore Size	0.5 µm, 1 µm, 1.2 µm, 2 µm, 2.5 µm, 5 µm 10 µm, 20 µm, 30 µm
Filter Media	Polypropylene
Support (Drainage) Layers	Polypropylene
Plastic Parts	Polypropylene
Size	
Size	5" 10" 20" 30"
Effective Filtration Area (Nominal)	2500 cm ² 5000 cm ² 10000 cm ² 15000 cm ²
Operational	
Max. Operating Temperature	80 °C @ ≤ 30 psi (2 Kg/cm ²)
Max. Differential Pressure	50 psi (3.5 Kg/cm ²) @ 25 °C

Water Flow Rates

ClariSure PP, 10" Cartridge Filters



Ordering Information

Type		Size		Pore Size		Adapter		Elastomer		Sterility		Pack Size	
	Code		Code		Code		Code		Code		Code		Code
ClariSure PP	CPPP	5"	53	0.5 µm	04	7P	A0	Silicone	SS	Non Sterile	1	1	01
		10"	54	1 µm	05	7P without fin	A1	Viton	SV				
		20"	55	1.2 µm	10	28 with fin	C0	EPDM	SE				
		30"	56	2 µm	15	BEO	B0	PTFE	ST*				
				2.5 µm	06	'O'	D0	FEP Encapsulated Viton	FV				
				5 µm	07								
				10 µm	08								
				20 µm	11								
				30 µm	19								

Example:

CPPP	55	04	A0	FV	1	01
------	----	----	----	----	---	----

*PTFE (ST) gasket seals are available with Adapter Code B0 only
For Adapters and Elastomers availability refer Page 91.

ClariSure PL Standard Cartridge Filters

Complies with
USFDA 21 CFR 210.3 (b)(6)

Meets and Exceeds
USFDA 21 CFR 177.1520

ClariSure PL cartridge filters are all polypropylene high retention efficiency (99.99%) pleated cartridge filters offering large filtration area.



Specifications

Construction

Pore Size 1 µm, 2.5 µm, 5 µm, 10 µm

Filter Media Polypropylene

Support (Drainage) Layers Polypropylene

Plastic Parts Polypropylene

Size

Size 5" 10" 20" 30"

Effective Filtration Area (Nominal) 2500 cm² 5000 cm² 10000 cm² 15000 cm²


Operational

Max. Operating Temperature 80 °C @ ≤ 30 psi (2 Kg/cm²)

Max. Differential Pressure 50 psi (3.5 Kg/cm²) @ 25 °C

Water Flow Rates

ClariSure PL, 10" Cartridge Filters



Ordering Information

Type		Size		Pore Size		Adapter		Elastomer		Sterility		Pack Size	
	Code		Code		Code		Code		Code		Code		Code
ClariSure PL	CPPL	5"	53	1 µm	05	7P	A0	Silicone	SS	Non Sterile	1	1	01
		10"	54	2.5 µm	06	7P without fin	A1	Viton	SV				
		20"	55	5 µm	07	28 with fin	C0	EPDM	SE				
		30"	56	10 µm	08	BEO	B0	PTFE	ST*				
						'O'	D0	FEP Encapsulated Viton	FV				

Example:

CPPL	54	07	B0	SV	1	01
------	----	----	----	----	---	----

*PTFE (ST) gasket seals are available with Adapter Code B0 only

For Adapters and Elastomers availability refer Page 91.

ClariSure DP Standard Cartridge Filters

Complies with
USFDA 21 CFR 210.3 (b)(6)

Meets and Exceeds
USFDA 21 CFR 177.1520

ClariSure DP cartridge filters are specially designed medium retention efficiency, multilayered, all polypropylene depth cartridge filters offering high dirt holding capacity as well as high flow rates due to their pleated configuration.

These filters are primarily used as pre-filters but can also be used as terminal filters for non-critical applications.



Specifications

Construction				
Pore Size	1 µm, 1.5 µm, 3 µm, 5 µm, 10 µm, 20 µm 30 µm, 40 µm			
Filter Media	Polypropylene			
Support (Drainage) Layers	Polypropylene			
Plastic Parts	Polypropylene			
Size				
Size	5"	10"	20"	30"
Effective Filtration Area (Nominal)	2000 cm ²	4000 cm ²	8000 cm ²	12000 cm ²
Operational				
Max. Operating Temperature	80 °C @ ≤ 30 psi (2 Kg/cm ²)			
Max. Differential Pressure	50 psi (3.5 Kg/cm ²) @ 25 °C			

Ordering Information

Type		Size		Pore Size		Adapter		Elastomer		Sterility		Pack Size	
	Code		Code		Code		Code		Code		Code		Code
ClariSure DP	CPDP	5"	53	1 µm	05	7P	A0	Silicone	SS	Non Sterile	1	1	01
		10"	54	1.5 µm	14	7P without fin	A1	Viton	SV				
		20"	55	3 µm	16	28 with fin	C0	EPDM	SE				
		30"	56	5 µm	07	BEO	B0	PTFE	ST*				
				10 µm	08	'O'	D0	FEP Encapsulated Viton	FV				
				20 µm	11								
				30 µm	19								
				40 µm	20								

Example:

CPDP	56	14	A0	SS	1	01
------	----	----	----	----	---	----

*PTFE (ST) gasket seals are available with Adapter Code B0 only
For Adapters and Elastomers availability refer Page 91.

Chemical Compatibility

Table below shows the chemical compatibility of various process filtration products with some commonly used solvents. All products were exposed to specified chemicals for 72 hours at 25 °C. Chemical compatibility data with specific reagents is available on request.

Reagents	AseptiSure Cartridge Filters					ClariSure Cartridge Filters								'O' Rings/Gasket Seals			
	HS/KS	HSR	WS	NS	TH/TF	GK	GS	GP	PA	PP	PL	DP	Silicone	Viton	EP	FEP Encapsulated Viton	
Solvents																	
Acetone	N	N	G	G	G	N	G	G	G	G	G	G	N	N	G	G	
Acetonitrile	G	G	G	G	G	G	G	G	G	G	G	G	G	N	G	G	
Benzene	G	G	N	G	G	F	F	F	F	F	F	F	N	G	N	G	
Benzyl Alcohol	N	N	N	G	G	N	G	G	G	G	G	G	G	G	G	G	
Benzyl Alcohol 2%	G	G	G	G	G	G	G	G	G	G	G	G	G	G	G	G	
Diethyl Ether	G	G	G	G	G	N	N	N	N	N	N	N	N	N	N	N	G
Dimethylformamide	N	N	G	G	G	N	G	G	G	G	G	G	G	N	N	G	
Ethanol	G	G	G	G	G	G	G	G	G	G	G	G	G	G	G	G	
Ethyl Acetate	G	G	G	G	G	G	G	G	G	G	G	G	N	N	G	G	
Ethylene Glycol	G	G	G	G	G	G	G	G	G	G	G	G	G	G	G	G	
Hexane	F	F	N	F	F	F	F	F	F	F	F	F	N	G	N	G	
Iso Propyl Alcohol	G	G	G	G	G	G	G	G	G	G	G	G	G	G	G	G	
Methanol	G	G	G	G	G	G	G	G	G	G	G	G	G	N	G	G	
Methylene Chloride	N	N	N	N	N	N	N	N	N	N	N	N	N	N	N	N	G
n-Butanol	G	G	G	G	G	G	G	G	G	G	G	G	N	G	G	G	
Peanut oil	G	G	G	G	G	G	G	G	G	G	G	G	G	G	G	G	
Tetrahydrofuran/Water (50:50)	N	N	N	N	G	N	N	N	N	N	N	N	N	N	N	N	G
Toluene	G	G	N	G	G	G	G	G	G	G	G	G	N	G	N	G	
Trichloroethylene	N	N	N	N	N	N	N	N	N	N	N	N	N	G	N	G	
Acids																	
Hydrochloric Acid 37%	G	G	G	N	G	G	G	G	G	G	G	G	N	G	N	G	
Hydrofluoric Acid 10%	G	G	NT	F	G	N	N	N	G	G	G	G	N	G	F	G	
Nitric Acid 67%	N	N	N	N	G	N	G	G	G	G	G	G	N	G	N	G	
Nitric Acid 7%	G	G	G	N	G	G	G	G	G	G	G	G	G	G	G	G	
Sulphuric Acid 10%	G	G	G	F	G	G	G	G	G	G	G	G	G	G	G	G	
Bases																	
Ammonium Hydroxide 25%	N	G	N	N	G	N	N	G	G	G	G	G	G	G	G	G	
Sodium Hydroxide 32%	N	G	N	N	G	N	N	G	G	G	G	G	G	G	G	G	
Potassium Hydroxide 32%	N	G	N	N	G	N	N	G	G	G	G	G	G	G	G	G	

G: Good F: Fair N: Not recommended NT : Not Tested

Chemical Compatibility

Table below shows the chemical compatibility of various process filtration products with some commonly used solvents. All products were exposed to specified chemicals for 72 hours at 25 °C. Chemical compatibility data with specific reagents is available on request.

Reagents	AseptiCap Capsule Filters				AseptiVent Capsule Filters	ClariCap Capsule Filters			Membrane Disc Filters	Microglassfiber Pre-filter Discs
	KL/ KS	KO/ KSO	NL/ NS	WS	TF	GK	GS	PP	NN	GF
Solvents										
Acetone	N	N	G	G	G	N	G	G	G	G
Acetonitrile	G	G	G	G	G	G	G	G	G	G
Benzene	G	G	G	N	G	F	F	F	G	G
Benzyl Alcohol	N	N	G	N	G	N	G	G	G	G
Benzyl Alcohol 2%	G	G	G	G	G	G	G	G	G	G
Diethyl Ether	G	G	G	G	G	N	N	N	G	G
Dimethylformamide	N	N	G	G	G	N	G	G	G	G
Ethanol	G	G	G	G	G	G	G	G	G	G
Ethyl Acetate	G	G	G	G	G	G	G	G	G	G
Ethylene Glycol	G	G	G	G	G	G	G	G	G	G
Hexane	N	N	N	N	N	N	N	N	G	G
Iso Propyl Alcohol	G	G	G	G	G	G	G	G	G	G
Methanol	G	G	G	G	G	G	G	G	G	G
Methylene Chloride	N	N	N	N	N	N	N	N	G	G
n-Butanol	G	G	G	G	G	G	G	G	G	G
Peanut oil	G	G	G	G	G	G	G	G	G	G
Tetrahydrofuran/Water (50:50)	N	N	N	N	G	N	N	N	G	G
Toluene	N	N	N	N	N	N	N	N	G	G
Trichloroethylene	N	N	N	N	N	N	N	N		G
Acids									N	
Hydrochloric Acid 37%	G	G	N	G	G	G	G	G	G	G
Hydrofluoric Acid 10%	G	G	G	NT	G	N	N	G	N	N
Nitric Acid 67%	N	N	N	N	G	N	G	G	G	G
Nitric Acid 7%	G	G	G	G	G	G	G	G	G	G
Sulphuric Acid 10%	G	G	G	G	G	G	G	G		G
Bases									G	
Ammonium Hydroxide 25%	N	G	N	N	G	N	G	G	G	G
Sodium Hydroxide 32%	N	G	N	N	G	N	G	G	N	F
Potassium Hydroxide 32%	N	G	N	N	G	N	G	G	N	G

G: Good F: Fair N: Not recommended NT: Not Tested

End Connection Availability Chart for Capsule Filters

Connections Available			
Inlet/ Outlet	25mm	37mm	50mm
¼" SHB I/O	X	√	√
¾" Sanitary Flange I/O	X	X	Outlet Only
Female Luer Lock	Inlet Only	X	X
Male Luer Slip	Outlet Only	X	X
1/8" Hose Barb I/O	√	X	X

End Connections	Small Capsule Filters			
	1"	2"	5"	8"
½" Hose Barb	√	√	√	√
Single Step ½" Hose Barb	X	√	√	√
¼" Stepped Hose Barb	√	√	√	√
1½" Sanitary Flange	√	√	√	√
¾" Sanitary Flange	√	√	√	√
½" MNPT	X	√	√	√
¼" MNPT	√	√	√	√
Quick Connector	√	√	√	√
Female Luer Lock	√	√	√	√
Male Luer Slip	Outlet Only	X	X	X
¾" Hose Barb	√	√	Outlet Only	X
⅜" Hose Barb	X	√	√	√

Bell is available with ¼" SHB outlet in 1" Capsule Filters only

Bell is available with ½" Hose Barb outlet in 1", 2", 5" and 8" Capsule Filters

End Connections	Large Capsule Filters						
	Inline				T-line		
	5"	10"	20"	30"	10"	20"	30"
Single Step ½" Hose Barb	√	√	√	√	X	X	X
1½" Sanitary Flange	√	√	√	√	√	√	√
¾" Sanitary Flange	√	√	X	X	X	X	X
⅜" Hose Barb	√	√	√	√	X	X	X
1" Hose Barb	X	√	√	√	X	X	X

Dimensions: Capsule Filters

Inlet/ Outlet	Inline Capsule Filters		
	25mm	37mm	50mm
¼" - ⅜" Stepped Hose Barb I/O	-	64 mm	79 mm
¼" Single Step Hose Barb I/O	38 mm	-	-
¾" Sanitary Flange I/O	-	-	51 mm
Female Luer Lock Inlet/ Male Luer Slip Outlet	23 mm	-	-
1/8" Hose Barb I/O	36 mm	-	-
Operational Radius	15 mm	23 mm	28 mm

End Connections	Small Capsule Filters			
	1"	2"	5"	8"
¼" SHB I/O	94 mm	122 mm	172 mm	223 mm
¾" Sanitary Flange Inlet I/O	85 mm	104 mm	155 mm	206 mm
Quick Connector	100 mm	113 mm	164 mm	218 mm
1½" Sanitary Flange I/O	92 mm	112 mm	164 mm	216 mm
½" Hose Barb I/O	90 mm	112 mm	162 mm	214 mm
½" Single Step Hose Barb I/O	-	115 mm	165 mm	218 mm
1½" Sanitary Flange Inlet ½" Single Step Hose Barb Outlet	-	112 mm	165 mm	216 mm
⅜" Hose Barb I/O	-	115 mm	167 mm	217 mm
Operational Radius	40 mm	65 mm	65 mm	65 mm

End Connections	Inline Capsule Filters				T-line Capsule Filters		
	5"	10"	20"	30"	10"	20"	30"
1½" Sanitary Flange I/O	205 mm	330 mm	600 mm	855 mm	340 mm	580 mm	840 mm
¾" Sanitary Flange I/O	214 mm	335 mm	x	x	x	x	x
½" Single Step Hose Barb I/O	218 mm	336 mm	630 mm	890 mm	x	x	x
1½" Sanitary Flange Inlet ½" Hose Barb Outlet	212 mm	334 mm	620 mm	870 mm	x	x	x
⅜" Hose Barb I/O	211 mm	332 mm	634 mm	885 mm	x	x	x
1" Hose Barb I/O	x	405 mm	635 mm	895 mm	x	x	x
Operational Radius	80 mm	80 mm	80 mm	80 mm	80 mm	80 mm	80 mm

Adapter and Elastomers Availability Chart for Cartridge Filters

Mini Cartridge Filters		
Adapters	2.5"	5"
4463	√	√
4463B	√	√
4440	√	√
Seal-K	√	√
Seal-O	X	√
Seal-M	√	√

Mini Cartridge Filters	
Adapters	Elastomer
	Silicone
4463	√
4463B	√
4440	√
Seal-K	X
Seal-O	√
Seal-M	√

Standard Cartridge Filters				
Adapters	5"	10"	20"	30"
7P	√	√	√	√
7P without Fin	√	√	√	√
28 with Fin	X	√	√	√
'O'	X	√	√	√

Standard Cartridge Filters				
Adapters	Elastomers			
	Silicone	Viton	EPDM	FEP Encapsulated Viton
7P	√	√	√	√
7P without Fin	√	√	√	√
28 with Fin	√	√	√	X
'O'	√	√	√	X

Ordering Information

Shipment details for customers outside India

Through Federal Express, UPS, or DHL courier (specify complete street address).

By air freight for large quantities (specify airport of discharge).

Goods usually reach destination within 5-10 days from date of shipment.

Membrane products are light weight and air freight charges usually vary between 3% to 10% of the value. Any duties/taxes in the country of destination are the responsibility of the consignee.

Shipment details for customers in India

The consignments can be sent through courier. Courier charges will be borne by the customer. Please specify the preferred courier and provide any form and instructions for octroi etc. that may be required for shipment.

How to order

Orders may be placed by phone/fax/email/mail directly to Sales.

Advanced Microdevices Pvt. Ltd.

20-21, Industrial Area, Ambala Cantt - 133 006, INDIA

Tel: +91-171-2699290, 2699471

Email: info@mdimembrane.com

Quality Policy

Quality is built into **mdi** products and services by not only adhering to well designed quality systems to consistently produce high quality, internationally acceptable products but also by striving to incorporate superior performance parameters into all our products and services and provide our customers with a unique performance advantage in their application. Our quality policy provides a glimpse of our commitment:

mdi strives to provide to its customers products and services of highest standards possible, consistently superior, and more satisfying than what is available anywhere else."

Stride Towards Excellence

At **mdi**, our mission is to constantly strive to achieve excellence in all our endeavors by establishing systems to create excellent products and services to fulfil the needs of our customers. To achieve this we

- ◆ Frequently compare our products with competing brands
- ◆ Simulate tests for functional use
- ◆ Develop easy-to-use innovative products

We are constantly working on improvements and welcome suggestions from our customers.

Guarantee

All **mdi** products are guaranteed and are backed by our

- ◆ Technical expertise and experience of over 45 years
- ◆ Validated **mdi** process' for consistency and repeatability
- ◆ Strict quality control and quality assurance regimen
- ◆ Certificate of Analysis accompanying all shipments

We have an unconditional replacement policy in case of any defects.





Worldwide Exports

INDIA Branch Offices

AHMEDABAD

Mobile : +91-9328257987

BANGALORE

Mobile : +91-9972587761

CHENNAI

Mobile : +91-9972587761

HYDERABAD

Mobile : +91-9391935423

MUMBAI

Tel : +91-22-40214436

Fax : +91-22-40214435

Mobile : +91-9820026178

USA Office

Tel : +1-717-412-0943

Fax : +1-717-695-9637

Other Literature Available

mdi Laboratory Filtration Product Guide

mdi Biotech Product Guide

mdi Diagnostic Product Guide

Advanced Microdevices Pvt. Ltd.

20-21, Industrial Area, Ambala Cantt -133006, INDIA

Tel : +91-171-2699290, 2699471

E-mail : info@mdimembrane.com

Website : www.mdimembrane.com

mdi Membrane Technologies INC.

75 Utlely Drive STE 103, Camp Hill, PA 17011

Phone: +1-717-412-0943, Fax: +1-717-695-9637

E-mail: rs@mdimembranetech.com

Website: www.mdimembranetech.com