

**Reference :** CHANGE NOTICE No. MB/09/06

**Subject :** Change in mdi Closed Sterility Test Device Stericheck-C and Stericheck-N (All types)

**Scope :** This change notification will affect the following parts/components of mdi Stericheck:

- Drain Plug - Improved grip

**Background:**

In our efforts to provide the user better products with more convenience and better performance, certain changes have been made in the parts/components of **mdi Stericheck-C and Stericheck-N**.

**Need for Change:**

As above.

**How Does It Affect the User:**

From a practical point of view the changes done will improve the product performance characteristics and ease-of-use under actual conditions of use. All performance specifications and product attributes remain unchanged.

- **Drain Plug:** In addition to the change in design, the size of wings has also been increased for better grip and fitment.



**Effective Date:**

**mdi** Closed Sterility Test Device **Stericheck-C** and **Stericheck-N** (all types) manufactured after November 12, 2009 will have said changes.

In case you have any queries, please feel free to contact our Technical Support Team at 'mdi@vsnl.com'.

**Quality Assurance**  
**Advanced Microdevices Pvt. Ltd.**