

Reference : CHANGE NOTICE No. MB/09/07

Subject : Change in mdi Closed Sterility Test Device Stericheck (All types)

Scope : This change notification will affect the following parts/components of mdi Stericheck:

- Needles - Improved thermal performance with better grip

Background:

In our efforts to provide the user better products with more convenience and better performance, certain changes have been made in the parts/components of **mdi Stericheck**.


Need for Change:

As above.

How Does It Affect the User:

From a practical point of view the changes done will improve the product performance characteristics and ease-of-use under actual conditions of use. All performance specifications and product attributes remain unchanged.

- **Needle:** The plastic polymer used now can withstand higher temperature and the needle design was improved for better grip.



Effective Date:

mdi Closed Sterility Test Device **Stericheck** (all types) manufactured after November 18, 2009 will have said changes.

In case you have any queries, please feel free to contact our Technical Support Team at 'mdi@vsnl.com'.



Quality Assurance

Advanced Microdevices Pvt. Ltd.