

Reference : CHANGE NOTICE No.: CN/15/03

Subject : Improved design of Protective Cap for end connection type Single Step 1/2" Hose barb (Code - Q) used with mdi Capsule Filters

Scope : This change notification will affect the following products supplied with end connection type Single Step 1/2" Hose barb (Code - Q):
mdi Gamma Sterilizable and Gamma Irradiated Capsule Filters with sizes 2", 5", 8", 10", 20" and 30" (All types)

Background:

In our efforts to provide the user better products, certain changes have been made in the design of Protective Cap for end connection type Single Step 1/2" Hose barb (Code - Q) used with mdi Capsule Filters as mentioned above.

Need for Change:

As above.

How Does It Affect the User:

From a practical point of view the changes done will improve the robustness of product during transportation and under actual use. All performance specifications and product attributes remain unchanged.

- Protective caps for Single step 1/2" Hose barb end connection are now inserted in the end connection, which were earlier assembled over the end connection. This will reduce the loosening/cracking of the protective caps during transportation and under actual use e.g. Gamma irradiation etc.



Effective Date:

mdi Capsule filters (types and sizes as mentioned above) manufactured after August 06, 2015 will have said changes.

In case you have any queries, please feel free to contact our Technical Support Team at 'info@mdimembrane.com'.



Quality Assurance
Advanced Microdevices Pvt. Ltd.