

Change Notice No.	Notification date	Implementation date
CN/A/22/002	June 15, 2022	July 15, 2022

Subject: Change in the Design of bottle cap adapters used in mdi Single Use Assemblies/Systems, and Single Use Components.

Scope: This change notification will affect the products with following catalog numbers: (X can be any alphabet or numeral):

- AXXXXXXXXXXXXXXXXX
- APXXXXXXXXXXXXXXXXX
- VBXXXXXXXXXXXXXXXXX

Background:

The design of Bottle cap adapter connector used in single use assembly/ system has been improved to make it more robust.

Need for Change:

There is a possibility that Bottle cap adapters may get damaged upon application of excessive force during use. In order to provide more robust bottle cap adapters, the wall thickness of the bottle cap adapters has been increased.


How Does It Affect the User:

From the user point of view, the changes done shall provide more robust bottle cap adapters. There is no impact on fitment of bottle cap adapters on the bottle. It shall have no impact on the performance specifications, product attributes, material of construction and traceability of the product.

Effective Date:

mdi Bottle cap adapters manufactured after the implementation date (**July 15, 2022**) will have the said changes.

In case you have any queries, please feel free to contact our Technical Support Team at 'info@mdimembrane.com'.



**Head - Quality Assurance (Single Use Systems)
Advanced Microdevices Pvt. Ltd.**