

Change Notice No.	Notification date	Implementation date
CN/A/23/004	September 18, 2023	December 18, 2023

**Subject:** Change in height of mdi 8" TC port used in mdi Single Use Assemblies.


**Scope:** This change will affect the Single Use Assemblies having mdi 8" TC port as a component.

The affected drawings are mentioned in annexure.


**Background:**

The existing design of 8" TC port in 3D bag is not user friendly, as the space between port and bag is not enough for free movement of the screw on TC clamp, during the clamping process.

In order to make the clamping process convenient and user friendly, the height of the TC port has been increased.



**Old Design**



**New Design**

**Need for Change:**

The change will make the clamping process of 8" TC port convenient and user friendly.

**How Does It Affect the User:**

From a user point of view, the change in height of the port will make it use friendly and it will not affect the other properties e.g fitment, dimension of the product under actual conditions of use. The change also does not affect the quality of the product. All performance specifications, product attributes and material of construction will remains the same.

Change Notice No.	Notification date	Implementation date
CN/A/23/004	September 18, 2023	December 18, 2023

**Effective Date:**

mdi Single Use Assemblies and Systems, released after the implementation date (December 18, 2023) will have the said changes.

In case you have any queries, please feel free to contact our Technical Support Team at 'info@mdimembrane.com'.



**Head - Quality Assurance (Single Use Systems)  
Advanced Microdevices Pvt. Ltd.**

**Annexure for affected drawing**

A2696XXXXXXXXC301	A3242XXXXXXXXC301	A4061XXXXXXXXC301
A550XXXXXXXXXC301	A3358XXXXXXXXC301	A4060XXXXXXXXC301
A3813XXXXXXXXC301	A3322XXXXXXXXC301	A4089XXXXXXXXC301
A2360XXXXXXXXC301	A3323XXXXXXXXC301	A3814XXXXXXXXC301
A2159XXXXXXXXC301	A3324XXXXXXXXC301	A4149XXXXXXXXC301
A2160XXXXXXXXC301	A3325XXXXXXXXC301	A4088XXXXXXXXC301
A2471XXXXXXXXC301	A3326XXXXXXXXC301	A4577XXXXXXXXC301
A2390XXXXXXXXC301	A3673XXXXXXXXC301	A4574XXXXXXXXC301
A2622XXXXXXXXC301	A3675XXXXXXXXC301	A4572XXXXXXXXC301
A2663XXXXXXXXC301	A3674XXXXXXXXC301	A4573XXXXXXXXC301
A2662XXXXXXXXC301	A3979XXXXXXXXC301	A4641XXXXXXXXC301
A2633XXXXXXXXC301	A3624XXXXXXXXC301	A3958XXXXXXXXC301
A2694XXXXXXXXC301	A3663XXXXXXXXC301	A4392XXXXXXXXC301
A2869XXXXXXXXC301	A3878XXXXXXXXC301	A4391XXXXXXXXC301
A2473XXXXXXXXC301	A3888XXXXXXXXC301	A4390XXXXXXXXC301
A2695XXXXXXXXC301	A3879XXXXXXXXC301	A4844XXXXXXXXC301
A2474XXXXXXXXC301	A4062XXXXXXXXC301	A5104XXXXXXXXC101
A3241XXXXXXXXC301	A4063XXXXXXXXC301	A4993XXXXXXXXC301
A5060XXXXXXXXC301	A5059XXXXXXXXC301	-