

Change Notice No.	Notification date	Implementation date
CN/A/25/003	March 03, 2025	April 03, 2025

**Subject : Change in position of 4" TC port used in AseptiMix MI Single Use Mixer Bags.**

**Scope : This change will affect AseptiMix MI Mixer Bag with 4" TC port. The affected drawing numbers are mentioned in Annexure 1.**

### **Background**

During packaging of AseptiMix MI Single Use Mixer Bags with 4" TC port, the mixer bags are packed in such a way that the blades of impeller are vertical to the port and directly come in contact to 4" TC port.

### **Need for Change**

Due to direct contact of 4" TC port with the impeller, it is sometimes observed that the impeller get damage when an excessive force is applied on the 4" TC port against the impeller. To mitigate the damage of impeller in AseptiMix MI Single Use Mixer Bags with 4" TC port, the position of port shall be shifted in such a way that it will not directly come in contact to the impeller.

### **How Does It Affect the User**

From the user point of view, only the location of TC port will be shifted and it does not impact the fitment of mixer bag. All product specifications including material of construction as well as other performance specifications shall remain unchanged.

### **Implementation of Change**

**mdi** AseptiMix MI Single Use Mixer Bag with 4" TC port manufacture after the implementation date will have the said changes in product.

In case you have any queries, please feel free to contact our Technical Support Team at 'info@mdimembrane.com'.



**Head - Quality Assurance Single Use Assembly  
Advanced Microdevices Pvt. Ltd.**

<b>Change Notice No.</b>	<b>Notification date</b>	<b>Implementation date</b>
<b>CN/A/25/003</b>	<b>March 03, 2025</b>	<b>April 03, 2025</b>

**Annexure 1**

**Affected Drawings**

A2463XXXXXXXXC301	A4711XXXXXXXXC301	A6458XXXXXXXXC301
A2474XXXXXXXXC301	A4759XXXXXXXXC301	A6830XXXXXXXXC301
A2710XXXXXXXXC301	A4780XXXXXXXXC301	A7004XXXXXXXXC301
A2804XXXXXXXXC301	A4782XXXXXXXXC301	A7039XXXXXXXXC301
A2965XXXXXXXXC301	A4783XXXXXXXXC301	A7197XXXXXXXXC301
A3301XXXXXXXXC301	A4784XXXXXXXXC301	A7198XXXXXXXXC301
A3917XXXXXXXXC301	A4785XXXXXXXXC301	A4393XXXXXXXXC301
A4028XXXXXXXXC301	A4805XXXXXXXXC301	A4394XXXXXXXXC301
A4029XXXXXXXXC301	A5574XXXXXXXXC301	A4395XXXXXXXXC301
A4030XXXXXXXXC301	A5690XXXXXXXXC301	A5923XXXXXXXXC301
A4031XXXXXXXXC301	A5703XXXXXXXXC301	A6086XXXXXXXXC301
A4230XXXXXXXXC301	A5909XXXXXXXXC301	-