

Change Notice No.	Notification date	Implementation date
CN/A/25/008	May 17, 2025	May 26, 2025

Subject : Change in design of mdi AseptiFit Straight Connector used as component in Single Use Assemblies.

**Scope : This change will affect the AseptiFit Straight Connector with catalog number mentioned below:
FPXA04QQXXXXXXXX**

Background:

In our efforts to improve the product design and flow rate, a change has been initiated for AseptiFit straight connector used as a component in Single Use Assemblies. Following observations were made for existing connector:

- As per customer's feedback the flow rate of straight connector with size 2 x 1/2" is not appropriate.
- Sometimes, whiteness observed in straight connector when an excessive force applied on it.

Need for Change:

In order to improve the flow rate and to strengthen the connector, it shall be improved with addition of ring design same as other AseptiFit Fittings.

No changes have been made/done to the material of construction, manufacturing process and traceability.

How Does It Affect the User:

This change is related to design of AseptiFit straight connector and it will impact the flow rate of the product. However, material of construction and performance specifications other than the flow rate shall remain unchanged. The pictures describing the change are shown below.



Existing



New

Change Notice No.	Notification date	Implementation date
CN/A/25/008	May 17, 2025	May 26, 2025

Implementation of Change:

mdi AseptiFit Straight Connector with part number as mentioned in the scope, manufactured after the implementation date, will be with improved design.

In case you have any queries, please feel free to contact our Technical Support Team at 'info@mdimembrane.com'.



Head - Quality Assurance Single Use Systems
Advanced Microdevices Pvt. Ltd.