

Change Notice No.	Notification date	Implementation date
CN/F/19/016	November 18, 2019	November 25, 2019

Subject : Improved design of vent/drain valves of mdi Large Capsule Filters

Scope : This change notification will affect the products with following catalog numbers: (# can be any alphabet or numeral):

- LB#####
- LD#####
- LG#####
- LK#####
- LN#####
- LO#####
- LP#####
- LT#####
- LV#####
- LW#####

Background:

In our efforts to provide the user, better products which are more robust and able to withstand rough handling, certain changes have been made in vent/drain valves of **mdi** Large Capsule Filters with catalog numbers as mentioned above.


Need for Change:

It was observed that existing vent/drain valves sometimes used to get damaged/broken due to rough handling during usage or transit. To overcome this problem/issue, the vent/drain valves of **mdi** Large Capsule Filters have been made more robust through improvement in design. The new design of vent/drain valves has been verified for critical product performance parameters including functionality (leakproof fitment) and ability to withstand rough handling. No change has been done in any other critical dimensions or materials of construction or manufacturing process of **mdi** Large Capsule Filters.

Change Notice No.	Notification date	Implementation date
CN/F/19/016	November 18, 2019	November 25, 2019

How Does It Affect the User:

From a practical point of view the changes done in **mdi** Large Capsule Filters will improve the ability of vent/drain valves to withstand rough handling. All performance specifications and product attributes will remain unchanged, except the improved design of vent/drain valves.



Implementation of Change:

mdi Large Capsule Filters having catalog numbers as mentioned above, manufactured after the implementation date will have the said changes. However, the available stock of **mdi** Large Capsule Filters with catalog numbers as mentioned above having vent/drain valves with old design will be received by you till the stock lasts.

In case you have any queries, please feel free to contact our Technical Support Team at 'info@mdimembrane.com'.

Head - Quality Assurance
Advanced Microdevices Pvt. Ltd.