

Change Notice No.	Notification date	Implementation date
CN/F/20/005	January 01, 2021	February 01, 2021

Subject : Rationalization of the design of capsule housing used in mdi Large Capsule Filter

Scope : This change will affect the products with following catalog numbers: (# can be any alphabet or numeral):

- LB##54#####T###
- LB##55#####
- LB##56#####
- LG##54#####T###
- LG##55#####
- LG##56#####
- LK##54#####T###
- LK##55#####
- LK##56#####
- LN##55#####
- LO##54#####T###
- LO##55#####
- LO##56#####
- LT##54#####T###
- LT##55#####
- LT##56#####
- LV##55#####
- LW##55#####

Background:

In our efforts to rationalize our products, certain changes have been made in the design of capsule housing used in **mdi** Large Capsule Filters with catalog numbers as mentioned above.


Need for Change:

The capsule housing used in **mdi** Large Capsule Filters (except T-line) of sizes 5" and 10" are having multiple ribs on the outer surface. But, the capsule housing used in **mdi** T-Line Large Capsule Filters and capsule filters having size greater than 10" are without ribs on the outer surface. To rationalize **mdi** Large Capsule Filters aesthetically, the design of capsule housing used in **mdi** T-Line Large Capsule Filters and capsule filters having size greater than 10" has been improved by providing ribs on the outer surface, similar to the capsule housing used in Large Capsule Filters of sizes 5" and 10". No change whatsoever has been done in inner cartridge element, critical for the performance of Large capsule filters. The length of respective Large capsule filters has also been re-defined marginally. The new design of capsule housing used in **mdi** T-Line Large Capsule Filters and capsule filters having size greater than 10" has been verified for all critical product performance parameters including safety related attributes. No change has been done in product attributes or product specifications or materials of construction.

Change Notice No.	Notification date	Implementation date
CN/F/20/005	January 01, 2021	February 01, 2021

How Does It Affect the User:

From a practical point of view the changes done will improve the appearance of **mdi** T-Line Large Capsule Filters and capsule filters having size greater than 10" by providing ribs on the capsule housing and re-defined length. All performance specifications and product attributes will remain unchanged.



Implementation of Change:

mdi T-Line Large Capsule Filters and capsule filters having size greater than 10" with catalog numbers as mentioned above, manufactured after the implementation date will have the said changes. However, the available stock of the same having old design of capsule housing will be supplied till the stocks last.

In case you have any queries, please feel free to contact our Technical Support Team at 'info@mdimembrane.com'.

**Head - Quality Assurance
Advanced Microdevices Pvt. Ltd.**