

Change Notice No.	Notification date	Implementation date
CN/F/22/006	October 13, 2022	October 20, 2022

Subject : Change in appearance of **mdi** Quick Connectors for better identification

Scope : This change notification will affect the products with following catalog numbers: (# can be any alphabet or numeral):

- QCCX#####X##1##
- QCCX#####X##2##
- QCMX#####X##1##
- QCMX#####X##2##

Background:

In our efforts to provide the users, product with better identification, certain changes have been done in the appearance of **mdi** quick connectors with catalog numbers as mentioned above. The change will have no impact whatsoever on the product performance or it's regulatory compliance.

Need for Change:

mdi Quick connectors are available in two types, autoclavable and gamma irradiatable. Currently, all types of **mdi** Quick Connectors are manufactured as transparent blue in colour. To improve the product identification, it has been decided to assign different colours to different types of **mdi** Quick Connectors, without changing the basic material of construction.

The colour of autoclavable Quick connectors has been changed from transparent Blue to transparent. However, the colour of gamma irradiatable Quick connectors will remain the same (i.e. Transparent blue).

No change has been done in any performance parameters or any critical dimensions or manufacturing process.

How Does It Affect the User:

From a user point of view, the change done will only improve the identification of **mdi** Quick connectors with catalog numbers as mentioned above, through change in appearance. All performance specifications, product attributes, traceability and basic material of construction will remain unchanged.

Change Notice No.	Notification date	Implementation date
CN/F/22/006	October 13, 2022	October 20, 2022

Implementation of Change:

mdi Quick connectors having catalog numbers as mentioned in scope, manufactured after the implementation date will have the said changes. However, the available stock of **mdi** Quick connectors with catalog numbers as mentioned above and transparent blue in colour will be received by you till the stocks last.

In case you have any queries, please feel free to contact our Technical Support Team at 'info@mdimembrane.com'.



Head - Quality Assurance
Advanced Microdevices Pvt. Ltd.