

Reference : CHANGE NOTICE No. CP/12/02

Subject : Modification in identification on 10" Large Capsule Filter (All types)

Scope : This change notification will affect the following products:

- mdi Large capsule filter type LKL
- mdi Large capsule filter type LKS
- mdi Large capsule filter type LGS
- mdi Large capsule filter type LNN
- mdi Large capsule filter type LTL
- mdi Large capsule filter type LGK

Background:

In our efforts to provide the user better products with better identification and traceability, certain changes have been made in the identification on mdi large capsule filter (all types) of 10" length.

Need for Change:

To improve the identification and traceability of the large capsule filters under different conditions of use.

How Does It Affect the User:

From a user point of view the changes done will improve the identification, traceability and ease-of-use under actual conditions of use. All performance specifications, product attributes and main content of the identification remain the same.

- **Identification on large capsule filter** – Laser etching method is now used to print the identification details on large capsule filter which is done using a thermal label earlier. This reduces detaching of thermal label during sterilization process.

OLD (Thermal Label)



NEW (Laser Etching)



Effective Date:

mdi Large Capsule Filters (all types) of 10" length manufactured after December 05, 2012 will be with said changes.

In case you have any queries, please feel free to contact our Technical Support Team at 'mdi@vsnl.com'.



Quality Assurance
Advanced Microdevices Pvt. Ltd.