

Reference : CHANGE NOTICE No. CR/11/01

Subject : Change in Water Intrusion Rate specification for 5" 0.2 µm PTFE Membrane Cartridge Filters type CPTF and CPTH

Scope : This change notification will affect the following products:

- **mdi** 0.2 µm PTFE membrane cartridge filter types CPTF and CPTH of length 5" (All adaptor types except with adaptor code 4463 and 4463B).

Background:

In our efforts to rationalize the product test specifications, the water intrusion specification for 5" cartridge filters types CPTF and CPTH (all adaptor types except adaptor code 4463 and 4463B) has been rationalized/changed.

Need for Change:

As above.

How Does It Affect the User:

From a practical point of view this change will only rationalize the water intrusion rate specification for 5" cartridge filters types mentioned above (all adaptor types except adaptor code 4463 and 4463B). The change will be reflected in Certificate of Quality for the product mentioned above. All performance specifications and product attributes including integrity of cartridge filters remain unchanged.

- The water intrusion rate tested at pressure 2.0 Kg/cm² for 5" cartridge filters types mentioned above with all adaptor types except adaptor code 4463 and 4463B will hence forth be ≤ 7.5 ml/10 min instead of ≤ 6.0 ml/10 min earlier.
- The water intrusion rate tested at pressure 2.0 Kg/cm² for 5" cartridge filters types mentioned above with adaptor code 4463 and 4463B remains unchanged i.e. ≤ 6.0 ml/10 min.

Effective Date:

mdi 0.2 µm PTFE cartridge filters of length 5" (adaptor types mentioned above) shipped/ordered after February 28, 2011 will be with said changes.

In case you have any queries, please feel free to contact our Technical Support Team at 'mdi@vsnl.com'.

