

Reference: CHANGE NOTICE No. CR/13/05

Subject : Modification in identification on Cartridge Filter

Scope : This change notification will affect the following products:

- **mdi** Cartridge filter (all product types and sizes)

Background:

In our efforts to rationalize production process, certain changes have been made in identification of serial number on cartridge filters (all product types and sizes).

Need for Change:

To rationalize our production process, the cartridge serial number has been modified from two digits to three digits.

How Does It Affect the User:

From a practical point of view the changes done will only modify the cartridge serial number. All performance specifications, product attributes, product packaging and traceability remain unchanged.

- **Identification on cartridge filter** – Cartridge serial number is now printed in three digits which was earlier printed in two digits.

IDENTIFICATION ON CARTRIDGE FILTERS		
	OLD	NEW
Printing on cartridge filters	CT0413-13	CT0413-013
Description	CT0413 is cartridge Lot No. and 13 is two digit cartridge serial No.	CT0413 is cartridge Lot No. and 013 is three digit cartridge serial No.

Effective Date:

mdi Cartridge Filters manufactured after November 15, 2013 will be with said changes.

In case you have any queries, please feel free to contact our Technical Support Team at 'mdi@vsnl.com'.



Quality Assurance
Advanced Microdevices Pvt. Ltd.