

Reference : CHANGE NOTICE No. MB/10/01

Subject : Change in Grid Printing of mdi Grid Marked Membrane Disc Filters

Scope : This change notification will affect the following products:

- GCN (All pore sizes)
- GCN – S (All pore sizes)
- RGCN (All pore sizes)
- RGCN – S (All pore sizes)

Background:

In our efforts to provide the user better products with more convenience and better performance, certain changes have been made in the **mdi** Grid marked Membrane Disc Filters.


Need for Change:

As above.

How Does It Affect the User:

From a practical point of view the changes done will improve the ease-of-use under actual conditions of use. All performance specifications and product attributes remain unchanged.

The intensity of colour of grid lines has been changed to lighter one to reduce the interference of grid lines with bacterial colonies. Also, the grid lines are now with dotted lines instead of solid lines earlier. These changes will improve the usability of the grid marked membrane disc filters.



Effective Date:

mdi Grid Marked Membrane Disc Filters (all types and all pore sizes) manufactured after December 03, 2009 will have said changes.

In case you have any queries, please feel free to contact our Technical Support Team at 'mdi@vsnl.com'.

