

Reference : CHANGE NOTICE No. MB/10/03

Subject : Change in mdi Closed Sterility Test Device type Stericheck-LVP (All membrane types)

Scope : This change notification will affect the following component of mdi Stericheck:

- Needle - New improved design

Background:

In our efforts to provide the user better products with more convenience, certain changes have been made in the needle of mdi Closed Sterility Test Device type **Stericheck-LVP**.


Need for Change:

As above.

How Does It Affect the User:

From a practical point of view this change will improve the needle grip. All performance specifications and product attributes remain unchanged.

- **Stericheck LVP needle:** A new integral design of needle and Y connection has been made which will eliminate the assembling and jointing and further reinforce the needle grip.



Effective Date:

mdi Closed Sterility Test Device type **Stericheck-LVP** (all membrane types) manufactured after December 27, 2010 will have said changes.

In case you have any queries, please feel free to contact our Technical Support Team at 'mdi@vsnl.com'.

Quality Assurance

Advanced Microdevices Pvt. Ltd.