

Change Notice No.	Notification date	Implementation date
CN/A/20/003	30 <sup>th</sup> October 2020	12 <sup>th</sup> November 2020

**Subject: Change in the capacity of Bags used in MDI Single Use Assemblies and AseptiBag Gold Assemblies.**

**Scope: This change notification will affect the products with the following catalog number (where X can be any alphabet or numeral):**

- ABXXXXXXXXXXXXX
- AXXXXXXXXXXXXXX

**Background:** In our efforts to provide the user more robust products, it has been decided to make change in the capacity of 20 liters and 50 liters bags used in **mdi** Single Use Systems and AseptiBag Gold Assemblies with catalog numbers as mentioned above.

**Need for Change:** Existing 20 liters and 50 liters bag used in MDI Single Use Assemblies and AseptiBag Gold Assemblies are having liquid holding capacity of 21 and 50 liters respectively. In order to meet the market requirements and to provide the safety margin in the bags, the liquid holding capacity of 20 and 50 Liters bags has been increased to 24 and 55 liters, respectively. Hence, to incorporate the same, the length of the 20 and 50 liters bags has been increased from 49 cm to 52 cm and from 77.5 cm to 86 cm respectively in MDI Single use Assemblies and AseptiBag Assemblies with catalog numbers as mentioned above.

**How Does It Affect the User:** From a practical point of view, this change only increases the length and thus volume holding capacity of the 20 liter and 50 liter Bags to be used in MDI Single Use Assemblies and AseptiBag Gold Assemblies. All performance specifications, product attributes, traceability and lot release criteria will remain unchanged and possess no impact on the product quality. All packaging and labeling attributes such as packaging box, pack size and content of the product package label remain same as before, only length of the 20 liter and 50 liter bags has been increased.

**Effective Date:** **mdi** AseptiBag Gold and bag assemblies as mentioned above, manufactured after the implementation date **(12<sup>th</sup> November 2020)** will have the said changes. However, the available stock of MDI Single use Assemblies and AseptiBag Gold assemblies having Bags with less liquid holding capacity will be received by you till the stock's last.

In case you have any queries, please feel free to contact our Technical Support Team at [info@mdimembrane.com](mailto:info@mdimembrane.com)

