

# Pre-loaded ASESS® Aseptic Sampling System with Sampling Port Cartridge Adapter





## Pre-loaded Disposable Sampling Port Cartridge Adapter

**mdi** pre-loaded disposable sampling port cartridge with adapter does away with tedious, time consuming assembly of ASESS® sampling system ports in the cartridge and subsequently in the multiport adapters. This provides an easy, faster, ready to use sampling system which fits directly on to the vessel port.

## Aseptic Sampling Applications

Application Areas	Fermentation & Cell Culture	Buffer/Media Prep	Purification	Fill & Finish	Cold Chain
Cell Viability, Density and Characterization	✓	-	-	-	✓
Fluid Chemistry and Gas Analysis	✓	✓	✓	✓	✓
pH Analysis	✓	✓	✓	✓	✓
IgG Analysis	✓	✓	✓	✓	✓
Osmolality	✓	✓	✓	✓	✓
Sterility and Bioburden Testing	✓	✓	✓	✓	✓
Endotoxin Testing	✓	✓	✓	✓	✓

## Special Features

- Closed system for maximum operator safety and zero contamination
- Easier, faster and completely disposable, ready-to-use solution for sampling application
- No tedious time consuming assemblies of sampling system ports in cartridges and multi-port adaptors

# Sampling Container Options



Sampling Bags



Tamper Proof Sampling Bottles





Centrifuge Tubes




Tamper Proof Sampling Syringes

## How to Use



# Product Realization



## **Advanced Microdevices Pvt. Ltd.**

20-21, Industrial Area, Ambala Cantt- 133 006, India  
Tel: +91 - 171-2699290, 2699471  
E-mail: info@mdimembrane.com  
Website: www.mdimembrane.com

## **mdi Membrane Technologies INC.**

75 Utley Drive STE 103, Camp Hill, PA 17011  
E-mail: rs@mdimembranetech.com  
Phone: +1-717-412-0943, Fax: +1-717-695-9637  
Website: www.mdimembranetech.com