

Product	Combs
Intended Use	Manufacture of Immunodiffusion Tests. Particularly suited for semi quantitative analysis. Combs with 8 or 12 legs fit into a normal ELISA plate. It allows 8 or 12 samples to be tested at the same time.
Comb Type	CCN-12 Leg, CCN-8 Leg, CCN-II 10 Leg

Physical Properties

Property	DCN-II		
Appearance	White		
Wetting with Water	Pad Size (mm)	Volume	Time
	4	1 µl	≤5 sec
	5	1 µl	≤5 sec
Adhesion	The membrane does not comes off the backing after immersion in water for 24 hours.		

Supply Format

Type	Pad Size (mm)	Comb Size (mm)	Leg Spacing (mm)	No. of Pads
CCN-12 Leg	5	104 X 35	3.7	1
CCN-8 Leg	5	69 X 35	3.7	1
CCN-II, 10 Leg	4	95 X 82	5	Max. 8

Storage Condition

Store and handle in accordance with all current regulations and standards.

Store in original packing at 1-35 °C. Avoid direct sunlight.

Caution: Material should be kept away from heat radiators.

Shelf Life

Two Year from the date of dispatch, stored under above conditions.

Ordering Information

Product	Description			
Type	Comb Length	Comb Width	No. of Membrane Pads	Quantity (Nos.)
CCN-12 Leg	As per requirement	As per requirement	As per requirement	As per requirement
CCN-8 Leg				
CCN-II, 10 Leg				

Example:

CCN-12 Leg	104 mm	35 mm	1	-- Nos.
------------	--------	-------	---	---------