

<b>Product</b>	Combs
<b>Intended Use</b>	Manufacture of Immunodiffusion Tests. Particularly suited for semi quantitative analysis. Combs with 8 or 12 legs fit into a normal ELISA plate. It allows 8 or 12 samples to be tested at the same time.
<b>Comb Type</b>	CCN-12 12 Leg, CCN-8 8 Leg, CCN-II 10 Leg

**Physical Properties**

Property	CCN-II		
Appearance	White		
Wetting with Water	Pad Size (mm)	Volume	Time
	4	1 µl	≤5 sec
	5	1 µl	≤5 sec
Adhesion	The membrane does not comes off the backing after immersion in water for 24 hours.		

**Supply Format**

Type	Pad Size (mm)	Comb Size (mm)	Leg Spacing (mm)	No. of Pads
CCN-12 12 Leg	5	104 X 35	3.7	1
CCN-8 8 Leg	5	69 X 35	3.7	1
CCN-II, 10 Leg	4	95 X 82	5	Max. 8

**Storage Condition**

Store and handle in accordance with all current regulations and standards.

Store in original packing at 15-20°C. Avoid direct sunlight.

**Caution:** Material should be kept away from heat radiators.

**Shelf Life**

One Year from the date of dispatch, stored under above conditions.

DST DCNXXX1505G

**DATASHEET**

Ordering Information				
Product	Description			
Type	Comb Length	Comb Width	No. of Membrane Pads	Quantity (Nos.)
CCN-12 12 Leg	As per requirement	As per requirement	As per requirement	As per requirement
CCN-8 8 Leg				
CCN-II, 10 Leg				

**Example:**

CCN-12 12 Leg	104 mm	35 mm	1	— Nos.
---------------	--------	-------	---	--------