

MDI cylindrical drum containers are robust, cost effective systems, specially designed and qualified for single use MDI *AseptiLiners* used for preparation, storage and transfer of media, buffer etc. in biopharmaceutical manufacturing. These plastic drums, constructed of high density polyethylene material, have smooth surfaces which are clean room compatible and are available with volumes of 50L, 100L, 200L and 500L.

The drums can be placed on specially designed MDI SS trolleys for safe and easy transfer.

MDI Liner Drums have applicability in large volume liquid handling operations such as storage, mixing and transfer.

Technical Specifications

Materials of construction

Drums : HDPE

Trolleys : Stainless Steel

Key Features

- Ergonomic design for ease of use
- High Density Polyethylene drum construction
- Multi use application for shipping and storage
- Wheeled trolley available for portability
- 50L to 500L size option availability

Ordering Information

For Cylindrical Drums

Size	Catalog Number	
	with Drain Port	without Drain Port
50 L	PEDXHXXXXXXXXX01	PEDXHXXXWDXXX01
100 L	PEDXIXXXXXXXXX01	PEDXIXXXWDXXX01
200 L	PEDXJXXXXXXXXX01	PEDXJXXXWDXXX01
500 L	PEDXKXXXXXXXXX01	PEDXKXXXWDXXX01



For Wheeled SS Trolleys

Size	Catalog Number
50 L	TSSXHXXXXXXXXX01
100 L	TSSXIXXXXXXXXX01
200 L	TSSXJXXXXXXXXX01
500 L	TSSXKXXXXXXXXX01

Liner Drum with Trolleys

Size	Catalog Number	
	with Drain Port	without Drain Port
50 L	PEDXHTXXXXXXXXX01	PEDXHTXXWDXXX01
100 L	PEDXITXXXXXXXXX01	PEDXITXXWDXXX01
200 L	PEDXJTXXXXXXXXX01	PEDXJTXXWDXXX01
500 L	PEDXKTXXXXXXXXX01	PEDXKTXXWDXXX01