

Product CoverTapes

Intended Use Manufacture of Lateral Flow Immunoassay Test Kit. The Clear Tape is used to cover the membrane so that the membrane does not dry up and at the same time the user can see the lines that develop. The Masking Taps are used to hold the pads and membrane together with the plastic backing material.

Type Masking Tape (MT1 and MT2)
ClearTape CT

Physical Properties

Property	MT1	MT2	CT
Appearance	White	White	Transparent
Thickness	180 - 280 µm		80 - 140 µm
Adhesion	The tape does not comes off the membrane plastic after immersion in water for 24 hours.		

Supply Format

	Masking Tape	Clear Tape
Format	Strip	Strip/Reel
Dimensions	Desired Dimensions	

Storage Condition Store and handle in accordance with all current regulations and standards.
Store in original packing at 15-30 °C. Avoid direct sunlight.
Caution: Material should be kept away from heat radiators.

Shelf Life One Year from the date of dispatch, stored under above conditions.

Ordering Information			
Product	Description		
Type	Strip/Reels Length	Strip/Reels Width	Quantity (Nos.)
*MT1	As per requirement	**As per requirement	As per requirement
*MT2	* MT1 and MT2 is only available in strip format ** Maximum strip/reels width is 5cm		
CT			

Example: ↓ ↓ ↓ ↓ ↓

MT2	30 cm	3 cm	-
-----	-------	------	---