



Membrane Technologies

# Steridisc Dispenser Model SD-01

**mdj** Steridisc Dispenser is an advanced, light weight table top instrument for dispensing individually packed membrane disc filters from the pre-sterilized reels.

This all-new fully automatic dispenser with electronic control is very easy to operate. The push-button system makes disc handling very convenient for the user. It delivers a single sterile membrane and simultaneously wraps the peeled-off protective packaging.

The stainless steel clips hold the free end of the packaging material (paper on one side and the plastic film on the other side), while releasing the membrane disc filter, absolutely untouched by hand.

A special feature of the Steridisc Dispenser is operation by foot paddle, which allows totally hands free operation for the microbiologist.



## Specifications

### Dimensions

Height	: 270 mm
Width	: 175 mm
Length	: 250 mm

**Weight** : 3 Kg

**Voltage** : 100-240V AC

**Frequency** : 50 Hz

**Dispense Speed** : ~0.5 Second per disc

**Storage** : It is recommended to be stored in a clean environment at 10° to 40° C

## Applications

- Membrane dispensing for sterility testing and microbiology

## Advantages

- Automatic single unit dispensing
- Less handling time
- No accidental contamination
- Easy loading of membrane roll
- Easy portability: Small and light
- Push-button or foot paddle operated dispensing
- Easy to clean and maintain

## Ordering Information

Type	
<b>mdj</b> Steridisc Dispenser	Model SD-01

To Order please Specify:		
Type	Model No.	Qty.