


Bovine Whole Blood Separator (Type FRB)

Point of Care Testing for Cattle

Property	FRB (0.45mm)
Appearance	White
Tare weight (mg/cm ²)	5.0 ± 1
Thickness	460 ± 65 µm




Added 20µl of bovine blood and 80µl of assay buffer

Canine Whole Blood Separator (Type FRC)

Point of Care Testing for Dogs

Property	FRC (0.35mm)	FRC (0.6mm)
Appearance	White	White
Tare weight (mg/cm ²)	7.2 ± 2	10.5 ± 3
Thickness	385 ± 65 µm	570 ± 80 µm

Competition **mdi** FRC (0.35 mm)



RBC Containment


Added 20µl of canine blood and 100µl of assay buffer

Equine Whole Blood Separator (Type FRE)

Point of Care Testing for Horse

Property	FRE (0.35mm)	FRE (0.6mm)
Appearance	White	White
Tare weight (mg/cm ²)	7.2 ± 2	10.5 ± 3
Thickness	385 ± 65 µm	570 ± 80 µm

mdi FRE (0.35 mm) **mdi** FRE (0.6 mm)



Added 30µl of equine blood and 60µl of assay buffer

Supply Format

Format	Strip/ Sheets/Reels
Dimensions	Sheets: Dimension as desired upto 25 cm x 30 cm Reels: 15 – 20 M Reels of desired width

Storage Condition

Store and handle in accordance with all current regulations and standards. Store in original packing at 1-35 °C. Avoid direct sunlight.

Caution: Material should be kept away from heat radiators.

Shelf Life

Two Year from the date of slitting, stored under above conditions.

Ordering Information

Product	Description			
Type	Thickness (mm)	Strip/Sheet/Reels Length	Strip/Sheet/Reels Width	Quantity (Nos.)
Type FRB (Bovine)	0.45	As per requirement	As per requirement	As per requirement (Numbers/Meters)
Type FRC (Canine)	0.35			
	0.6			
Type FRE (Equine)	0.35			
	0.6			

Example:

FRC	0.6	25 cm	30 cm	___Nos.
-----	-----	-------	-------	---------

Please contact for any animal blood related query.

mdi Membrane Technologies INC.

75 Utlely Drive STE 103
Camp Hill, Pennsylvania 17011
United States of America
Website: www.mdimembranetech.com

ADVANCED MICRODEVICES PVT. LTD.

20-21, Industrial Area, Ambala Cantt- 133 006, India
Tel: +91 - 171-2699290, 2699471
E-mail: diagnostics@mdimembrane.com
Website: www.mdimembrane.com