

Certificate of Quality

The **Medi G-M blood Genomic DNA microprep kits** have been manufactured in **mdi** facility in compliance with **ISO 9001** regulations using **validated production processes**.

Medi G-M blood Genomic DNA microprep kit

Catalog No. : BRGKXXXXXXXXP0250
Type : BRGK
Lot Number : BR1-3073-199

LOT RELEASE CRITERIA

Binding capacity of column : $\geq 10 \mu\text{g}$

Typical Yield :

Leucocyte count	Yields
1 X 10 ⁶	2.3 μg
6.4 X 10 ⁶	6.7 μg
14.8 X 10 ⁶	10.4 μg

DNA Purity (A₂₆₀ /A₂₈₀) : 1.7 – 1.9

DNase Free : Yes

Volume of blood = 200 μl

Elution volume = 60 - 200 μl

VALIDATED FOR

- Downstream Processing
 - Restriction Digestion
 - PCR
 - Sequencing
 - Cloning
- Shelf Life
 - Minimum one year at Room Temperature
 - More than two years at 2-8°C
- Purification of small as well as large gDNA fragments
- Consistent intra lot (column to column) and inter lot genomic DNA recovery

CUSTOMER SUPPORT

mdi offers its unique interdisciplinary skills to provide solutions to specific problems. Please contact our factory or the local application specialist.



Head of Quality Control

Issue Date: 18-Dec-13

Advanced Microdevices Pvt. Ltd.

21, Industrial Area, Ambala Cantt, INDIA,

Tel: +91-171-2699290/ 2699274

Fax: +91-171-2699221/ 2699008

Website: www.mdimembrane.com

Email: mdi@vsnl.com

An ISO 9001 Company