

mdi BioKart™ is a specially designed stainless steel trolley with adjustable holder stand for 2D single use bags. These can hold 5L - 20L AseptiBag™ Gold 2D bag assemblies and are fully adjustable to a maximum height of 72 inches. The trolley has separate capsule filter holder to hold mdi large capsule filters of different sizes.

These trolleys support storage bag assemblies placed on the head and hang down from a stainless steel hook. The bottom of the trolley incorporates SS tray for holding accessories, tubing etc. The trolley is equipped with two front locking wheels for safety, when not in motion.

BioKart™ trolleys are made from high quality stainless steel 316L which are easy to clean and disinfect. The mirror polished contact surface is smooth to ensure safety of the bags.

Key Features

- Ergonomic access
- Enhanced mobility with strong wheels
- Lockable wheels for safety
- Convenient hanging
- Adjustable height to handle up to 20 liter bags

Technical Specifications

Materials of construction

SS316L

Dimensions

Height: 70 inch (179 cm)

Length: 25 inch (62.5 cm)

Width: 18 inch (46 cm)



Wheels with locking system



Clamp for height adjustment

Ordering Information

To order, please specify:

BKSSTAXXXXXX01