

Product	CoverTapes
Intended Use	Manufacture of Lateral Flow Immunoassay Test Kit. The Clear Tape is used to cover the membrane so that the membrane does not dry up and at the same time the user can see the lines that develop. The Masking Taps are used to hold the pads and membrane together with the plastic backing material.
Type	Masking Tape (MT1 and MT2) Clear TapeCT

Physical Properties	Property	MT1	MT2	CT
	Appearance	White	White	Transparent
	Thickness	230±50 µm		110±30 µm
	Adhesion	The tape does not comes off the membrane plastic after immersion in water for 24 hours.		

Supply Format		Masking Tape	Clear Tape
	Format	Strip	Strip/Reel
	Dimensions	Desired Dimensions	

Storage Condition Store and handle in accordance with all current regulations and standards.
Store in original packing at 1-35 °C. Avoid direct sunlight.
Caution: Material should be kept away from heat radiators.

Shelf Life Two Year from the date of dispatch, stored under above conditions.

Ordering Information			
Product	Description		
Type	Strip/Reels Length	Strip/Reels Width	Quantity (Nos.)
*MT1	As per requirement	**As per requirement	As per requirement
*MT2	* MT1 and MT2 is only available in strip format ** Maximum strip/reels width is 5cm		
CT			

Example: ↓ ↓ ↓ ↓ ↓

MT2	30 cm	3 cm	-
-----	-------	------	---