

Certificate of Quality

The M – Funnel Nxt has been manufactured in a **mdi** facility in compliance with **ISO 9001** regulations using **validated production processes**.

M – Funnel Nxt

Catalog No.	: FXCN0902XXXX212
Type	: M - Funnel Nxt
Pore Size	: 0.45 µm
Lot Number	: 9965L
Ster. No.	: 969
Date of Sterilization	: 2025 - 12
Expiry Date	: 2028 - 12

SPECIFICATION

Membrane	: Cellulose Nitrate Membrane with Ink grid
Plastic Components	
Funnel Cover	: Acrylic
Funnel & Base	: Polypropylene
Plug	: LDPE
Filter Diameter	: 47 mm
Funnel capacity	: 100 ml
Sterilization	: Pre sterilized with EO

LOT RELEASE CRITERIA

Typical Water Flow Rate	: 250 ml/min at Vacuum of 250 mm of Hg
Microbial Recovery of <i>E. coli</i>	: Exhibits greater than 70% recovery for <i>E. coli</i> (ATCC 8739).

VALIDATED FOR

Sterility	: The sterilization process has been validated to assure Sterility Assurance Level (SAL) of 10^{-6} in accordance with ISO 11135.
Microbial Recovery Test	: Passes Test as per ISO 7704: Water Quality – Evaluation of membrane filter used for microbiological analysis.
Inhibitory Effect of Ink Grids	: Passes test as per ASTM D 4200-82: Standard test method for evaluating inhibitory effect of ink grids on membrane filters.
Package Integrity Test	: Passes as per ASTM F1929-98: Standard test method for detecting seal leaks in porous medical packaging by Dye Penetration Test.
EO Residual Level	: ≤ 6 ppm

T. No.: COQ/FUN/016-00



Head of Quality Assurance
Issue Date: 08-Dec-25

Advanced Microdevices Pvt. Ltd.
Jawahargarh Road,
Village-Tepla, Ambala, INDIA.
Tel: +91-171-2699290/2699274
Website: www.mdimembrane.com
Email: info@mdimembrane.com

An ISO 9001 Company