

# Certificate of Quality

The Absorbent Pad Type AP has been manufactured in **mdi** facility in compliance with ISO-9001 regulations using validated production processes.

## ABSORBENT PAD

**CATALON NO.** : P080XXX021X260X  
**TYPE** : AP080  
**LOT NO.** : AQ305034D-125

### PHYSICAL TESTS

	RESULTS	SPECIFICATION
Thickness (µm)	: 733	800±150
Tare weight (mg/cm <sup>2</sup> )	: 16.0	16±6

### FUNCTIONAL TESTS

Water Holding Capacity (µl/cm <sup>2</sup> )	: 67.2	82±30
--	--------	-------

### DIMENSIONS

Length (mm)	: 260	260±1.0
Width (mm)	: 21	21±0.2

**SLITTING/PACKAGING DATE** : 20-Apr-2024

**BEST BEFORE** : 19-Apr-2026

### STORAGE CONDITIONS

Temperature : 1-35°C

Note : Absorbent pad Type AP should be stored in a cool, dry place in original packing and should be kept away from direct sun light and heat radiators.

T.No.: COQ/ABP/002-01



Head of Quality Assurance

Issue Date: 20-Apr-24

**Advanced Microdevices Pvt. Ltd.**  
21, Industrial Area, Ambala Cantt, INDIA,  
Tel: +91-171-2699290/ 2699274  
Website: [www.mdimembrane.com](http://www.mdimembrane.com)  
Email: [info@mdimembrane.com](mailto:info@mdimembrane.com)

**An ISO 9001 Company**