

# Certificate of Quality

The Conjugate Release Matrix Type PT has been manufactured in **mdi** facility in compliance with ISO-9001 regulations using validated production processes.

## CONJUGATE RELEASE MATRIX

**CATALOG NO** : PTR5XXX014X300X  
**TYPE** : PT-R5  
**LOT NO.** : PV020314D

### PHYSICAL TESTS

	RESULTS	SPECIFICATION
Thickness ( $\mu\text{m}$ )	: 393	405 $\pm$ 80
Tare weight ( $\text{mg}/\text{cm}^2$ )	: 9.3	9.5 $\pm$ 2.5

### FUNCTIONAL TESTS

H <sub>2</sub> O Wicking Rate, 4cm distance (sec)	: 28	26 $\pm$ 8
Accelerated Gold Release Test 45°C/24 hours	: Pass	No visible mark at point of application

### DIMENSIONS

Length (mm)	: 300	300 $\pm$ 1.0
Width (mm)	: 14	14 $\pm$ 0.2

**SLITTING/PACKAGING DATE** : 19-Apr-2024

**BEST BEFORE** : 18-Apr-2026

### STORAGE CONDITIONS

Temperature : 1-35°C

Note : Conjugate Release Matrix Type PT should be stored in a cool, dry place in original packing and should be kept away from direct sun light and heat radiators.

T.No.: COQ/RLM/002-01



Head of Quality Assurance

Issue Date: 19-Apr-24

**Advanced Microdevices Pvt. Ltd.**  
21, Industrial Area, Ambala Cantt, INDIA,  
Tel: +91-171-2699290/ 2699274  
Website: [www.mdimembrane.com](http://www.mdimembrane.com)  
Email: [info@mdimembrane.com](mailto:info@mdimembrane.com)

**An ISO 9001 Company**