

# Certificate of Quality

The **SVP Dilutor** has been manufactured in a **mdi** facility in compliance with **ISO 9001** regulations.

## SVP DILUTOR

Catalog No. : SVPDXXXXXXXXX209  
Type : **SVP Dilutor**  
Lot No. : SD5557G  
Ster. No. : 012  
Date of Sterilization : 2017 - 07  
Expiry Date : 2020 - 07

### MATERIALS

Tubes : PVC

### DIMENSIONS

Tube Length : 850mm

Sterilization : Pre-sterilized with EO gas using validated sterilization process with well controlled Relative Humidity, Temperature, Pressure and Exposure Time.

### LOT RELEASE TESTS

Each lot has been sampled, tested and released for the following characteristics:

- o Integrity of joints
- o Sterility

### VALIDATED FOR

Maximum Operating Temperature : 35°C Continuous  
Residual EO : < 6 ppm

### PRECAUTIONS

1. **mdi** SVP Dilutor is a pre-sterilized device and its sterility will be compromised in case the pack is damaged or tampered.  
Do not use the device in such cases.
2. During handling, avoid contamination of inlet.

### CUSTOMER SUPPORT

**mdi** offers its unique interdisciplinary skills to provide solutions to specific problems. Please contact our factory or the local application specialist.



Head of Quality Assurance

Issue Date: 24-Jul-17

**Advanced Microdevices Pvt. Ltd.**  
21, Industrial Area, Ambala Cantt, INDIA,  
Tel: +91-171-2699290/ 2699274  
Website: [www.mdimembrane.com](http://www.mdimembrane.com)  
Email: [info@mdimembrane.com](mailto:info@mdimembrane.com)

**An ISO 9001 Company**