

# Certificate of Quality

The E – Funnel has been manufactured in a **mdi** facility in compliance with **ISO 9001** regulations using **validated production processes**.

## E - Funnel

Catalog No. : FMEN0902XXXXX212  
Type : E - Funnel  
Pore Size : 0.45 µm  
Lot Number : 5597G  
Ster. No. : 045  
Expiry Date : 2020 - 07

### SPECIFICATION

Membrane : Edge Hydrophobic Cellulose Nitrate Membrane  
Plastic Components  
    Funnel Cover : Acrylic  
    Funnel & Base : Polypropylene  
    Plug : LDPE  
Filter Diameter : 47 mm  
Funnel capacity : 100 ml  
Sterilization : Pre sterilized with EO

### LOT RELEASE CRITERIA

**Water Flow Rate** :  $\geq$  200 ml/min at Vacuum of 250 mm of Hg  
**Sterility** : Samples passed the sterility test in accordance with U.S Pharmacopoeia.

### VALIDATED FOR

**Biosafety** : Passes biological test for plastics Class VI as per USP and complies with Code of Federal Regulations FDA 21CFR 177.1520(a)1(i).

### CUSTOMER SUPPORT

**mdi** offers its unique interdisciplinary skills to provide solutions to specific problems. Please contact our factory or the local application specialist.



Head of Quality Assurance

Issue Date: 24-Jul-17

**Advanced Microdevices Pvt. Ltd.**  
21, Industrial Area, Ambala Cantt, INDIA,  
Tel: +91-171-2699290/ 2699274  
Website: [www.mdimembrane.com](http://www.mdimembrane.com)  
Email: [info@mdimembrane.com](mailto:info@mdimembrane.com)

**An ISO 9001 Company**