

# Certificate of Quality

The Sample Pad Type GFB has been manufactured in **mdi** facility in compliance with ISO-9001 regulations using validated production processes.

## SAMPLE PAD

**CATALOG NO** : PGF4D35014X260X  
**TYPE** : GFB-R4(0.35)  
**LOT NO.** : GL244444D-101

### PHYSICAL TESTS

|   | RESULTS | SPECIFICATION |
|---|---------|---------------|
| Thickness ( $\mu\text{m}$ )             | : 379   | 370 $\pm$ 70  |
| Tare weight ( $\text{mg}/\text{cm}^2$ ) | : 5.1   | 5.5 $\pm$ 1.5 |

### FUNCTIONAL TESTS

|   |        |  |
|---|--------|--|
| H <sub>2</sub> O Wicking Rate, 4cm distance (sec) | : 52   | 50 $\pm$ 20                                |
| Accelerated Gold Release Test<br>45 °C/24 hours   | : Pass | No visible mark at<br>point of application |

### DIMENSIONS

|             |       |               |
|-------------|-------|---------------|
| Length (mm) | : 260 | 260 $\pm$ 1.0 |
| Width (mm)  | : 14  | 14 $\pm$ 0.2  |

### SLITTING/PACKAGING DATE

: 16-Apr-2024

### BEST BEFORE

: 15-Apr-2026

### STORAGE CONDITIONS

Temperature : 1-35°C

Note : Sample Pad Type GFB should be stored in a cool, dry place in original packing and should be kept away from direct sun light and heat radiators.

T.No.: COQ/SMP/002-01



Head of Quality Assurance

Issue Date: 17-Apr-24

**Advanced Microdevices Pvt. Ltd.**  
21, Industrial Area, Ambala Cantt, INDIA,  
Tel: +91-171-2699290/ 2699274  
Website: [www.mdimembrane.com](http://www.mdimembrane.com)  
Email: [info@mdimembrane.com](mailto:info@mdimembrane.com)

**An ISO 9001 Company**