

Certificate of Quality

The Microglassfibre Disc Filters have been manufactured in a **mdi** facility in compliance with **ISO 9001** regulations using **validated production processes**.

Microglassfibre Disc Filters

Catalog No. : GFSX1615XXXX103
Type : GFS
Pore Size : 2.0 µm
Diameter : 142 mm
Lot Number : GH887877G

LOT TEST DATA

Thickness : 330 - 380 µm

Water Flow Rate : ≥ 900 ml/min/cm² @ 0.70 Kg/cm²

Heat stability : Maintains integrity after autoclaving at 121 °C of 30 minute.

EXPIRY : 3 years when stored in a cool and dry place.



Head of Quality Assurance

Issue Date: 24-Jul-17

Advanced Microdevices Pvt. Ltd.
21, Industrial Area, Ambala Cantt, INDIA,
Tel: +91-171-2699290/ 2699274
Website: www.mdimembrane.com
Email: info@mdimembrane.com

An ISO 9001 Company