

# Certificate of Quality

The Masking Tape Type MT has been manufactured in **mdi** facility in compliance with ISO-9001 regulations using validated production processes.

## MASKING TAPE

**CATALOG NO.** : TMT1XXX036X300X  
**TYPE** : MT1  
**LOT NO.** : MT308924C-066

### PHYSICAL TESTS

Backing Film

: B.O.P.P

Thickness of Backing Film ( $\mu\text{m}$ )

: 30

Thickness of Glue ( $\mu\text{m}$ )  
(excluding backing film)

: 101

Tare weight of Glue ( $\text{gm}/\text{cm}^2$ )  
(excluding backing film)

: 11.1

### DIMENSIONS

Size (mm)

: 36x300

### SLITTING/PACKAGING DATE

: 27-Mar-2024

### BEST BEFORE

: 26-Mar-2026

### STORAGE CONDITIONS

Temperature

: 1-35°C

Note : Masking Tape Type MT should be stored in a cool, dry place in original packing and should be kept away from direct sun light and heat radiators.

T.No.: COQ/CVT/002-01



Head of Quality Assurance

Issue Date: 27-Mar-24

**Advanced Microdevices Pvt. Ltd.**  
21, Industrial Area, Ambala Cantt, INDIA,  
Tel: +91-171-2699290/ 2699274  
Website: [www.mdimembrane.com](http://www.mdimembrane.com)  
Email: [info@mdimembrane.com](mailto:info@mdimembrane.com)

**An ISO 9001 Company**