

Certificate of Quality

The Nitrocellulose Membrane Type CLW has been manufactured in a **mdi** facility in compliance with ISO-9001 regulations using validated production processes.

NITROCELLULOSE MEMBRANE

TYPE : CLW-040-SH34
PORE SIZE : 0.45 µm
LOT NO. : LB316163

A. PHYSICAL TEST DATA

	VALUES	SPECIFICATION
Thickness (µm)	: 496	500 ± 70 µm
Bubble point (psi)	: 39.2	36 ± 6 psi

B. FUNCTIONAL TEST DATA

Wettability	: Passed Test	
Protein binding with green LED(3mm device)	: 1.158	≥0.9
Time to absorb 30µl buffer in 3mm hole (sec.)	: 26	23 ± 5

DIMENSIONS

Size (mm)	: 240x240	240 x 240 (±5) mm
-----------	-----------	-------------------



Head of Quality Control

Advanced Microdevices Pvt. Ltd.

21, Industrial Area, Ambala Cantt, INDIA,
 Tel: +91-171-2699290/ 2699274/ 2699471
 Fax: +91-171-2699221/ 2699008
 Website: www.mdimembrane.com
 Email: mdi@vsnl.com