

## **Certificate of Quality**

The Microglassfibre Disc Filters have been manufactured in a **mdi** facility in compliance with **ISO 9001** regulations using **validated production processes**.

## **Microglassfibre Disc Filters**

Catalog No. : GF2X1914XXXX103

Type : GF2

Pore Size : 1.5 μm

Diameter : 293 mm

Lot Number : GG857577G

## **LOT TEST DATA**

Thickness :  $400 - 500 \mu m$ 

Water Flow Rate :≥ 1100 ml/min/cm² @ 0.70 Kg/cm²

**Heat stability** : Maintains integrity after autoclaving at 121 °C of 30 minute.

**EXPIRY** : 3 years when stored in a cool and dry place.

Om

Head of Quality Assurance Issue Date: 24-Jul-17 **Advanced Microdevices Pvt. Ltd.** 

21, Industrial Area, Ambala Cantt, INDIA, Tel: +91-171-2699290/ 2699274 Website: www.mdimembrane.com Email: info@mdimembrane.com

An ISO 9001 Company